

V79 COLONY FORMING ASSAY

Experiment Name : ^{100% labeled cells} $^{131}\text{IUdR}$ toxicity (crossed dose); Exp. #: 1.1; Investigator: A. Bishayee
 Date: 01/29/98

1. Set the ~~rocker~~ ^{roller} at 37°C incubator with 5% CO₂, set the Coulter Counter, wash cells (from two 150 cm² flusk, subcultured 1:2, 24h before) with PBS, trypsinize cells, resuspend in 10 ml MEMB for each flusk, pool, vortex, pass five times through 3 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml isotone (Coulter balanced electrolyte solution)
2. Dilute to ~40,00,000 cells/ml in MEMB (final volume 11 ml) [Actual count : 4518666.6 cells/ml cells/ml)
3. Transfer 1 ml of cell suspension into ten 12 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall
4. Roll the tubes for 3-4 h at 37°C, 5% CO₂ Date/Time: 01/29/98; 2-30 p.m.
5. Prepare MEMB containing radioactivity in hood
 3.5 µl $^{131}\text{IUdR}$ (prepared on 01/29/98) + 3 ml MEMB
6. After 3-4 h, remove test tubes from roller and add MEMB with or without radioactivity according to Table below. Date/Time: 01/29/98; 6-15 p.m.

| Tube # | $^{131}\text{IUdR}$ uCi/ml | Cells in MEMB (ml) | MEMB (ml) | MEMB+ $^{131}\text{IUdR}$ (ml) [2.0 uCi/ml] |
|--------|-------------------------------|--------------------------|--------------|---|
| 1 | 0 | 1.0 | 1.0 | 0 |
| 2 | 0 | 1.0 | 1.0 | 0 |
| 3 | 0.010 | 1.0 | 0.990 | 0.010 |
| 4 | 0.025 | 1.0 | 0.975 | 0.025 |
| 5 | 0.050 | 1.0 | 0.950 | 0.050 |
| 6 | 0.100 | 1.0 | 0.900 | 0.100 |
| 7 | 0.250 | 1.0 | 0.750 | 0.250 |
| 8 | 0.500 | 1.0 | 0.500 | 0.500 |
| 9 | 1.000 | 1.0 | 0.000 | 1.000 |

7. Return test tubes to roller for 12 h. Date/Time: 01/29/98; 6-30 p.m.

8. While test tubes are rolling label 40 (4x10) gamma-tubes (13 X 100 mm VWR glass test tube)
9. After ~12 h incubation period, remove tubes and centrifuge at 2000 rpm at 4°C for 10 min
(precooled centrifuge). Date/Time: 01/30/98; 9-00 a.m.
10. During centrifugation, obtain ice
12. Remove buckets from centrifuge and carefully remove 100 µl of supernatant and place in pre-labeled gamma-tube.
13. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
14. Centrifuge tubes for 10 min at 2000 rpm, 4°C
15. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
16. Centrifuge tubes for 10 min at 2000 rpm, 4°C
17. Decant supernatant, click tubes, vortex, resuspend in 10 ml MEMA
18. Centrifuge tubes for 10 min at 2000 rpm, 4°C
19. Decant supernatant, click tubes, resuspend in 200 ul ice cold MEMA ,transfer the cell suspension in polypropylene microcentrifuge tubes with attached caps (Helena Plastics, 400 ul) using pipet tips
20. Again add 200 ul ice cold MEMA, resuspend and transfer the cell suspensions in the same polypropylene microcentrifuge tubes (Total volume ~400 ul)
21. Centrifuge tubes for 5 min at 1000 rpm, 4°C
22. Transfer tubes at 10°C for 72 h. Date/Time: 01/30/98; 11-30 a.m.
23. Transfer 10 ul supernatant in three sets of tubes containing small pieces of tissue paper from 100 ul supernatant removed earlier and count them for radioactivity
Date/Time: 01/30/98; 11-45 a.m.
24. After 72 h, carefully remove the supernatant from the top, resuspend pellet in 200 ul wash MEMA and transfer the content to ten 12 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using pasteur pipet
Date/Time: 02/02/98; 9-45 a.m.
25. Again add 200 ul wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions in 12 ml tubes
26. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)
27. Labeling and preparation of dilution tubes and colony dishes
 - load 69 60 mm petri dishes with 4 ml MEMA
 - load 36 T-tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.
28. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
29. Centrifuge tubes for 10 min at 2000 rpm, 4°C
30. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
31. Centrifuge tubes for 10 min at 2000 rpm, 4°C

32. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 3 cc syringe with 21 gauge needle
33. Determine cell concentration by transferring 100 μ l to Coulter cup
34. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5 and transfer 0.5 ml to tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3 and vortex tube X.3 and transfer 0.5 ml to tube X.2. Keep tubes on ice.
35. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.
36. Transfer 300 μ l of cell suspension (in triplicate) to gamma tubes for each tube
37. Incubate petridishes for 1 week
38. Count gamma tubes for radioactivity Date/Time : 02/02/98; 2-45 p.m.
39. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol.
Stain colonies with 0.05% crystal violet
40. Count colonies. There must be between 25 and 250 colonies for the flask to be a valid data point.

Expt # 1-1

Take two 150cm² flasks containing 80-90% confluent cells
Suspend in 10 ml MEMB for each flask.

01/29/98

Initial cell count = 1111, 1127, 1151 (Counting time = 1.14 Sec)
Avg. cell count = 1129.6 (Manometer correct = 50 μ)
Cell conc. = $1129.6 \times 4000 = 4518666.6$ cells/ml

For dilution,

vol. of original cell suspension taken = $\frac{44000000}{\quad}$

Take \quad ml of cell suspension + \quad ml MEMB = 11 ml.

After dilution,

Final count =
Avg. Count =
Cell conc. =

943
968
975
956

Expt # 101

01/29/98

Preparation of ^{131}I UDR in MEMB (2 $\mu\text{Ci}/\text{ml}$)

Stock ^{131}I UDR 1.71 $\mu\text{Ci}/\text{ml}$

Prepare 3 ml of 2 $\mu\text{Ci}/\text{ml}$ \Rightarrow 6 μCi required

From Stock 1.71 μCi \Rightarrow 1 μl

6 μCi \Rightarrow $\frac{1}{1.71} \times 6$

\Rightarrow 3.50 μl .

- i) Pipet 3.50 μl ^{131}I UDR from Stock in hood
- ii) Keep it at RT for ~2-3 h
- iii) Add 3 ml of MEMB

01/27/98

~~wt of chloramid T tube = 5.19567
 wt of " + chl. T = 5.18300 g

 chl. T = 0.0126 = 12.67 mg~~

~~Ph. buffer = $\frac{12.67}{2} = 6.33 \text{ ml}$~~

~~wt of Pot. metabisulfite tube = 5.19136 g
 " " " + Pot. = 5.24000

 Pot. metabisulfite = 0.04864~~

~~dH₂O = $\frac{48.64}{10} = 4.86 \text{ ml}$~~

~~wt of 2'-deoxyurine tube = 5.25464
 wt of 2'-deoxy = 5.26438

 2'-Deoxy = 0.00974~~

~~Ph. buffer = $\frac{9.74}{2} = 4.88 \text{ ml}$~~

0.02

2 — 1 ml
 1 — $\frac{1}{2}$
 0.02

TABLE-1

Expt. #: 101

Date/Time: 01/30/98; 11-45 a.m.

| Tube # | Medium count for 10 ul (cpm) | Avg. cpm | dpm [cpm/0.142] | μ Ci/ml (A) on counting [dpm/22200] | μ Ci/ml (A ₀) on addition [A/e ^{-λt}] |
|--------|---------------------------------|----------|--------------------|---|---|
| 1 | 3, 2, 4 | 0 | 0 | 0 | 0 |
| 2 | 4, 5, 3 | 0 | 0 | 0 | 0 |
| 3 | 17, 23, 20 | 20 | 140.84 | 0.0063 | 0.006 0.0067 |
| 4 | 48, 36, 47 | 43.66 | 307.51 | 0.0138 | 0.012 0.0146 |
| 5 | 73, 89, 68 | 76.66 | 539.90 | 0.0243 | 0.022 0.0258 |
| 6 | 163, 172, 139 | 158 | 1112.67 | 0.0501 | 0.047 0.0532 |
| 7 | 378, 384, 406 | 389.33 | 2741.78 | 0.1235 | 0.1161 0.1313 |
| 8 | 763, 847, 818 | 809.33 | 5699.53 | 0.2567 | 0.2414 0.2730 |
| 9 | 1587, 1607, 1625 | 1606.33 | 11312.20 | 0.5095 | 0.4743 0.5420 |
| 10 | | | | | |

0.01
0.02
0.05
0.1
0.25
0.5
1

$$\begin{aligned}
 & e^{-\lambda t} \\
 = & e^{-\frac{\ln 2}{t_{1/2}} \times 17.25 \text{ h}} \\
 = & e^{-\frac{0.693}{193.2} \times 17.25} \\
 = & e^{-0.0618} \\
 = & \cancel{1.063} \quad 0.940
 \end{aligned}$$

$$A_0 = \frac{Ae}{1.063} = 0.940$$

TABLE-2

Expt. # : 1-1

Date/Time : 02/02/98; 2-45 p.m.

| Tube # | Radioactivity for 300 ul cell suspension (cpm) | Avg. cpm | dpm [cpm/0.142] | μ Ci/ml (A) on counting [dpm/666000] | μ Ci/ml (A ₀) after 12 h incubation [A _t /e ^{-λt}] |
|--------|--|----------|-----------------|--------------------------------------|---|
| 1 | 0, 2, 1 | 0 | 0 | 0 | 0 |
| 2 | 2, 2, 0 | 0 | 0 | 0 | 0 |
| 3 | 105, 106, 105 | 105.33 | 741.78 | 0.00113 | 0.00149 |
| 4 | 250, 261, 278 | 263 | 1852.11 | 0.00278 | 0.00367 |
| 5 | 433, 475, 433 | 447 | 3147.88 | 0.00472 | 0.00624 |
| 6 | 1064, 1108, 1091 | 1087.66 | 7659.6 | 0.01150 | 0.01521 |
| 7 | 2562, 2552, 2472 | 2528.66 | 17807.51 | 0.02673 | 0.03535 |
| 8 | 4694, 4883, 4924 | 4833.66 | 34039.85 | 0.05111 | 0.06760 |
| 9 | 9588, 9790, 9829 | 9735.66 | 68561.03 | 0.10294 | 0.13616 |
| 10 | | | | | |

$$A_0 = A_t / e^{-\lambda t}$$

$$= A_t / e^{-\frac{\ln 2}{195.2} \times 77.75}$$

$$= A_t / e^{-\frac{0.693 \times 77.75}{195.2}}$$

$$= A_t / e^{-0.278}$$

$$= A_t / 0.756$$

$$= 77.75 \text{ h} \left\{ \begin{array}{l} 02/02/98 : 2-45 \text{ p.m.} \\ 01/30/98 : 9-00 \text{ a.m.} \end{array} \right.$$

TABLE-3

Expt. # : 1.1

Date/Time : 02/02/98; 11-00 a.m.

| Tube # | Coulter count for 100 ul cell suspension (Manometer select : 50 μ l) | Avg. count | Cells/ml [Avg. count x 400 ρ] | pCi/cell [$\frac{\mu\text{Ci}/\text{cell} \times 10^6}{\text{Cells/ml}}$] |
|--------|---|------------|---|--|
| 1 | 290, 310, 324 | 308 | 1232000 | 0 |
| 2 | 514, 536, 510 | 520 | 2080000 | 0 |
| 3 | 539, 535, 487 | 520.33 | 2081333.3 | 0.0007 |
| 4 | 484, 464, 469 | 472.33 | 1889333.3 | 0.00194 |
| 5 | 438, 436, 422 | 432 | 1728000 | 0.00361 |
| 6 | 499, 486, 503 | 496 | 1984000 | 0.00766 |
| 7 | 517, 537, 493 | 515.66 | 2062666.6 | 0.01713 |
| 8 | 439, 464, 467 | 456.66 | 1826666.6 | 0.03700 |
| 9 | 442, 429, 440 | 437 | 1748000 | 0.07789 |
| 10 | | | | |

uci/ml

0.01

0.025

0.05

0.1

0.25

0.5

1.0

0.00071

0.5 -

1.0 -

Expt # : 1.1 Date : 02/09/98

Colony Counts and Survival Fraction


| Tube.dilution | Colony 1 | Colony 2 | Colony 3 | Avg Colony | SF |
|---------------|----------|----------|----------|------------|--------|
| 1.2 | 44 | 37 | 44 | 41.66* | — |
| 2.2 | 118 | 107 | 113 | 112.66 | — |
| 3.2 | 73 | 91 | 93 | 85.66 | 0.7603 |
| 4.2 | 71 | 78 | 66 | 71.66 | 0.6360 |
| 5.2 | 69 | 71 | 68 | 69.33 | 0.6153 |
| 6.2 | 60 | 61 | 62 | 61 | 0.5414 |
| 7.2 | 55 | 45 | 60 | 53.33 | 0.4733 |
| 8.2 | 44 | 54 | 53 | 50.33 | 0.4467 |
| 9.2 | 28 | 25 | 31 | 28 | 0.2485 |
| | | | | | |

* not considered.

Expt #1.1

130. 100R tonicity

NATIONAL
12 183
MADE IN U.S.A.



Semi-Logarithmic
3 Cycles x 10 to the right


pci/cell

Expt # 1.1

131I dR toxicity

NATIONAL
12 183

SF



Semi-Logarithmic
3 Cycles x 10 to the inch

mCi