



## V79 COLONY FORMING ASSAY

**Experiment Name :** <sup>3</sup>HTdR toxicity (cluster, 50% labeling)  
**Investigator:** R Howell

**Exp. # :** 4  
**Date:** 6/28/01

**Serum/Lot #s:** Gibco Lot # 1023298 for MEMA and MEMB

1. Set the rocker-roller at 37°C incubator with 5% CO<sub>2</sub>, set the Coulter Counter, wash cells (from one 225 cm<sup>2</sup> flasks, subcultured 1:2 the day before) with PBS-PS, trypsinize cells with 2 ml trypsin 3 min at 37°C, resuspend in 10 ml MEMB, pool, pass five times through 10 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml Isotone II (Coulter balanced electrolyte solution).

2. Dilute to ~2,000,000 cells/ml in MEMB [Actual count :  $2.20 \times 10^6$  cells/ml]

$\left. \begin{array}{l} 2584 \\ 2507 \\ 2602 \end{array} \right\}$	$\left. \begin{array}{l} 100 \text{ ml cells} \\ 100 \text{ ml counter} \\ \text{setting} \end{array} \right\}$	$5.13 \times 10^6 / \text{ml}$ $18 \text{ ml} = 7 \frac{36 \times 10^6}{5.13 \times 10^6} = 7.0 \text{ ml cells} + 11 \text{ ml MEMB}$	$\left. \begin{array}{l} 5562 \\ 5411 \\ 5562 \end{array} \right\} 2.20 \times$
---	---	---	---

3. Transfer 1 ml of cell suspension into 20 14 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall

4. Keep the tubes in the roller for 3-4 h at 37°C, 5% CO<sub>2</sub>

**Date/Time:** 2:00pm

5. Prepare MEMB containing radioactivity in hood

$\mu\text{l } ^3\text{HTdR (Stock : } 1 \mu\text{Ci}/\mu\text{l on } 5/30/01) + \text{ ml MEMB}$

Manufacturer: NEN NET-0272 Lot #: 3106-427 Calibration:

6. After 3-4 h, remove first set of 10 test tubes from roller and add MEMB with or without radioactivity according to Table below. Also add 1 ml MEMB to each of second set of tubes.

Tube #	<sup>3</sup> HTdR μCi/ml	Cells in MEMB (ml)	MEMB (ml)	MEMB+ <sup>3</sup> HTdR [160μCi/ml] (ml)
1	0	1.0	1.0	0
2	0	1.0	1.0	0
3				
4				
5	10	1.0	0.975	0.125
6	20	1.0	0.75	0.25
7	30	1.0	0.625	0.375
8	40	1.0	0.5	0.5
9	60	1.0	0.25	0.75
10	80	1.0	0	1.0

3ml

Didn't have this  
much. Had 350-400ml  
Ringed stock viral  
to get all of it.  
Concentrations will be  
lower than expected.

3.2 ml  
at 160 μCi/ml

512 μl H-3

7. Return test tubes to roller for 12-14 h.
8. Next day, while test tubes are in roller label 10 gamma-tubes (13 X 100 mm VWR glass test tube)
9. After ~12-14 h incubation period, remove all tubes and centrifuge at 2000 rpm at 4°C for 10 min (precooled-centrifuge).
10. Remove buckets from centrifuge and carefully remove 150 μl of supernatant from tubes containing radioactivity and place in pre-labeled gamma-tubes.
11. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
12. Centrifuge tubes for 10 min at 2000 rpm, 4°C
13. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
14. Centrifuge tubes for 10 min at 2000 rpm, 4°C
15. ~~Decant supernatant, click tubes, vortex, resuspend in 10 ml of wash MEMA~~
16. ~~Centrifuge tubes for 10 min at 2000 rpm, 4°C~~
17. Decant supernatant, click tubes, vortex, resuspend unlabeled cells in 5 ml MEMA.
18. Transfer unlabeled cells to the corresponding tubes containing 2,000,000 labeled cells.
19. Wash unlabeled cell tube with 5 ml MEMA and transfer to corresponding labeled cell tube.
20. ~~Syringe the pooled cells 5 times with 5 ml syringe with 21 G needle.~~
21. Centrifuge tubes for 10 min at 2000 rpm, 4°C.
22. Decant supernatant completely, click tubes, vortex.
23. Transfer the cell suspension in polypropylene microcentrifuge tubes with attached caps (Helena Plastics, 400 μl) using 200 μl pipette tip.
24. Again add 200 μl ice cold MEMA, resuspend and transfer the cell suspensions in the same

Date/Time: 06/28/01; 6-00pm

Date/Time: 8:00 am 6/29/01

polypropylene microcentrifuge tubes (Total volume ~400  $\mu$ l)

25. Centrifuge tubes for 5 min at 1000 rpm, 4°C

26. Transfer tubes at 10.5°C for 72 h. (Lebline shaker w/o shaking) **Date/Time:** 10:10 am 6/29/01

27. Transfer 10  $\mu$ l supernatant in three sets of 7 ml scintillation vials and add 6 ml liquid scintillation cocktail (Ecoscint) from 150  $\mu$ l supernatant removed earlier and count them for radioactivity **Date/Time:** 10:00 am 6/29/01

28. After 72 h, carefully remove the supernatant from the top, resuspend pellet in 200  $\mu$ l wash MEMA and transfer the content to ten 12 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using Pasteur pipette **Date/Time:** 10:10 am 7/2/01

29. Again add 200  $\mu$ l wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions in 12 ml tubes

30. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)

31. Labeling and preparation of dilution tubes and colony dishes

- load 60 mm tissue culture dishes with 4 ml MEMA

- load 40 sterile tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.

32. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

33. Centrifuge tubes for 10 min at 2000 rpm, 4°C

34. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

35. Centrifuge tubes for 10 min at 2000 rpm, 4°C

36. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 5 cc syringe with 21 gauge needle

37. Determine cell concentration by transferring 100  $\mu$ l to Coulter cup

38. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5, transfer 0.5 ml into dilution tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3, vortex tube X.3 and transfer 0.5 ml to tube X.2 and vortex. Keep tubes on ice.

39. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.

40. Transfer 200  $\mu$ l of cell suspension (in triplicate) to 7  $\mu$ l scintillation vial containing 6 ml cocktail (Ecolume)

41. Incubate tissue culture dishes for 1 week

42. Count vials for radioactivity

**Date/Time:** 7/3/01 17:02

43. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol. Stain colonies with 0.05% crystal violet

44. Count colonies. There must be between 25 and 250 colonies for the dish to be a valid data point.

→ Transfer  $2 \times 10^6$  cells to T75 flask w/ 12.5 ml MEMA (same lot) for HPRT expression

7/9/01 RWH

Tube	Dilution	Colony Counts
<del>1</del>	200	180, 182, 207
2	200	178, 134, 179
5	200	125, 109, 133
6	200	81, 81, 91
7	200	113, 121, 105
8	200	117, 100, 104
9	200	113, 89, 112
10	200	72, 78, 72

Rogev's exp. #

6

V79, H1dR,

TABLE-3

4 6/29/01

Expt. #: ~~██████████~~

Date/Time: July 2, 2001

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 400] - Bckgr.	<del>pCi/cell [uCi/ml] x 10<sup>6</sup> Cells/ml]</del>
1	9779, 9910, 9659		3 908 666	0.51
2	9996, 10014, 10072		4 006 533	0.50
3	—NO SAMPLE			→
4	—NO SAMPLE			→
5	9369, 9064, 8939		3 644 533	0.56
6	9687, 9505, 9647		3 840 800	0.53
7	9080, 9335, 9210		3 678 933	0.54
8	9613, 9596, 9431		3 814 267	0.53
9	8609, 8538, 8421		3 412 667	0.59
10	8744, 8699, 8716		3 483 467	0.57

Bckgr - 11  
Mode - 500  $\mu$ l

↑  
2x10<sup>6</sup>/T75 flask  
for Hprt expression  
(day 0)

Exp-4 R Howell

USER: 6 ID:H3 HOWELL      PRESET TIME: 1.00      TUE 03 JUL 2001 17:02  
 SAMPLE REPEAT: 1 CYCLE REPEAT: 1 SCR:N      RS232:N  
 H 1 ADC:N DCF:N RCM:N  
 CHANNEL 1-LL: 0 UL: 400 2SIGMA: 2.00 BKG SUB: 0.00 BKG 2SIG: 0.00 LSR: 0  
 DATA CALC: CPM, UNKNOWN REPLICATES: 1      NORM FACTOR:0 1.00000  
 HALF LIFE(DAYS):N

SAM	POS	CH	CPM	2SIG%	TIME	EL TIME	AVG H#	ERR
1	**	1	12.00	57.74	1.00	1.48	86.0	
2	**	2	13.00	55.47	1.00	3.06	85.0	
3	**	3	11.00	60.30	1.00	4.63	85.0	
4	**	4	16.00	50.00	1.00	6.21	85.0	
5	**	5	11.00	60.30	1.00	7.78	85.0	
6	**	6	13.00	55.47	1.00	9.34	84.0	
7	**	7	184813.33	1.20	0.15	10.05	86.0	
8	**	8	186820.00	1.19	0.15	10.77	84.0	
9	**	9	182226.66	1.21	0.15	11.48	85.0	
10	**	10	362040.00	1.05	0.10	12.20	84.0	
11	**	11	355233.31	0.87	0.15	12.93	86.0	
12	**	12	357940.00	0.86	0.15	13.65	85.0	
13	**	13	550046.62	0.70	0.15	14.39	86.0	
14	**	14	541413.31	0.70	0.15	15.13	85.0	
15	**	15	543220.00	0.70	0.15	15.88	85.0	
16	**	16	718720.00	0.61	0.15	16.63	84.0	
17	**	17	723886.62	0.61	0.15	17.39	84.0	
18	**	18	724306.62	0.61	0.15	18.15	85.0	
19	**	1	1086426.62	0.50	0.15	18.99	85.0	
20	**	2	1101986.62	0.49	0.15	19.78	84.0	
21	**	3	1099560.00	0.49	0.15	20.58	85.0	
22	**	4	1445173.25	0.43	0.15	21.39	84.0	
23	**	5	1426913.25	0.43	0.15	22.22	83.0	
24	**	6	1413553.25	0.43	0.15	23.04	85.0	
25	**	7	11.00	60.30	1.00	24.67	102.0	
26	**	8	10.00	63.25	1.00	26.23	101.0	
27	**	9	6.00	81.65	1.00	27.80	103.0	
28	**	10	5.00	89.44	1.00	29.41	102.0	
29	**	11	7.00	75.59	1.00	30.97	101.0	
30	**	12	12.00	57.74	1.00	32.53	103.0	
31	**	13	22020.00	1.91	0.50	33.59	104.0	
32	**	14	24328.89	1.91	0.45	34.60	102.0	
33	**	15	24695.55	1.90	0.45	35.61	104.0	
34	**	16	51160.00	1.98	0.20	36.36	104.0	
35	**	17	51735.00	1.97	0.20	37.11	105.0	
36	**	18	50395.00	1.99	0.20	37.92	103.0	
37	**	1	70573.33	1.94	0.15	38.67	102.0	
38	**	2	72293.33	1.92	0.15	39.37	102.0	
39	**	3	62770.00	1.79	0.20	40.12	103.0	
40	**	4	98573.33	1.64	0.15	40.82	102.0	
41	**	5	94066.66	1.68	0.15	41.52	103.0	
42	**	6	88240.00	1.74	0.15	42.23	102.0	
43	**	7	136593.33	1.40	0.15	42.95	103.0	
44	**	8	145240.00	1.36	0.15	43.66	105.0	
45	**	9	145073.33	1.36	0.15	44.37	103.0	
46	**	10	187586.66	1.19	0.15	45.07	103.0	

1M  
2M  
5M  
6M  
7M  
8M  
9M  
10M  
1C  
2C  
5C  
6C  
7C  
8C  
9C

30 ML medium

200 ML cells

SAM	POS	CH	CPM	2SIG%	TIME	EL TIME	AVG H#	ERR
●	**	-11	1	183713.33	1.20	0.15	45.79 104.0	} 100
48	**	-12	1	185440.00	1.20	0.15	46.51 104.0	
49	**	-13	1	28762.86	1.99	0.35	47.47 0.0	



Parameters

Date	6/28/2001
Experiment No.	4
Investigator	R. Howell
Cell Line	V79 - Anupam's 10/00 freeze
Modifier	None
Radionuclide	H-3
Half-life (days)	4500.45
Radiation Yield	1
Radiochemical	3HTdR
Manufacturer/Lot	NEN NET-027Z Lot 3106-427
Original Calibration Date/Time	5/30/2001 12:00
Present Calibration Date/Time	6/28/2001 0:00
Fraction of Cells Labeled	0.5
Liquid Scintillation Cocktail	Ecoscint
Volume of LSC Cocktail (ml)	6
Volume/Type Counting Vial	7 ml
Model of Counter	Beckman
Counting Efficiency	0.65
Activity Added (Date/Time)	6/28/2001 18:00
Cells Washed (Date/Time)	6/29/2001 8:00
Medium Tubes Counted (Date/Time)	7/3/2001 17:02
Cell Tubes Counted (Date/Time)	7/3/2001 17:02
Vol. Supernatant Counted (µl)	30
Vol. Suspension Counted Cell Activity (µl)	200
Vol. Suspension Coultter (µl)	100
Coultter Manometer Volume (µl)	500
Average Coultter Background Counts	11
Coultter Calibration Parameter	400
Hemocytometer Counting (Yes or No)?	No
Original Activity Concentration (MBq/ml)	37
Time Elapsed Since Original Calibration (d)	9
Present Activity Concentration (MBq/ml)	36.95
I-125=59.408, H-3=4500.45, Po-210=138.376, I-131=8.02	
I-125=1.47, H-3=1.0, Po-210=1.0, I-131=8.02	
Time Elapsed Between Add and Wash (hr)	14.00
Time Elapsed Between Add and Count (hr)	59.00
Time Elapsed Between Wash and Count (hr)	59.00
Background	
Coultter 1	11
Coultter 2	11
Coultter 3	11

MediumActivity

Experiment: 4  
Date: 6/28/2001

Tube #	1st	2nd	3rd	Medium count (CPM)	CPM Average	CPM corrected for control	DPM CPM/(y e)	At $\mu\text{Ci/ml}$ on counting	Ao $\mu\text{Ci/ml}$ at addition	Ao kBq/ml at addition
1	12	13	11	13	0	0	0	0	0	0
2	16	11	13	13	0	0	0	0	0	0
3				#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
4				#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
5	184813	186820	182226	184620	184607	284011	4.2644	4.2660	157.8435	
6	362040	355233	357940	358404	358392	551372	8.2789	8.2820	306.4336	
7	550046	541413	543220	544893	544880	838277	12.5867	12.5915	465.8860	
8	718720	723886	724306	722304	722291	1111217	16.6849	16.6913	617.5768	
9	1086426	1101986	1099560	1095991	1095978	1686120	25.3171	25.3267	937.0880	
10	1445173	1426913	1413553	1428546	1428534	2197744	32.9992	33.0117	1221.4312	

CellSuspension

Experiment: 4  
Date: 06/28/01

Tube #	Suspension count (CPM)			CPM Average	CPM corrected for control	DPM CPM/(y e)	A <sub>i</sub> μCi/ml on counting	A <sub>o</sub> μCi/ml after uptake	A <sub>o</sub> kBq/ml after uptake
	1st	2nd	3rd						
1	11	10	6	9	0	0	0.0000	0	0.0000
2	5	7	12	0	0	0	0.0000	0	0.0000
3				#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
4				#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
5	22020	24328	24695	23681	23673	36419	0.08203	0.08206	3.0361
6	51160	51735	50395	51097	51088	78597	0.17702	0.17709	6.5522
7	70573	72293	62770	68545	68537	105441	0.23748	0.23757	8.7901
8	98573	94066	88240	93626	93618	144027	0.32439	0.32451	12.0068
9	136593	145240	145073	142302	142294	218913	0.49305	0.49323	18.2497
10	187586	183713	185440	185580	185571	285494	0.64300	0.64325	23.8002

CoulterSurvival

Experiment: 4  
 Date/Time: #####



Tube #	Coulter count			Average Cells/ml	Hemocytometer Count in Grid			
	1st	2nd	3rd		1st	2nd	3rd	4th
1	9779	9910	9659	9783	3908667			
2	9996	10014	10072	10027	4006533			
3				#DIV/0!	#DIV/0!			
4				#DIV/0!	#DIV/0!			
5	9369	9004	8934	9102	3636533			
6	9687	9505	9647	9613	3840800			
7	9080	9335	9210	9208	3678933			
8	9613	9596	9431	9547	3814267			
9	8609	8598	8421	8543	3412667			
10	8744	8699	8716	8720	3483467			

Tube #	Predicted # Cells Seeded	Actual # Cells Seeded	Colony count			Average	PE (%)	SF	SF
			1st	2nd	3rd			Uncorrected	Corrected
1	200	391	180	182	207	185	46.796	1.00	1.0000
2	200	401	178	179					
3	200	#DIV/0!				#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
4	200	#DIV/0!				#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
5	200	364	125	109	133	122	33.640	0.6605	0.7189
6	200	384	81	81	91	84	21.957	0.4554	0.4692
7	200	368	113	121	105	113	30.715	0.6102	0.6564
8	200	381	117	100	104	107	28.053	0.5778	0.5995
9	200	341	113	89	112	105	30.670	0.5652	0.6554
10	200	348	72	78	72	74	21.243	0.3996	0.4540

Summary

Experiment: 4  
 Date/Time: 6/28/2001

Tube #	Activity Conc. (kBq/ml)	Activity/Cell (mBq/cell)	Survival Uncorrected	Survival Corrected
1	0.000	0.000	1.0000	1.0000
2	0.000	0.000		
3	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
4	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
5	157.844	1.670	0.6605	0.7189
6	306.434	3.412	0.4554	0.4692
7	465.886	4.779	0.6102	0.6564
8	617.577	6.296	0.5778	0.5995
9	937.088	10.695	0.5652	0.6554
10	1221.431	13.665	0.3996	0.4540



7/31/01

Scoring of Plating Efficiency  
for mutation study

Exp # 4

1	180	133	167	
2	140	140	179	
3				
4				
5	143	149	169	
6	<del>140</del>	180	175	162
7	141	147	162	
8	123	136	124	
9	112	140		
10	106	105	105	

HPRT expression → HPRT selection  
 Eugene's exp #4, V79, HDR, 50%  
 Day 0 | 3 | 6 | 9 | 11 | 13 | Chall  
 Date 7/2 | 7/5 | 7/8 | 7/11 | 7/13 | 7/15 | 7/17  
 Fix, count = 27 July

Colony Counts - Mutant Selection  
Exp. 4

1	1	2	0	2	2	1	0	1	2	1
2	0	1	1	0	0	0	0	0	0	0
3										
4										
5	1	2	2	0	1	0	0	0	0	1
6	1	1	4	1	2	4	2	1	0	1
7	1	1	1	2	1	2	2	3	2	
8	3	3	4	3	0	0	0	2	0	4
9	0	2	3	0	1	2	2	6	3	
10	3	5	5	6	3	4	0	0	3	