

V79 COLONY FORMING ASSAY

Experiment Name : ³HTdR toxicity (cluster, 100% labeling)

Exp. #: 5

Investigator: R. Howell

Date: 7/16/01

Serum/Lot #s: Gibco Lot# 1023288
FCS

1. Set the rocker-roller at 37°C incubator with 5% CO₂, set the Coulter Counter, wash cells (from one 225 cm² flasks, subcultured 1:2 the day before) with PBS-PS, trypsinize cells with 2 ml trypsin 3 min at 37°C, resuspend in 10 ml MEMB, pool, pass five times through 10 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml Isotone II (Coulter balanced electrolyte solution).

2. Dilute to ~4,000,000 cells/ml in MEMB [Actual count : 4.05×10^6 cells/ml)

$$\left. \begin{array}{l} 10135 \\ 10089 \\ 10152 \end{array} \right\} 10125 \times 400 = 4.05 \times 10^6 \text{ cells/ml}$$

3. Transfer 1 ml of cell suspension into 10 14 ml tubes (Falcon polypropylene test tube, 17x100 mm) labeled 1-10 both on cap and wall.

4. Keep the tubes in the roller for 3-4 h at 37°C, 5% CO₂

Date/Time: 2:15 pm

5. Prepare MEMB containing radioactivity in hood

7/16/01

160 μ l ³HTdR (Stock : 1 μ Ci/ μ l on 7/5/01) + ml MEMB

Manufacturer: NEN NET-027E Lot #:

Calibration: 7/5/01

6. After 3-4 h, remove tubes from roller and add MEMB with radioactivity according to Table below.

Tube #	³ HTdR μCi/ml	Cells in MEMB (ml)	MEMB (ml)	MEMB+ ³ HTdR [40 μCi/ml] (ml)
1	0	1.0	1.0	0
2	0	1.0	1.0	0
3	1	1.0	0.95	0.05
4	2	1.0	0.90	0.10
5	3	1.0	0.85	0.15
6	4	1.0	0.8	0.20
7	5	1.0	0.75	0.25
8	10	1.0	0.5	0.5
9	15	1.0	0.25	0.75
10	20	1.0	0	1.0

Make
4 ml
160 μl H-3
+ 3.84 ml MEMB

8 ml 3 ml

7. Return test tubes to roller for 12-14 h. **Date/Time:** 07/16/01; 6-30 pm
8. Next day, while tubes are in roller label 10 gamma-tubes (13 X 100 mm VWR glass tube)
9. After ~12-14 h incubation period, remove all tubes and centrifuge at 2000 rpm at ~~4°C~~ for 10 min (precooled centrifuge). **Date/Time:**
10. Remove buckets from centrifuge and carefully remove 150 μl of supernatant from tubes containing radioactivity and place in pre-labeled gamma-tubes.
11. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
12. Centrifuge tubes for 10 min at 2000 rpm
13. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
14. Centrifuge tubes for 10 min at 2000 rpm
17. Decant supernatant, click tubes, vortex, resuspend in 10 ml MEMA.
18. Centrifuge tubes for 10 min at 2000 rpm.
22. Decant supernatant completely, click tubes, vortex.
23. Transfer the cell suspension in polypropylene microcentrifuge tubes with attached caps (Helena Plastics, 400 μl) using 200 μl pipette tip.
24. Again add 200 μl MEMA, resuspend and transfer the cell suspensions in the same polypropylene microcentrifuge tubes (Total volume ~400 μl)
25. Centrifuge tubes for 5 min at 1000 rpm, 4°C
26. Transfer tubes at 10.5°C for 72 h. **Date/Time:** 10:30am 7/17/01
27. Transfer 30 μl supernatant in three sets of 7 ml scintillation vials and add 6 ml liquid scintillation cocktail (Ecoscint) from 150 ul supernatant removed earlier and count them for radioactivity **Date/Time:** 10:30am 7/17/01

28. After 72 h, carefully remove the supernatant from the top, resuspend pellet in 200 μ l wash MEMA and transfer the content to ten 12 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using Pasteur pipette

Date/Time: 7/20/01

29. Again add 200 μ l wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions in 12 ml tubes

30. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (*precooled centrifuge*)

31. Labeling and preparation of dilution tubes and colony dishes

- load 60 mm tissue culture dishes with 4 ml MEMA

- load 40 sterile tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.

32. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

33. Centrifuge tubes for 10 min at 2000 rpm, 4°C

34. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

35. Centrifuge tubes for 10 min at 2000 rpm, 4°C

36. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 5 cc syringe with 21 gauge needle

37. Determine cell concentration by transferring 100 μ l to Coulter cup

38. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5, transfer 0.5 ml into dilution tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3, vortex tube X.3 and transfer 0.5 ml to tube X.2 and vortex. Keep tubes on ice.

39. Transfer 2×10^6 cells to T75 flask for mutant expression.

39. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.

40. Transfer 200 μ l of cell suspension (in triplicate) to 7 ml scintillation vial containing 6 ml cocktail (Ecolume)

41. Incubate tissue culture dishes for 1 week

42. Count vials for radioactivity

Date/Time:

43. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol. Stain colonies with 0.05% crystal violet

44. Count colonies. There must be between 25 and 250 colonies for the dish to be a valid data point.

For hprt
expression

TABLE-3

100% cluster, HDR, 179 cells

Expt. #: 5

July 20, 2001. Date/Time: ~12:00 am

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count - Backgr x 400] - Backgr Total in 2ml	pCi/cell [uCi/ml x 10 ⁶ Cells/ml] 2 x 10 ⁶ cells (ml)
1	4827 4878 4582	1896800	3793600	1.05
2	5348 5384 5312	2130400	4260800	0.96
3	4752 4848 4672	1894133	3788267	1.06
4	5120 4944 4934	1990933	3981867	1.00
5	5457 5625 5562	2210400	4420800	0.91
6	5265 5084 4725	2061000 2001733	4122000	0.97
7	4929 4878 4753	1932533	3865067	1.04
8	5115 4965 5094	2014400	4028800	0.95
9	4615 4711 4583	1,839733	3679466	1.08
10	4662 4596 4647	1845200	3690400	1.08

Backgr - 22.

Mode - 500 μ m

7/23/2001

Exp. 5

Replating of cells for expression. Trypsinize, syringe, count 100µl, transfer 2×10^6 cells to new flasks. T75 flasks

Vol (ml)	ml cells	Coulter counts	# of cells/ml	total $\times 10^6$	2×10^6 cells
6	6	8660 8632	3,46	20.74	0.58
6	6	8221 8653	3,37	20.23	0.59
6	6	8321 8571	3.38	20.26	0.59
6	6	7231 7155	2.88	17.25	0.69
6	6	7769 7770	3,11	18.63	0.64
<hr/>					
7	10.8	3920, 3990, 3970	3960	1.58×10^6	1.26
8	10.5	3793, 3876, 3873		1.56×10^6	1.28ml
9	7.7	4388, 4288, 4308	4328	1.73×10^6	1.16ml
<hr/>					
10	7.8	contaminated	4328	1.73×10^6	1.16ml
<hr/>					
		4728, 5001, 4888	4872	1.95×10^6 /ml	1.03ml

NOTE:

①. Sample 1-5 - Belgor-4
 Mode-500

②. Would not be surprise if any of 1-5 get contamination, because of syringe piston

Be sure to mention that you saw this ad on *Science Careers*

▲ PAGE TOP



Copyright © 2001 by the American Association for the Advancement of Science.

TABLE-3

Expt. #: 5

Date/Time: 7/24/01

only manom.
vol. 100µL

Tube #	Coulter count for 100 ul cell suspension Manometer vol. 100µL	Avg. count	Cells/ml [Avg. count x 400] x 5	pCi/cell $\frac{\mu\text{Ci/ml} \times 10^6}{\text{Cells/ml}}$ ml cells for 2×10^6
1	1533, 1634, 1624	1597	$3.19 \times 10^6/\text{ml}$	0.63
2	1642, 1653, 1705	1667	$3.33 \times 10^6/\text{ml}$	0.60
3	1246, 1154, 1201	1200	$2.40 \times 10^6/\text{ml}$	0.83
4	1699, 1627, 1667	1664	3.33×10^6	0.60
5	991, 994, 905	963	1.93×10^6	1.04
6	1257, 1323, 1224	1268	2.54×10^6	0.79
7	1397, 1350, 1405	1384	2.77×10^6	0.72
8	1485, 1383, 1367		2.82×10^6	0.71
9				
10	1381, 1331, 1364		2.75×10^6	0.73

Replating of cells for expression. Remove media from T75, wash 1X with PBS-PS, add 2 ml trypsin, set at room temp > 5 min. Resuspend in 8 ml wash MEMA, transfer to 12 ml tube, centrifuge 2000 rpm, decant, resuspend in 5 ml MEMA. Syringe 5X, 100µL to Coulter, coulter manometer volume 100µL. Transfer 2×10^6 cells to T150 flask containing 20 ml MEMA.

Split for ^{mutation} expression on 7/30/01

TABLE-3

Expt. #: 5

Date/Time: 7/30/01

Tube #	Coulter count for 100 ul cell suspension <i>500 ul manometer</i>	Avg. count	Cells/ml [Avg. count x 400]	pcf/cell <i>[uCi/ml x 10⁶]</i> Cells/ml <i>2 x 10⁶ cells</i>
1	6039, 5960, 6050	6016	2.41×10^6 /ml	0.83mR
2	7639, 7858, 7775	7757	3.10×10^6	0.65mR
3	<i>contaminated</i>			
4				0.75
5				0.75
6				0.75
7				0.75
8				0.75
9				
10	7333, 7150			0.75

Seed into T75 flasks

7/31/01

Colony Counting
for Survival

4

Tube	Dilution	Colony Counts
1	200	108 114 114
2	200	85 81
3	200	63 68 83
4	200	94 84 73
5	200	66 108 67
6	200	61 93 73
7	200	56 55 52
8	200	47 38 37
9	200	46 38 54
10	200	40 30 30

Exp. 5

USER: 6 ID:H3 HOWELL PRESET TIME: 1.00 MON 23 JUL 2001 13:22
 SAMPLE REPEAT: 1 CYCLE REPEAT: 1 SCR:N RS232:N
 1 AQC:N OCF:N RCM:N
 CHANNEL 1-LL: 0 UL: 400 2SIGMA: 2.00 BKG SUB: 0.00 BKG 2SIG: 0.00 LSR: 0
 DATA CALC: CPM, UNKNOWN REPLICATES: 1 NORM FACTOR: 1.00000
 HALF LIFE(DAYS):N

SAM	POS	CH	CPM	2SIG%	TIME	EL TIME	AVG H#	ERR
1	**	1	9.00	66.67	1.00	1.42	82.0	
2	**	2	8.00	70.71	1.00	3.00	82.0	
3	**	3	7.00	75.59	1.00	4.57	83.0	
4	**	4	5.00	89.44	1.00	6.13	81.0	
5	**	5	16.00	50.00	1.00	7.69	81.0	
6	**	6	8.00	70.71	1.00	9.26	82.0	
7	**	7	28785.71	1.99	0.35	10.16	82.0	
8	**	8	26352.50	1.95	0.40	11.17	81.0	
9	**	9	28608.57	2.00	0.35	12.07	83.0	
10	**	10	57790.00	1.86	0.20	12.82	82.0	
11	**	11	55970.00	1.89	0.20	13.57	83.0	
12	**	12	57645.00	1.86	0.20	14.32	83.0	
13	**	13	85580.00	1.77	0.15	15.02	84.0	
14	**	14	83353.33	1.79	0.15	15.72	83.0	
15	**	15	82146.66	1.80	0.15	16.43	82.0	
16	**	16	115919.99	1.52	0.15	17.13	83.0	
17	**	17	116026.66	1.52	0.15	17.83	83.0	
18	**	18	117066.66	1.51	0.15	18.54	83.0	
19	**	1	142806.66	1.37	0.15	19.29	82.0	
20	**	2	134140.00	1.41	0.15	20.01	81.0	
21	**	3	142500.00	1.37	0.15	20.72	82.0	
22	**	4	346220.00	0.88	0.15	21.44	81.0	
23	**	5	330266.66	0.90	0.15	22.16	82.0	
24	**	6	336800.00	0.89	0.15	22.89	80.0	
25	**	7	427873.31	0.79	0.15	23.63	82.0	
26	**	8	409359.97	0.81	0.15	24.35	82.0	
27	**	9	416326.66	0.80	0.15	25.09	82.0	
28	**	10	525273.31	0.71	0.15	25.83	82.0	
29	**	11	548100.00	0.70	0.15	26.58	83.0	
30	**	12	548013.31	0.70	0.15	27.32	83.0	
31	**	13	7.00	75.59	1.00	29.88	102.0	
32	**	14	9.00	66.67	1.00	30.44	101.0	
33	**	15	17.00	48.51	1.00	32.01	99.0	
34	**	16	7.00	75.59	1.00	33.57	101.0	
35	**	17	6.00	81.65	1.00	35.13	102.0	
36	**	18	7.00	75.59	1.00	36.71	102.0	
37	**	1	7.00	75.59	1.00	38.32	74.0	
38	**	2	2515.00	3.99	1.00	39.89	101.0	
39	**	3	2487.00	4.01	1.00	41.47	100.0	
40	**	4	5132.00	2.79	1.00	43.03	100.0	
41	**	5	4865.00	2.87	1.00	44.61	102.0	
42	**	6	4606.00	2.95	1.00	46.18	100.0	
43	**	7	9053.00	2.10	1.00	47.76	101.0	
44	**	8	8805.00	2.13	1.00	49.32	100.0	
45	**	9	8504.00	2.17	1.00	50.90	101.0	
46	**	10	11281.11	1.98	0.90	52.37	101.0	

1M
2M
3M
4M
5M
6M
7M
8M
9M
10M
1C
2C
3C
4C
5C

SAM	POS	CH	CPM	2SIGZ	TIME	EL TIME	AVG H#	ERR
	**	-11	1	11345.55	1.98	0.90	53.83	101.0
	**	-12	1	10826.32	1.97	0.95	55.36	102.0
49	**	-13					55.83	101
50	**	-14	1	13424.00	1.99	0.75	57.15	102.0
51	**	-15	1	13333.33	2.00	0.75	58.47	101.0
52	**	-16	1	41116.00	1.97	0.25	59.32	103.0
53	**	-17	1	39310.00	1.84	0.30	60.17	101.0
54	**	-18	1	36616.66	1.91	0.30	61.04	102.0
55	**	- 1	1	27582.50	1.90	0.40	62.10	94.0
56	**	- 2	1	20047.27	1.90	0.55	63.21	95.0
57	**	- 3	1	21732.00	1.92	0.50	64.27	96.0
58	**	- 4					64.73	101
59	**	- 5	1	50000.00	1.79	0.25	65.52	102.0
60	**	- 6	1	53260.00	1.94	0.20	66.28	101.0

6C
7C
8C
9C
10C

Control std. 29323 cpm

6-IG.

1. $(2.5 \text{ mg} + 10 \text{ ml MEMA}) \times 2 \text{ sample}$
Stock - $\boxed{5 \text{ mg} / 20 \text{ ml}}$
2. Dilute powder in each bottle (10 ml)
3. Syringe each one thru filter (0.2)
4. Collect all 6-IG into 50 ml conical tube
5. If it is your stock, sterile 6-IG solution.

Final concentration = $5 \mu\text{g} / \text{ml}$

- a) $10 \text{ ml} / \text{P}100 \Rightarrow 10 \text{ ml} \times 5 \mu\text{g} \Rightarrow 50 \mu\text{g} / \text{P}100 (10 \text{ ml})$
- b) you have already 7 ul in each P100
- c) $3 \text{ ml} + 50 \mu\text{g} \times \# \text{ of P}100\text{'s} \rightarrow \text{X}$

Example:

~~10 sample~~

$$\text{Need } 250 \text{ ml @ } \frac{50}{3} \text{ kg/ml} = 4167 \text{ kg}$$

$$\text{Stock} = \frac{5 \text{ mg}}{20 \text{ ml}} = 0.25 \frac{\text{mg}}{\text{ml}} = \frac{250 \text{ kg}}{\text{ml}}$$

$$\text{Need } 16.7 \text{ ml stock} + 233.3 \text{ mg MEMA}$$

MCA #1 - Canberra S100 -

Tag Number: 1028 Readout: Sat 14 Jul 2001 @ 15:04:45
Report Group: Full Acquire Started: Sat 14 Jul 2001 @ 14:48:03
Channel Size: 4096

Elapsed Live Time: 10.00 min.
Elapsed True Time: 10.03 min.
Dead Time: 0.27 %

ADC: Can1510,LLD=0.04,4k
Detector: Can CAM300 PIPS
Geometry: Am-241 disk against PIPS
Sample Description: Am-241 disk (Tracerlab)

MCA Mode: PHA+
Preset Conditions: Live Time = 10.00 min.

Area Statistics: Background Channels = 4
 % Error Sigma = 1.65

$$\frac{x(5.13 \times 10^6 \frac{\text{cells}}{\text{ml}})}{4.5 + x} = 2 \times 10^6 / \text{ml}$$

$$5.13 \times 10^6 x = 9 \times 10^6 + 2 \times 10^6 x$$

$$3.13 \times 10^6 x = 9 \times 10^6 = 2.87$$

1.3

8/2/01

Exp. 5

Processing for mutant selectionCoulter counts (100µl manometer, 100 µl cells)

				ml cells to add to 4.5 ml to get 2×10^6 /ml
1)	2503	2631	2542	
			$2559 \times 400 \times 5 = 5.12 \times 10^6$ /ml	2.9
2)	2411	2478	2399	
			$2429 \times 400 \times 5 = 4.86 \times 10^6$ /ml	3.1
4)	2616 5269	2806	2650	
			$2691 \times 400 \times 5 = 5.38 \times 10^6$ /ml	2.7
5)	2692	2730	2659	
			$2694 \times 400 \times 5 = 5.38 \times 10^6$ /ml	2.7
6)	2224	2289	2213	
			$2242 \times 400 \times 5 = 4.48 \times 10^6$ /ml	3.6
7)	2206	2335	2294	
			$2278 \times 400 \times 5 = 4.56 \times 10^6$ /ml	3.5
8)	2497	2317	2410	
			$2379 \times 400 \times 5 = 4.76 \times 10^6$ /ml	3.3
10)	2087	2138	2051	
			$2092 \times 400 \times 5 = 4.18 \times 10^6$ /ml	4.1

Preload: 80 100 mm dishes w/ 7 ml MEMA and put in incubator
 24 60 mm dishes w/ 4 ml MEMA "
 40 12 ml round bottom with 4.5 ml wash MEMA

Insufficient volume of stock.
 dilute with 4 ml instead of 4.5

Procedure:

Trypsinize flasks, ~~resusp~~ with 2 ml trypsin after 10 ml PBS PS wash.

Add 10 ml wash MEMA to flask & transfer to conical 15 ml tube. Centrifuge 4-5 min 2000 rpm. Decant, resusp in 5 ml MEMB. Syringe 5X with 21g, counter count 100µl, transfer aliquot to 4.5 ml wash MEMA so final

conc. is 2×10^6 cells/ml. Transfer 10 x 100 µl aliquots to 100 mm dishes. Perform 4 serial dilutions 0.5 → 4.5 ml and transfer 3 x 1 ml of final dilution to P60's.

$$\frac{x(4.18 \times 10^6)}{4+x} = 2 \times 10^6 \text{ /ml}$$

$$x = \frac{8 \times 10^6}{2.18 \times 10^6} = 3.7$$

Exp. 5Mutant Colonies

- 1) 3, 2, 2, 1, 3, 1, 1, 2, 1, 4
- 2) 3, 2, 5, 4, 8, 8, 4, 5, 9
- 4) 6, 12, 9, 7, 8, 10, 6, 7, 9, 10
- 5) 7, 6, 9, 6, 5, 6, 9, 6, 6, 4
- 6) 8, 8, 3, 4, 6, 2, 5, 5, 4, 4
- 7) 7, 7, 10, 9, 6, 14, 5, 8, 9, 7
- 8) 15, 4, 8, 5, 8, 9, 7, 9, 10, 9
- 10) 9, 3, 8, 12, 8, 8, 9, 10, 6, 6

Plating Efficiency

- 1) 191, 193, 199
- 2) 272, 208, ~~280~~ 241
- 4) 271, 281, 290
- 5) 234, 223, 275
- 6) 221, 238, 225
- 7) 260, 239, 230
- 8) 208, 225, 211
- 10) 245, 248, 227

Exp 5 Selection
Stain 8/13
8/12

Exp 5 Plating Eff. 8/12
Stain 8/9

Parameters

Date	7/16/01
Experiment No.	5
Investigator	R. Howell
Cell Line	V79 Frozen 10/00
Modifier	None
Radionuclide	H-3
Half-life (days)	4500.45
Radiation Yield	1
Radiochemical	H3TdR
Manufacturer/Lot	NEN NET-027Z Lot 3106-431
Original Calibration Date/Time	7/5/01 0:00
Present Calibration Date/Time	7/16/01 0:00
Fraction of Cells Labeled	1
	Original Activity Concentration (MBq/ml) 37
	Time Elapsed Since Original Calibration (d) 11
	Present Activity Concentration (MBq/ml) 36.94
Liquid Scintillation Cocktail	Ecolume
Volume of LSC Cocktail (ml)	6
Volume/Type Counting Vial	7 ml
Model of Counter	Beckman 5500
Counting Efficiency	0.65
Activity Added (Date/Time)	7/16/01 18:30
Cells Washed (Date/Time)	7/17/01 10:30
Medium Tubes Counted (Date/Time)	7/23/01 13:22
Cell Tubes Counted (Date/Time)	7/23/01 13:22
Vol. Supernatant Counted (μl)	30
Vol. Suspension Counted Cell Activity (μl)	200
	Time Elapsed Between Add and Wash (hr) 14.00
	Time Elapsed Between Add and Count (hr) 163.00
	Time Elapsed Between Wash and Count (hr) 163.00
Vol. Suspension Coultter (μl)	100
Coultter Manometer Volume (μl)	500
Average Coultter Background Counts	22
Coultter Calibration Parameter	400
Hemocytometer Counting (Yes or No)?	No
	Background
	Coultter 1 22
	Coultter 2 22
	Coultter 3 22

MediumActivity

Experiment: 5
Date: 7/16/01

Tube #	1st	2nd	3rd	CPM Average	CPM corrected for control	DPM CPM(y e)	At $\mu\text{Ci/ml on counting}$	Ao $\mu\text{Ci/ml at addition}$	Ao $\text{kBq/ml at addition}$
1	9	8	7	9	0	0	0	0	0
2	5	16	8	8	0	0	0	0	0
3	28785	26352	28608	27915	27906	42933	0.6446	0.6453	23.8764
4	57790	55970	57645	57135	57126	87886	1.3196	1.3210	48.8769
5	85580	83353	82146	83693	83684	128745	1.9331	1.9351	71.5998
6	115919	116026	117066	116337	116328	178966	2.6872	2.6900	99.5298
7	142806	134140	142500	139815	139807	215087	3.2295	3.2329	119.6178
8	346220	330266	336800	337762	337753	519620	7.8021	7.8103	288.9800
9	427873	409359	416326	417853	417844	642837	9.6522	9.6623	357.5052
10	525273	548100	548013	540462	540453	831466	12.4845	12.4975	462.4091

CellSuspension

Experiment: 5
Date: 07/16/01

Tube #	Suspension count (CPM)			CPM Average	CPM corrected for control	DPM CPM/(y e)	A _i μCi/ml on counting	A _o μCi/ml after uptake	A _o kBq/ml after uptake
	1st	2nd	3rd						
1	7	9	17	9	0	0	0.0000	0	0.0000
2	7	6	7	0	0	0	0.0000	0	0.0000
3	2515	2487		2501	2492	3834	0.00864	0.00864	0.3198
4	5132	4865	4606	4868	4859	7475	0.01684	0.01685	0.6236
5	9053	8805	8504	8787	8779	13505	0.03042	0.03045	1.1266
6	11281	11345	10826	11151	11142	17141	0.03861	0.03865	1.4299
7	13424	13333		13379	13370	20569	0.04633	0.04637	1.7159
8	41116	39310	36616	39014	39005	60008	0.13515	0.13529	5.0059
9	27582	20047	21732	23120	23112	35556	0.08008	0.08017	2.9661
10	50000	53260		51630	51621	79417	0.17887	0.17905	6.6250

CoulterSurvival

Experiment: 5
 Date/Time: 7/16/01



Tube #	Coulter count			Average Cells/ml		Hemocytometer Count in Grid			
	1st	2nd	3rd			1st	2nd	3rd	4th
1	4827	4878	4587	4764	1896800				
2	5348	5384	5312	5348	2130400				
3	4752	4898	4672	4774	1900800				
4	5120	4944	4934	4999	1990933				
5	5457	5625	5562	5548	2210400				
6	5265	5084		5175	2061000				
7	4929	4878	4753	4853	1932533				
8	5115	4965	5094	5058	2014400				
9	4615	4711	4583	4636	1845733				
10	4662	4596	4647	4635	1845200				

Tube #	Predicted # Cells Seeded	Actual # Cells Seeded	Colony count			Average	PE (%)	SF Uncorrected	SF Corrected
			1st	2nd	3rd				
1	200	190	108	114	114	100	49.861	1.00	1.0000
2	200	213	85	81					
3	200	190	63	68	83	71	37.528	0.7105	0.7527
4	200	199	94	84	73	84	42.024	0.8333	0.8428
5	200	221	66	108	67	80	36.343	0.8001	0.7289
6	200	206	61	93	73	76	36.714	0.7537	0.7363
7	200	193	56	55	52	54	28.115	0.5412	0.5639
8	200	201	47	38	37	41	20.188	0.4050	0.4049
9	200	185	46	38	54	46	24.922	0.4582	0.4998
10	200	185	40	30	30	33	18.065	0.3320	0.3623

Summary



Experiment: 7/16/01
 Date/Time:

Tube #	Activity Conc. (kBq/ml)	Activity/Cell (mBq/cell)	Survival Uncorrected	Survival Corrected
1	0.000	0.000	1.0000	1.0000
2	0.000	0.000	0.7105	0.7527
3	23.876	0.168	0.8333	0.8428
4	48.877	0.313	0.8001	0.7289
5	71.600	0.510	0.7537	0.7363
6	99.530	0.694	0.5412	0.5639
7	119.618	0.888	0.4050	0.4049
8	288.980	2.485	0.4582	0.4998
9	357.505	1.607	0.3320	0.3623
10	462.409	3.590		



Experiment:
Date/Time: 7/16/01

Tube #	Activity Conc. (kBq/ml)	Activity/Cell (mBq/cell)	Survival Uncorrected	Survival Corrected
1	0.000	0.000	1.0000	1.0000
2	0.000	0.000	0.7105	0.7527
3	23.876	0.168	0.8333	0.8428
4	48.877	0.313	0.8001	0.7289
5	71.600	0.510	0.7537	0.7363
6	99.530	0.694	0.5412	0.5639
7	119.618	0.888	0.4050	0.4049
8	288.980	2.485	0.4582	0.4998
9	357.505	1.607	0.3320	0.3623
10	462.409	3.590		



HPRT mutant selection

Date: 20-Aug-01

Experiment: V79, HTdR, 100 % cluster, Roger's exp. #5

Sample #	Coulter H	count		Coulter bckgr.	Coulter Mode (μ)	# of cells per ml	Cell susp. volume for 200 000 cells (ml)
1	2503	2631	2542	5	100	1021467	0.196
2	2411	2478	2399	5	100	969733	0.206
3		N O S A M P L I		5			
4	2616	2806	2650	5	100		
5	2692	2730	2659	5	100	1075467	0.186
6	2224	2289	2213	5	100	894800	0.224
7	2206	2335	2294	5	100	909333	0.220
8	2497	2317	2410	5	100	961200	0.208
9				5	100	#DIV/0!	#DIV/0!
10	2087	2138	2051	5	100	834800	0.240

Protocol:

1. Wash T75 1-2x with PBS (no Ca++, no Mg++)
2. Trypsinize the cells (2 ml trypsin /T75, 2-3 min, RT)
3. Add 10 ml wash medium (wash MEMA) / T75 flask
4. Resuspend cells and transfere them in 15 ml conical tube.
5. Spin the cells down: 4-5 min, 2K rpm)
6. Aspirate supernatant, click tube to disperse pelet
7. Add 5 ml regular medium (MEMB/FCS10%) and resuspend the cells
8. Siringe the cells using 5 ml singe & 21G needle.
9. Count the cells using coulter couter: 100 μ l cell susp. + 20 ml Isotone.
10. Plate 2 x 10(5) cells / P100 dish in 7 ml MEMA/FCS10% x 10 dishes/dose point
11. Transfere dishes into incubator and let the cells to attach (2-4 hrs)
12. Take 2 x 10(5) cells suspension, make serial dilution (2 x 10x) to obtain 200 cells/ml.
13. Plate 200 cells / P 60 x 3 dishes / dose point for Plating Efficiency.
14. Repaet staep 11 and 12 for each sampling point.
15. After 2-4 hrs add 3 ml MEMA + 6-TO (60 μ g 6-TO) into each P 100.
(Final concentartion for 6-TO = 6 μ g/ml) 6-TO Stock sol. = 60 μ g/ 3 ml MEMA)
16. Keep the dishes in standard culture condition for 8-10 days.
17. Wash HPRT- colonies 1-2 times with PBS, fix them with MetOH and stain.
18. Count colonies.

Table 1. Plating efficiency for HPRT challenge (Roger's exp. #5)

Set	HTdR (mBq/cell)	Number of cells originally plated	S u r v i v a l		Abs PE	Abs PE Avg.	+/- Std
			Numer of colonies	Avg. # of col./ plate			
1	0	200	191	194	0.96	0.97	0.02
		200	193		0.97		
		200	199		1.00		
2	0	200	272	240	1.36	1.20	0.16
		200	208		1.04		
		200	241		1.21		
3	0.1685	200		#DIV/0!	0.00	0.00	0.00
		200			0.00		
		200			0.00		
4	0.313	200	271	281	1.36	1.40	0.05
		200	281		1.41		
		200	290		1.45		
5	0.5095	200	234	244	1.17	1.22	0.14
		200	223		1.12		
		200	275		1.38		
6	0.694	200	221	228	1.11	1.14	0.04
		200	238		1.19		
		200	225		1.13		
7	0.888	200	260	243	1.30	1.22	0.08
		200	239		1.20		
		200	230		1.15		
8	2.485	200	208	215	1.04	1.07	0.05
		200	225		1.13		
		200	211		1.06		
9	1.607	200		#DIV/0!	0.00	0.00	0.00
		200			0.00		
		200					
10	3.5905	200	245	240	1.23	1.20	0.06
		200	248		1.24		
		200	227		1.14		

Table 2. HPRT mutant frequency for V79 cells exposed to HTdR in 50% cluster model (Roger's exp #5)

Sample #	Dose (mBq/cell) ($\times 10^5$)	HPRT mutants per plate										SUM	No of dishes plated	Total # of cells plated	PE for challenge (in Tab.1)	HPRT per plated cell	mutants per 100000 survivors	#DIV/0!		
		1st	2nd	3rd	4th	5th	6th	7th	8th	9th	10th									
1	0.00	2	3	2	2	1	3	1	1	2	1	4	20	10	2000000	0.97	0.000010	0.000010	1.0	
2	0.00	2	3	2	5	4	8	4	5	9		48	9	1800000	1.20	0.000027	0.000022	2.2		
1+2	0.00	2										68	19	3800000	1.09	0.000018	0.000016	1.6	0.0	
3															0	0.00	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
4	0.31	2	6	12	9	7	8	10	6	7	9	10	84	10	2000000	1.40	0.000042	0.000030	3.0	1.3
5	0.51	2	7	6	9	6	5	6	9	6	6	4	64	10	2000000	1.22	0.000032	0.000026	2.6	1.0
6	1.39	2	8	8	3	4	6	2	5	5	4	4	49	10	2000000	1.14	0.000025	0.000021	2.1	0.5
7	0.89	2	7	7	10	9	6	14	5	8	9	7	82	10	2000000	1.22	0.000041	0.000034	3.4	1.7
8	2.49	2	15	4	8	5	8	9	7	9	10	9	84	10	2000000	1.07	0.000042	0.000039	3.9	2.3
9															0	0.00	#DIV/0!	#DIV/0!	#DIV/0!	#DIV/0!
10	3.59	2	9	3	8	12	8	8	9	10	6	6	79	10	2000000	1.20	0.000040	0.000033	3.3	1.6

Statistic

t-test
(p)

#DIV/0!	p<0.05
#DIV/0!	p<0.05
#DIV/0!	p<0.05
#DIV/0!	p<0.05
#DIV/0!	p<0.01
#DIV/0!	p<0.01
#DIV/0!	p<0.01

