

## HPRT mutant selection

Date: 27-Sep-01

Experiment: V79, HTdR, 100 % cluster, Roger's exp. #6

Sample #	Coulter H	count			Coulter bckgr.	Coulter Mode ( $\mu$ l)	# of cells per ml	Cell susp. volume for 200 000 cells (ml)
1	2545	2567	2587	5	100	5122667	0.039	
2	2452	2471	2458	5	100	4910667	0.041	
3	2429	2309	2446	5	100	4779333	0.042	
4	2352	2353	2288	5	100	4652000	0.043	
5	2350	2377	2367	5	100	4719333	0.042	
6	2184	2096	2178	5	100	4295333	0.047	
7	2133	2242	2163	5	100	4348667	0.046	
8	2491	2436	2396	5	100	4872000	0.041	
9	2286	2298	2245	5	100	4542667	0.044	
10	2289	2281	2433	5	100	4658667	0.043	

## Protocol:

1. Wash T75 1-2x with PBS (no Ca<sup>++</sup>, no Mg<sup>++</sup>)
2. Trypsinize the cells (2 ml trypsin /T75, 2-3 min, RT)
3. Add 10 ml wash medium (wash MEMA) / T75 flask
4. Resuspend cells and transfere them in 15 ml conical tube.
5. Spin the cells down: 4-5 min, 2K rpm)
6. Aspirate supernatant, click tube to disperse pelet
7. Add 5 ml regular medium (MEMB/FCS10%) and resuspend the cells
8. Siringe the cells using 5 ml singe & 21G needle.
9. Count the cells using coulter couter: 100  $\mu$ l cell susp. + 20 ml Isotone.
10. Plate 2 x 10(5) cells / P100 dish in 7 ml MEMA/FCS10% x 10 dishes/dose point
11. Transfere dishes into incubator and let the cells to attach (2-4 hrs)
12. Take 2 x 10(5) cells suspension, make serial dilution (2 x 10x) to obtain 200 cells/ml.
13. Plate 200 cells / P 60 x 3 dishes / dose point for Plating Efficiency.
14. Repaet staep 11 and 12 for each sampling point.
15. After 2-4 hrs add 3 ml MEMA + 6-TO (60  $\mu$ g 6-TO) into each P 100.  
(Final concentartion for 6-TO = 6  $\mu$ g/ml) 6-TO Stock sol. = 60  $\mu$ g/ 3 ml MEMA)
16. Keep the dishes in standard culture condition for 8-10 days.
17. Wash HPRT- colonies 1-2 times with PBS, fix them with MetOH and stain.
18. Count colonies.