

13. The tubes were irradiated using Mark I irradiator (^{137}Cs gamma-ray), one tube at a time, while placing onto a separate Rainin pipet tip box containing ice as per the Table below

Tube #	Total Dose (Rad)	Dose rate (Rad/min)	Time (min)	Attenuat.
1	0	0	0	0
2	0	0	0	0
3	200 100	101.4	0.98 0.98	X-10
4	400 200	101.4	1.97 1.97	X-10
5	600 300	101.4	2.95 2.95	X-10
6	800 400	101.4	3.94 3.94	X-10
7	1000 500	101.4 101.4	4.93 4.93	X-10
8	1500 750	169.6	7.42 7.42	X-5
9	2000 1000	169.6	5.89	X-5
10	2500 1250	169.6	7.37	X-5

14. After irradiation, carefully remove the supernatant from the top, resuspend pellet in 200 μl wash MEMA and transfer the content to ten 12 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using pasteur pipet
15. Again add 200 μl wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions in 12 ml tubes
16. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)
17. Labeling and preparation of dilution tubes and colony dishes
- load 57 60 mm petri dishes with 4 ml MEMA
 - load 39 T-tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.
18. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
19. Centrifuge tubes for 10 min at 2000 rpm, 4°C
20. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 3 cc syringe with 21 gauge needle
21. Determine cell concentration by transferring 100 μl to Coulter cup
22. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5 and transfer 0.5 ml to tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3 and vortex tube X.3 and transfer 0.5

ml to tube X.2. Keep tubes on ice.

23. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.
24. Incubate petridishes for 1 week
25. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol.
Stain colonies with 0.05% crystal violet
26. Count colonies. There must be between 25 and 250 colonies for the flask to be a valid data point.

Expt # 1

05/07/98

Monomer Select = 50 ml.

80-90% confluent: suspend in 5 ml
MEMB for each flask

Initial cell count = 1296, 1347, 1372

Avg cell count = 1338.3

Cell conc. = 1338.3×4000
= 5353333 Cells/ml

For dilution

$$\text{Vol. of cell suspension required} = \frac{44000000}{5353333} = 8.2 \text{ ml.}$$

Take 8.2 ml of cells + 2.8 ml MEMB = 11 ml

After dilution,

Final Cell count = 1047, 1012, 977

Avg. Cell count = 1012

Cell conc. = 1012×4000

= 40,48,000 Cells/ml

1.2

2.2

3.2

4.2

5.2, 5.3

6.2, 6.3

7.2, 7.3, 7.4

8.2, 8.3, 8.4

9.2, 9.3, 9.4

~~10.2~~, 10.3, 10.4

05/11/98

Expt #1

Tube #	Coulter count for 100 μ l cell suspension (MS = 50 μ l)
1	511, 506, 469
2	460, 417, 418
3	432, 452, 463
4	501, 487, 482
5	428, 420, 444
6	422, 407, 427
7	398, 402, 419
8	414, 408, 441
9	403, 405, 413
10	455, 429, 427

05/18/98

Expt #1

Tube #

of colonies

Avg. colonies
for X.2

SF

1.2	65, 60, 62	}	64.16	
2.2	59, 67, 72		61.00	
3.2	64, 58, 61		50.00	0.8950
4.2	52, 47, 62		53.66	0.83
5.2	42, 40, 38		40	0.62
6.2 ✓	41, 31, 32		33.77	0.52
7.2 ✓	28, 30, 27		28.33	0.44
8.2 ✓	25, 29, 24		26	0.40
9.3	88, 84, 84		85	0.13
10.3	65, 62, 60		62	0.097


4x10⁶ cells cluster

D₃₇ = 560 r

47

46 6010
SF

K+E SEMI-LOGARITHMIC 4 CYCLES X 70 DIVISIONS
KEUFFEL & ESSER CO. MADE IN U.S.A.



Chronic