

### V79 COLONY FORMING ASSAY

Experiment Name : <sup>137</sup>Cs toxicity (10<sup>6</sup> Cells cluster, Chronic, 10.5° C); Exp. # : 1;  
 Investigator: A. Bishayee Date: 10/13/98

1. Set the Coulter Counter, wash cells (from two 150 cm<sup>2</sup> flusk, subcultured 1:2, 24h before) with PBS, trypsinize cells, resuspend in 7 ml MEMB for each flusk, pool, vortex, pass five times through 3 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml isotone (Coulter balanced electrolyte solution)
2. Dilute to ~4,000,000 cells/ml in MEMB (final volume 11 ml) [Actual count : 4,770,000 cells/ml]
3. Transfer 1 ml of cell suspension into ten 6 ml tubes (Falcon plastic test tube, 12x175 mm) labeled 1-10 both on cap and wall
4. Add 2 ml MEMA, vortex and centrifuge at 2000 rpm at 4°C for 10 min (precooled centrifuge).
5. Decant the supernatant, click tubes and vortex
6. Transfer the cell suspension in prelabeled (1-5 and C) polypropylene microcentrifuge tubes (Helena Plastics, 400 ul) with attached caps using 200 ul pipet tips
7. Again add 200 ul ice cold MEMA, vortex to resuspend and transfer the cell suspensions in the same polypropylene microcentrifuge tubes. (Total volume in each tube ~400 ul)
8. Centrifuge tubes for 5 min at 1000 rpm, 4°C
9. Place tubes at different distances (positions #1-5) from <sup>137</sup>Cs source (placed inside a 6 ml plastic centrifuge tube in a styrofoam platform) at 10°C for 72 h. Place 2 control tubes away from the source inside the 10.5° C incubator. Date/Time: 10/13/98; 7-00 p.m.
10. After 72h of chronic irradiation, remove the tubes from incubator, carefully remove the supernatant from the top, resuspend pellet in 200 ul wash MEMA and transfer the content (using pasteur pipet) to 12 ml tubes (Falcon plastic test tube, 17x100 mm, labeled both on cap and wall) containing 10 ml wash MEMA Date/Time: 10/16/98; 5-00 p.m.
11. Again add 200 ul wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions in 12 ml tubes
12. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)
13. Labeling and preparation of dilution tubes and colony dishes
  - load 60 mm petri dishes with 4 ml MEMA
  - load test tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.
14. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
15. Centrifuge tubes for 10 min at 2000 rpm, 4°C

16. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 3 cc syringe with 21 gauge needle
17. Determine cell concentration by transferring 100  $\mu$ l to Coulter cup
18. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5 and transfer 0.5 ml to tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3 and vortex tube X.3 and transfer 0.5 ml to tube X.2. Keep tubes on ice.
19. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.
20. Incubate petridishes for 1 week
21. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol. Stain colonies with 0.05% crystal violet
22. Count colonies. There must be between 25 and 250 colonies for the <sup>dish</sup>flask to be a valid data point.

CS-137 in clusters  
(4 x 10<sup>6</sup> cells)

Exp #1

10/13/98

50

Take 2x 150 cm<sup>2</sup> flask (80-90% confluent)

Total vol. = 16 ml

Initial cell count = 1407, 1397, 1422

Avg. cell count = 1408.6

Cell conc. = 1408.6 x 4000 = 5,634,666 Cells/ml

we need <sup>8</sup> ml of 4,000,000 Cells/ml = ~~32,000,000 Cells~~  
~~32 ml~~ = 32,000,000

Vol. required =  $\frac{32 \times 4,000,000}{5,634,666}$

= 5.67 ml.

Take 6 ml cells + 2 ml MEMB = 8 ml

Final cell count = 11921, 11909, 11945

Avg. count = 11925

Cell conc. = 11925 x 400

= 4,770,000 cells/ml

= 4.7 million/ml.

6/30 p.m.

10/13/98

SETUP

UNIT:

LOCATION:

ENERGY:

SSD:

GANTRY ANGLE:

COLLIMATION:

FIELD SIZE:

EXPOSURE:

Calibration of Cs-137  
Cluster Irradiation

10mCi Cs-137  
1/76

On 4:30 10/19

Read 8:52 10/20

16.37h

Calibration of Low  
Dose Rate Cs-137 Irradiation  
of V79 Clusters irradiated  
at 10.5°C

Source NEN 10mCi Cs-137  
pellet calibrated on 1/76.

\* 16:30 OCT 19/98  
ZERO SENSORS

Source in

A1 TOTAL: 1.024 mV	HIGH RESOLUTION
A2 TOTAL: 1.137 mV	HIGH RESOLUTION
A3 TOTAL: 986 mV	HIGH RESOLUTION
A4 TOTAL: 881 mV	HIGH RESOLUTION
A5 TOTAL: 1.173 mV	HIGH RESOLUTION

A1 08:52 OCT 20/98

DOSE: 298 RAD  
TOTAL: 1.322 mV

18.2  $\frac{\text{rad}}{\text{h}}$  position #1

A2 08:52 OCT 20/98

DOSE: 161.3 RAD  
TOTAL: 1.298 mV

9.86  $\frac{\text{rad}}{\text{h}}$  #2

A3 08:52 OCT 20/98

DOSE: 107.6 RAD  
TOTAL: 1.094 mV

6.57  $\frac{\text{rad}}{\text{h}}$  #3

A4 08:52 OCT 20/98

DOSE: 99.9 RAD  
TOTAL: 981 mV

6.10  $\frac{\text{rad}}{\text{h}}$  #4

A5 08:52 OCT 20/98

DOSE: 55.6 RAD  
TOTAL: 1.229 mV

3.40  $\frac{\text{rad}}{\text{h}}$  #5

THOMSON & NIELSEN  
MULTI DOSIMETER  
SYSTEM

REV 2.0

CAL FACTOR SETTINGS

A1: 1.00mV/RAD	C1: 1.00mV/RAD
A2: 1.00mV/RAD	C2: 1.00mV/RAD
A3: 1.00mV/RAD	C3: 1.00mV/RAD
A4: 1.00mV/RAD	C4: 1.00mV/RAD
A5: 1.00mV/RAD	C5: 1.00mV/RAD
B1: 3.15mV/RAD	D1: 1.00mV/RAD
B2: 3.15mV/RAD	D2: 1.00mV/RAD
B3: 3.15mV/RAD	D3: 1.00mV/RAD
B4: 3.15mV/RAD	D4: 1.00mV/RAD
B5: 3.15mV/RAD	D5: 1.00mV/RAD

NOTES

Calibration of  $^{137}\text{Cs}$  irradiator in  $10^5^\circ\text{C}$  incubator

12:00 OCT 13/98  
ZERO SENSORS

4.75 h irradiation

A1 TOTAL: 983 mV HIGH RESOLUTION  
A2 TOTAL: 1.122 mV HIGH RESOLUTION  
A3 TOTAL: 963 mV HIGH RESOLUTION  
A4 TOTAL: 865 mV HIGH RESOLUTION  
A5 TOTAL: 1.179 mV HIGH RESOLUTION

1

A1 16:46 OCT 13/98

DOSE: 37.7 RAD  
TOTAL: 1.021 mV

7.9 rad/h

2

570

A2 16:46 OCT 13/98

DOSE: 23.8 RAD  
TOTAL: 1.146 mV

5.0 rad/h

3

360

A3 16:46 OCT 13/98

DOSE: 25.1 RAD  
TOTAL: 988 mV

5.3 rad/h

4

A4 16:46 OCT 13/98

DOSE: 20.4 RAD  
TOTAL: 885 mV

4.3 rad/h

310

A5 16:46 OCT 13/98

DOSE: 8.5 RAD  
TOTAL: 1.187 mV

1.8 rad/h

5

130

100 $\mu$ l Cell Count

NOTES

MS = 50 $\mu$ l

C	485, 502, 498
1	447, 462, 444
2	411, 399, 436
3	428, 409, 425
4	389, 389, 374
5	495, 472, 458

NOTES



TABLE-2

Expt # : 1

Date : 10/23/98 ; 10-00 A.M.

Tube.dilution	Colony 1	Colony 2	Colony 3	Avg Colony	SF
1-2	73	83	75	77	—
2-2	24	26	29	26.33	0.342
3-2	39	37	43	39.66	0.515
4-2	58	49	67	58	0.7532
5-2	61	60	68	63	0.818
6-2	68	69	78	71.66	0.931

Chronic CS-137 toxicity :  $4 \times 10^6$  Cells cluster, 10.5°C for ~70h

10/23/98


57

### Summary of Results

Tube # (Position #)	Total Dose (Rad)	Dose rate (Rad/h)	Time (h)	SF
C	0	0	0	-
C	0	0	0	-
1	1274	18.2	70	0.342
2	690	9.8	70	0.515
3	460	6.5	70	0.753
4	427	6.1	70	0.818
5	238	3.4	70	0.931
6				
7				
8				

Chronic CS-137 :  $4 \times 10^6$  cells cluster at  $10^5$ °C for ~70h ;  $D_{37} = 1360.2$

SF



# Check out these WilsonJones differences!

**Cure GAPiTIS!** Conventional rings can become spread or misaligned—we call it GAPiTIS! Cure it with DublLock! Rings stay closed and aligned—sheets won't tear or hang-up on rings.

**Locking Rings!** WilsonJones patented ring mechanism actually locks the outside rings. Our 3 position trigger lets you lock, unlock, and open rings—sheets won't fall out!

**EasyFlow Sheet Lifters!** We've got a new angle on sheet lifters! The patented angle shape of WilsonJones EasyFlow lifter won't catch or tear papers stored in pockets, won't take up valuable ring capacity like conventional sheet lifters! Available in 1½, 2 & 3 inch capacities.

**Extra-Wide Covers!** WilsonJones DublLock binders feature extra-wide covers that keep sheet protectors and index tabs inside your binder—where they belong!


**Identify and Organize!** WilsonJones unique color coding labeling system lets you file instantly! Just write or type in titles—labels are pica spaced—and insert!

**Available in 12 great colors!** All WilsonJones traditional and contemporary colors are heavy duty suede grain virgin vinyl—for longer wear with no scuffing or fingerprints!


**The  
WilsonJones®  
Difference**




Cure GAPiTIS! DublLock



Sheets won't fall out



EasyFlow Sheet Lifters



Identify and organize

*Color coding made easy! Type or write in titles.  
Tear along perforations and insert into label holders.*

354 Line