

Investigator: A. Bishaye

**Granulocyte Macrophage Colony Forming Unit (GM-CFU) Assay**

Experiment # 3 (Survival #1) Source of Irradiation: Radionuclide injection  
 Mice Sex, Strain, Age: SW, F, 5-6wk  
 Type of Irradiation: Chronic  
 Animals per group: 3

**Aims:**

- i) To determine the bone marrow response with dose after Sn-117m administration.

**Summary of Results:**

**Brief Procedure:**

- 1) Inject animals in groups of 3 with desired activity of Sn-117m intravenously through lateral tail vein.
- 2) Sacrifice each mouse on optimal day by cervical dislocation and sterilize using 70% EtOH and immediately move it into laminar flow hood.
- 3) Remove both femurs carefully using sterile instruments and clean the attached tissue thoroughly.
- 4) Flush the bone marrow with 2% Horse Serum in Dulbecco's Modified Eagles Medium (2% HS-DMEM) using 21G needle and syringe.
- 5) Separate the mononuclear cells by density gradient procedures using Histopaque.
- 6) Plate the desired number of cells (cell suspension) in mixture of 60% HS-DMEM and 1.7 ml 0.6% bacto agar solution in the presence of 9.2 U (New Sigma Unit) of GM-CFS.
- 7) Keep the plated petri dishes for 20 min. in laminar flow hood and move them into incubator with 5% CO<sub>2</sub> and 95% air, at 37°C.
- 8) Count the granulocyte macrophage colonies on 7th day.

**Sn-117m injections:**

Group # Probe #	Date of injection	Activity injected ( <del>Bq</del> )	Date Sacrificed	# of days	Remarks if any
C1	8/28/98	0	9/4/98	7	0.2 ml 0.9% Nacl injected, control
<del>2</del>					
1	8/28/98	5.95 <sup>58</sup>	9/4/98	7	
2	"	8.0 <sup>47</sup>	"	"	
3	"	223 <sup>52</sup>	"	"	
4	"	31.5 <sup>468</sup>	"	"	
5	"				
6					
7					
8					

**Preparing Media and Agar:**

**Culture Medium (Double Strength, 2X):** 13.37g (1 pack) of DMEM powder (Gibco, Cat # 12100-038) + 490 ml deionized water + 16  $\mu$ l of L-asparagine (Gibco Cat # 12416-012) at a concentration of 5  $\mu$ g/ $\mu$ l + 150.4  $\mu$ l of DEAE dextran (mol. wt. =  $2 \times 10^6$ , intrinsic viscosity = 0.7) at a concentration of 1  $\mu$ g/ $\mu$ l (Sigma Cat # D-9885) + 10 ml of penstrep (Gibco Cat # 600-5070, 5,000 units/ml pen, 5,000  $\mu$ g/ml streptomycin) + 3.7 g of NaHCO<sub>3</sub> (Gibco Cat # 11810-025).

**Culture Medium (2X) with 60% Horse Serum :** Add 60 % Horse Serum to 2X DMEM

**Wash Medium:** i) Mix equal amounts of culture medium and sterile deionized water.  
ii) Add 2% HS

**Agar:** Prepare 0.6% agar by adding 0.6 g Difco Bacto agar (Difco Cat # 0140-15-4) to 100 ml boiling deionized water. Autoclave on liquid cycle for 20 min.

**Comments If any:**

**Flushing Bone marrow:**

- 1) Remove both femurs from each mouse and place them in a test tube containing wash medium kept in ice, if the femur can not be flushed immediately.
- 2) Flush the marrow from each femur by aspirating 3 ml of Wash Medium through the femur 5 times with a 21G needle/3 ml syringe in a 50 ml conical centrifuge tube. Follow with two flushes with 1 ml of fresh Wash Medium.
- 3) Spin the cells at 1200 rpm for 5 minutes at 4°C, decant, break up the cell pellet, resuspend the cells in 5 ml of cold Wash Medium, and vortex the cell suspension.

**Comments If any:**

One of the bones of Group 4 was broken, so not used in this study.

**Counting the Cells:**

Add 10 $\mu$ l of cell suspension to 20 ml of Isotone II in a coulter cup and count the cells using coulter counter. Calculate total # of cells in each group.

**Coulter Counter Parameters:**Current(I)=500  $\mu$ A

Attenuation= 4

Full Scale = 1

Alarm Threshold = off

 $T_L$  = 2.7

Preset Gain = 1

 $T_U$  = 99.9

Stirrer control = off

Manometer Select = 500  $\mu$ l

*Multiplication Factor to get total # of cells in 5 ml = 20,000 x Coulter count*

Group #	Coulter Count without ZG	Avg	Total # of cells	Coulter Count with 5 drops ZG	Avg	Total # of cells
C1	NOT PERFORMED					
C2						
1						
2						
4						
5						
6						
7						
8						

Comments If any:

**Separating Mononuclear cells and washing the cells:**

- 1) Transfer 3.5 ml of Histopaque (Sigma Cat #H8889) into fresh 15 ml tubes (1 tube per group).
- 2) Layer the cell suspension carefully on top of Histopaque and centrifuge at 1500 rpm, 4°C, for 30 minutes.
- 3) Using a Pasteur pipette transfer the mononuclear cells into fresh tubes.
- 4) Dilute the cell suspension to 15 ml by adding cold Wash Medium into each tube and spin them at 1200 rpm, 4°C, for 5 min.
- 5) Decant the supernatant, break the pellet, and add 15 ml cold Wash Medium, and spin them again at 1200 rpm, 4°C, for 5 min. Repeat this procedure 2 more times.
- 6) After 3rd wash break the pellet and resuspend in 2 ml Culture medium (2x DMEM) with 60% HS and keep the tubes in dry bath at 37°C.
- 7) Add 20  $\mu$ l of cell suspension to 20 ml of Isotone II in a coulter cup and determine total # of cells in each group using coulter counter

Coulter Counter Parameters: Same as above

*Multiplication Factor to get # of cells/ml = 2000 x Coulter Count*

Group #	Coulter Count without ZG	Avg	# cells per ml	Coulter Count with 5 drops ZG	Avg	# cells per ml
C1				5110, 5125, 5213	5149	10,298,666
C2				<del>5007, 5107, 5123</del>		
1				5007, 5107, 5123		
2				4813, 4918, 4975		
3				4439, 4556, 4579		
4				4322, 4429, 4372		
5						
6						
7						
8						

Average # of cells per  $\mu\text{l} = \text{ml}$  for control = 10,298,666

<b>DILUTIONS</b>
------------------

**GROUP 1****Dilution A:** (1.0x10<sup>6</sup> cells /ml, Total volume 3.4 ml)

$$\text{Vol of final cell suspension required} = 3400000 / \text{cells per ml} = \frac{3400000}{10298666} = 0.330 \text{ ml}$$

1.7 ml Agar + 1.370 ml Medium + 0.330 ml Cell Suspension

**Dilution B:** (3.0x10<sup>5</sup> cells /ml, Total volume 3.4 ml)

$$\text{Vol of final cell suspension required} = \frac{1020000}{340000} / \text{cells per ml} = \frac{1020000}{10298666} = 0.099 \text{ ml}$$

1.7 ml Agar + 1.60 ml Medium + 0.099 ml Cell Suspension

**GROUP 2****Dilution A:** (1.0x10<sup>6</sup> cells /ml, Total volume 3.4 ml)

$$\text{Vol of final cell suspension required} = 3400000 / \text{cells per ml} =$$

1.7 ml Agar +        ml Medium +        ml Cell Suspension

**Dilution B:** (3.0x10<sup>5</sup> cells /ml, Total volume 3.4 ml)

$$\text{Vol of final cell suspension required} = 3400000 / \text{cells per ml} =$$

1.7 ml Agar +        ml Medium +        ml Cell Suspension

**GROUP 3****Dilution A:** (1.0x10<sup>6</sup> cells /ml, Total volume 3.4 ml)

$$\text{Vol of final cell suspension required} = 3400000 / \text{cells per ml} =$$

1.7 ml Agar +        ml Medium +        ml Cell Suspension

**Dilution B:** (3.0x10<sup>5</sup> cells /ml, Total volume 3.4 ml)

$$\text{Vol of final cell suspension required} = 3400000 / \text{cells per ml} =$$

1.7 ml Agar +        ml Medium +        ml Cell Suspension

Plating the Cells:

Culture Medium: Maintain four 13mm tubes each containing 4.5 ml of Culture medium in dry bath at 37°C.

Horse Serum: Maintain five 13mm tubes each containing 4.5 ml of Horse Serum in dry bath at 37°C.

Agar: Maintain five 16mm tubes each containing 6.5 ml of Agar in dry bath at 37°C.

- 1) Warm up dilution tubes (2 or 3 per group) to 37°C in dry bath.
- 2) Warm up Agar (30 ml) and 60% HS in 2x DMEM (30 ml) in separate tubes to 37°C.
- 3) Mark the Six-well plates (3 wells for each dilution for each group) containing 20 ul of stock GM-CSF (9.2 U) in each well.
- 4) Mix 1.7ml agar + x ml of 2x DMEM with 60% HS + y ml cell suspension + 0.02 ml GMCSF (x + y = 1.7 ml) in a dilution tube.
- 5) Add 1 ml of mixture 4 to each well, mix properly and let it gel for about 30 minutes.
- 6) Repeat steps 4 and 5 for each dilution.
- 7) Repeat steps 1 to 6 for each group.
- 8) Incubate the cells in an incubator at 37°C and 5% CO<sub>2</sub>, 95% air for 7 days.
- 9) On 8th day of incubation count colonies and determine the survival fractions.

Counting the Colonies: (Inverted at 40X or dissecting at 35X)

Group #	Dose / Activity # injected	# of cells plated	# CFU-GM counted	Avg	SF
C1	0	$3 \times 10^5$	144, 131, 137	137.33	
C2					
1	3.85	$3 \times 10^5$	122, 131, 136	129.66	0.9441
2	8	$3 \times 10^5$	123, 115, 109	115.66	0.8421
3	22.8	$3 \times 10^5$	111, 105, 117	111	0.8082
4	31.5	$3 \times 10^5$	98, 93, 96	95.66	0.6966
5					
6					
7					
8					