

Investigator: A. Bishayee

Granulocyte Macrophage-Colony Forming Unit (GM-CFU) Assay

Experiment # 9 Serial #2 C-157 Source of Irradiation: External

Mice Sex, Strain, Age: SW, F, 5-6 wk

Irradiation

Type of Irradiation: Chronic External

Animals per group: 3

Aims: To determine the bone marrow GM-CFC response to chronic external ¹³⁷Cs irradiation (dose rate reduction half-life corresponds to T_{1/2} of Sr-117m = 223h)

Summary of Results:

Brief Procedure:

- 1) Irradiate animals in groups of 3 with desired initial dose rates (cGy/h) listed in Table 1.
- 2) Sacrifice each mouse on optimal day by cervical dislocation and sterilize using 70% EtOH and immediately move it into laminar flow hood.
- 3) Remove both femurs carefully using sterile instruments and clean the attached tissue thoroughly.
- 4) Flush the bone marrow with 2% Horse Serum in Dulbecco's Modified Eagles Medium (2% HS-DMEM) using 21G needle and syringe.
- 5) Separate the mononuclear cells by density gradient procedures using Histopaque.
- 6) Plate the desired number of in mixture of 60% HS-DMEM and 0.6% bacto agar solution in the presence of 9.2 U (New Sigma unit) GM-CFS.
- 7) Keep the plated well plates for 20 min. in laminar flow hood and move them into incubator with 5% CO₂ and 95% air, at 37°C.
- 8) Count the granulocyte macrophage colonies on 7th day.

Group # Probe #	Starting date	Initial Dose Rate R/h	Date Sacrificed	# of days	Remarks if any
C1	10/7/98	0	10/14/98	7	} unirradiated control
C2	10/7/98	0	10/14/98	7	
C3					
1	10/7/98	3.0	10/14/98	7	Cage #1
2	10/7/98	1.506	10/14/98	7	Cage #2
3	10/7/98	0.825	10/14/98	7	Cage #3
4	10/7/98	0.456	10/14/98	7	Cage #4
5					
6					
7					

Preparing Media and Agar:

Culture Medium (Double Strength): 13.37g (1 pack) of D-MEM powder (Gibco, Cat # 12100-046) + 490 ml deionized water + 16 μ l of L-asparagine (Gibco Cat # 12416-012) at a concentration of 5 μ g/ μ l + 150.4 μ l of DEAE dextran (mol. wt. = 2×10^6 , intrinsic viscosity = 0.7) at a concentration of 1 μ g/ μ l (Sigma Cat # D-9885) + 10 ml of penstrep (Gibco Cat # 600-5070, 5,000 units/ml pen, 5,000 μ g/ml streptomycin) + 3.7 g of NaHCO₃ (Gibco Cat # 11810-025).

Wash Medium: i) Mix equal amounts of culture medium and sterile deionized water.
ii) Add 2% HS

Agar: Prepare 0.6% agar by adding 0.6 g Difco Bacto agar (Difco Cat # 0140-15-4) to 100 ml boiling deionized water. Autoclave on liquid cycle for 20 min.

Comments If any:**Flushing Bone marrow:**

- 1) Remove both femurs from each mouse and place them in a test tube containing wash medium kept in ice, if the femur can not be flushed immediately.
- 2) Flush the marrow from each femur by aspirating 3 ml of Wash Medium through the femur 5 times with a 21G needle/3 ml syringe in a 50 ml conical centrifuge tube. Follow with two flushes with 1 ml of fresh Wash Medium.
- 3) Spin the cells at 1200 rpm for 5 minutes at 4°C, decant, break up the cell pellet, resuspend the cells in 5 ml of cold Wash Medium, and vortex the cell suspension.

Comments If any:**Separating Mononuclear cells and washing the cells:**

- 1) Transfer 3.5 ml of Histopaque (Sigma Cat #H8889) into fresh 15 ml tubes (1 tube per group).
- 2) Layer the cell suspension carefully on top of Histopaque and centrifuge at 1500 rpm, 4°C, for 30 minutes.
- 3) Using a Pasteur pipette transfer the mononuclear cells into fresh tubes.
- 4) Dilute the cell suspension to 15 ml by adding cold Wash Medium into each tube and spin them at 1200 rpm, 4°C, for 5 min.
- 5) Decant the supernatant, break the pellet, and add 15 ml cold Wash Medium, and spin them again at 1200 rpm, 4°C, for 5 min. Repeat this procedure 2 more times.
- 6) After 3rd wash break the pellet and resuspend in 2 ml Culture medium (2x DMEM) with 60% HS and keep the tubes in dry bath at 37°C.
- 7) Add 20 μ l of cell suspension to 20 ml of Isotone II in a coulter cup and determine total # of cells in each group using coulter counter.

Coulter Counter Parameters:Current(I)=500 μ A

Full Scale = 1

T_L = 2.7T_u = 99.9

Attenuation= 4

Alarm Threshold = off

Preset Gain = 1

Stirrer control = off

Coulter Counter Parameters: *Same as above*

Multiplication Factor to get # of cells/ μ l = $2000 \times$ Coulter Count

Group #	Coulter Count without ZG	Avg	# cells per μ l	Coulter Count with 5 drops ZG	Avg	# cells per μ l
C1				3830, 3833, 3740	} 1014	9629666
C2				5729, 5881, 5876		
C3						
1				2568, 2444, 2384		
2				5224, 5234, 5053		
3				3412, 3363, 3374		
4				3116, 3086, 3028		
5						
6						
7						

Average # of cells per $\mu\text{l} = \mu\text{l} = 9,629,666$

(4)

DILUTIONS

Dilution A: (1.0×10^6 cells / ml, Total volume 3.4 ml)

1.7 ml Agar + 1.34 μl Medium + 0.353 μl Cell Suspension

Dilution B: (3.0×10^5 cells / ml, Total volume 3.4 ml)

1.7 ml Agar + 1.6 μl Medium + 0.106 μl Cell Suspension

Dilution C: (1.0×10^5 cells / ml, Total volume 3.4 ml)

1.7 ml Agar + μl Medium + μl Cell Suspension

Plating the Cells:

Culture Medium: Maintain four 13mm tubes each containing 4.5 ml of Culture medium in dry bath at 37°C.

Horse Serum: Maintain five 13mm tubes each containing 4.5 ml of Horse Serum in dry bath at 37°C.

Agar: Maintain five 16mm tubes each containing 6.5 ml of Agar in dry bath at 37°C.

- 1) Warm up dilution tubes (one per group) to 37°C in dry bath.
- 2) Warm up Agar (30 ml) and 60% HS in 2x DMEM (30 ml) in separate tubes to 37°C.
- 3) Mark the well plates (3 wells for each dilution for each group)
- 4) Mix 1.7ml agar + x ml of 2x DMEM with 60% HS + y ml cell suspension + 0.02 ml GMCSF (x + y = 1.7 ml) in a dilution tube.
- 5) Add 1 ml of mixture 4 to each well, mix properly and let it gel for about 15 minutes.
- 6) Repeat steps 4 and 5 for each dilution.
- 7) Repeat steps 1 to 6 for each group.
- 8) Incubate the cells in an incubator at 37°C and 5% CO₂, 95% air for 7 days.
- 9) On 8th day of incubation count colonies and determine the survival fractions.

Comments If any:

Counting the Colonies: (Inverted at 40X or dissecting at 35X)

10/21/98

⑤

Group #	Dose (p) Activity # injected	# of cells plated	# CFU-GM counted	Avg	SF
C1	0	3×10^5	112, 114, 112	126.16	
C2	0	3×10^5	131, 141, 147		
C3					
1	390.5	1×10^6	30, 42, 28	10.10	0.0800
2	196	3×10^5	36, 45, 42	41	0.3249
3	107	3×10^5	52, 47, 39	46	0.3646
4	59	3×10^5	66, 78, 85	76	0.6050
5					
6					
7					