

Granulocyte Macrophage-Colony Forming Unit (GM-CFU) Assay

Experiment # 2 (Acute # 2)  
 Mice Age, Sex, Strain: 10-11 wks, F, SW  
 Type of Irradiation: Acute  
 Animals per group: 3

Source of Irradiation: <sup>137</sup>Cs HDR1 (Mark I)  
 Date In Irradiator: 08/11/98  
 Date Out Irradiator: 08/11/98  
 Date Sacrificed: 08/11/98

Aim: *To study GM-CFU survival against acute irradiation with three dose points*

Results:Brief Procedure :

- 1) Irradiate the mice acutely with desired doses.
- 2) Sacrifice each mouse by cervical dislocation and sterilize using 70% EtOH and move it into laminar flow hood.
- 3) Remove both femurs carefully using sterile instruments and clean the attached flesh thoroughly.
- 4) Flush the bone marrow with 2% Horse Serum in Dulbecco's Modified Eagles Medium (2% HS-DMEM) using 21G needle and syringe.
- 5) Separate the mononuclear cells by density gradient procedures using Histopaque.
- 6) Plate the desired number of cells (cell suspension) in mixture of 60% HS-DMEM and 1.7 ml 0.6% bacto agar solution in the presence of 9.2 U GM-CFS in six-well plates.
- 7) Keep the plates for 20 min. in laminar flow hood and move them into incubator with 5% CO<sub>2</sub> and 95% air, at 37°C.
- 8) Count the granulocyte macrophage colonies on 7th day.

High Dose Rate Irradiation:

Group #	Attenuator Used	Turn Table Position	# of Sources Used	Dose Rate (R/min)	Irradiation Time (min)	Total Dose (R)	Comments if any
C1	—					0	unmarked
C2							
1	X-10	#3	2	101.4	0.98	100	marked on head
2	X-10	#3	2	101.4	1.97	200	" " Body
3	X-10	#3	2	101.4	3.94	400	" " Tail
4							
5							
6							

**Preparing Media and Agar:**

**Culture Medium (Double Strength, 2X):** 13.37g (1 pack) of D-MEM powder (Gibco, Cat # 12100-038) + 490 ml deionized water + 16  $\mu$ l of L-asparagine (Gibco Cat # 12416-012) at a concentration of 5  $\mu$ g/ $\mu$ l + 150.4  $\mu$ l of DEAE dextran (mol. wt. =  $2 \times 10^6$ , intrinsic viscosity = 0.7) at a concentration of 1  $\mu$ g/ $\mu$ l (Sigma Cat # D-9885) + 10 ml of penstrep (Gibco Cat # 600-5070, 5,000 units/ml pen, 5,000  $\mu$ g/ml streptomycin) + 3.7 g of  $\text{NaHCO}_3$  (Gibco Cat # 11810-025).

✓ **Culture Medium (2X) with 60 % Horse Serum:** Add 60 % Horse Serum in 2X DMEM

**Wash Medium:** i) Mix equal amounts of culture medium and sterile deionized water.  
ii) Add 2% HS

**Agar:** Prepare 0.6% agar by adding 0.6 g Difco Bacto agar (Difco Cat # 0140-15-4) to 100 ml boiling deionized water. Autoclave on liquid cycle for 20 min.

**Comments If any:**

**Flushing Bone marrow:**

- 1) Remove both femurs and tibias from each mouse and place them in a test tube containing wash medium kept in ice if the femur can not be flushed immediately.
- 2) Flush the marrow from each femur by aspirating 3 ml of Wash Medium through the femur 10 times with a 21G needle/3 ml syringe in a 50 ml conical centrifuge tube. Follow with two flushes with 0.5 ml of fresh Wash Medium.
- 3) Spin the cells at 1200 rpm for 5 minutes at 4°C, decant, break up pellet, resuspend the cells in 5 ml of cold Wash Medium, and vortex the cell suspension.

**Comments If any:**

**Counting the Cells:**

Add 10 $\mu$ l of cell suspension to 20 ml of Isotone II in a coulter cup and count the cells using coulter counter. Calculate total # of cells in each group.

**Coulter Counter Parameters:**Current(I)=500  $\mu$ A

Full Scale = 1

T<sub>L</sub> = 2.7T<sub>u</sub> = 99.9Manometer Select = 500  $\mu$ l

Attenuation = 4

Alarm Threshold = off

Preset Gain = 1

Stirrer control = off

*Multiplication Factor to get total # of cells in 5 ml = 20,000 x Coulter count*

Group #	Coulter Count without ZG	Avg	Total # of cells	Coulter Count with 2 drops ZG	Avg	Total # of cells
C1	NOT PERFORMED					
C2						
1						
2						
3						
4						
5						
6						

Comments if any:

**Separating Mononuclear cells and washing the cells:**

- 1) Transfer 3.5 ml of Histopaque (Sigma Cat #H8889) into fresh 15 ml tubes (1 tube per group).
- 2) Layer the cell suspension carefully on top of Histopaque and centrifuge at 1500 rpm, 4°C, for 30 minutes.
- 3) Using a Pasteur pipette transfer the mononuclear cells into fresh tubes.
- 4) Dilute the cell suspension to 15 ml by adding cold Wash Medium into each tube and spin them at 1200 rpm, 4°C, for 5 min.
- 5) Decant the supernatant, break the pellet, and add 15 ml cold Wash Medium, and spin them again at 1200 rpm, 4°C, for 5 min. Repeat this procedure 2 more times.
- 6) After 3rd wash break the pellet and resuspend in 2 ml Culture medium (2x DMEM) with 60% HS and keep the tubes in dry bath at 37°C.
- 7) Add 20  $\mu$ l of cell suspension to 20 ml of Isotone II in a coulter cup and determine total # of cells in each group using coulter counter.

Exp. # \_\_\_\_\_  
Coulter Counter Parameters: *Same as above*

*Multiplication Factor to get # of cells/ml = 2000x Coulter Count*

Group #	Coulter Count without ZG	Avg	# cells per $\mu$ l	Coulter Count with 2 drops ZG	Avg	# cells per $\mu$ l
C1						
C2						
1						
2						
3						
4						
5						
6						

*Please See the separate sheet*

Average # of cells per  $\mu$ l =

DILUTIONS

Culture Medium: Maintain four 13mm tubes each containing 4.5 ml of Culture medium in dry bath at 37°C.

Horse Serum : Maintain five 13 mm tubes each containing 4.5 ml of Horse Serum in dry bath at 37°C.

Agar: Maintain five 16 mm tubes each containing 6.5 ml of Agar in dry bath at 37°C.

- 1) Warm up dilution tubes (three per group) to 37°C in dry bath.
- 2) Warm up Agar (30 ml) and 60% HS in 2x DMEM (30 ml) in separate tubes to 37°C.
- 3) Mark the six-well plates (three wells for each dilution) containing 20 ul of GM-CSF (9.2 U) in each well.
- 4) Mix 1.7ml agar + x ml of 2x DMEM with 60% HS + y ml cell suspension + 0.02 ml GMCSF (x + y = 1.7 ml) in a dilution tube.
- 5) Add 1 ml of mixture 4 to each well, mix properly and let it gel for about 15 minutes.
- 6) Repeat steps 4 and 5 for each dilution.
- 7) Repeat steps 1 to 6 for each group.
- 8) Incubate the cells in an incubator at 37°C and 5% CO<sub>2</sub>, 95% air for 7 days.
- 9) On 8th day of incubation count colonies and determine the survival fractions.

Comments If any:

Counting the Colonies: (Inverted at 40X or dissecting at 35X)

Group #	Dose (rads)	# of cells plated	# CFU-GM counted	Avg	SF
C1	0	3x10 <sup>5</sup>	134, 133, 128	131.66	
C2					
1	100	3x10 <sup>5</sup>	52, 47, 40	46.33	0.3519
2	200	3x10 <sup>5</sup>	18, 15, 13	15.33	0.1164
3	400	1x10 <sup>6</sup>	9, 10, 11	10	0.022
4					
5					

Cell count

Group	Coulter's count for 20ul Cells (MS = 500ul)	Avg. count	Cells/ml	Total # of cells	Cells/penar
1	3382, 3234, 3202	3272	6,545,333	13090666	2,181,777
2	2717, 2600, 2620	2645	5,291,333	10592666	1,763,777
3	3810, 3745, 3719	3758	7,516,000	15032000	2,505,333
4	4507, 4409, 4343	4419	8,839,333	1,7678666	2,946,444

## Dilution

### Group 1

Dilution A : 3.4 ml of 1,000,000 cells/ml = 3,400,000 cells

$$\text{Vol. required} = \frac{3,400,000}{6545333} = 0.519 \text{ ml}$$

For dilution,

1.7 ml Agar + 1.18 ml ~~2x~~ DMEM + 0.519 <sup>ml</sup> ~~ml~~ cells.

Dilution B : 3.4 ml of 300,000 cells/ml = 1,020,000 cells

$$\text{Vol. required} = \frac{1,020,000}{6545333} = 0.155 \text{ ml}$$

For dilution,

1.7 ml Agar + 1.54 ml ~~2x~~ DMEM + 0.155 <sup>ml</sup> ~~ml~~ cells

## Group 2

### Dilution A

$$\text{Vol. required} = \frac{3400000}{5291333} = 0.642 \text{ ml}$$

1.7 ml Agar + 1.05 2x DMEM + 0.642 ml Cells

### Dilution B

$$\text{Vol. required} = \frac{1,020,000}{5291333} = 0.192 \text{ ml}$$

1.7 ml Agar + 1.50 ml 2x DMEM + 0.192 ml Cells

## Group 3

### Dilution A

$$\text{Vol. required} = \frac{3400000}{7516000} = 0.452 \text{ ml}$$

1.7 ml Agar + 1.247 ml 2x DMEM + 0.452 ml Cells

### Dilution B

$$\text{Vol. required} = \frac{1,020,000}{7516000} = 0.135 \text{ ml}$$

1.7 ml Agar + 1.564 ml 2x DMEM + 0.135 ml Cells



Group 4

Dilution A      vol. required =  $\frac{3400000}{8839333} = 0.384 \text{ ml}$


1.7 ml Agar + 1.315 ml 2x DMEM + 0.384 ml Cells

Dilution B      vol. required =  $\frac{1,020,000}{3400000} = 0.115 \text{ ml}$   
8839333

1.7 ml Agar + 1.58 ml 2x DMEM + 0.115 ml Cells

EXPE #2

D37 = 100 R



CITADEL® NO. 652 - SEMI-LOG 2 X 10

330

Dose (Rad)

500