

Investigator: A. BISHAYEE

Granulocyte Macrophage-Colony Forming Unit (GM-CFU) Assay

Experiment # 1 (Acute #1) Source of Irradiation: ¹³⁷CS HDR (Mark-I)
 Mice Age, Sex, Strain: 8-9 wks, F, SW Date In Irradiator: 07/29/98
 Type of Irradiation: Acute Date Out Irradiator: 07/29/98
 Animals per group: 3 Date Sacrificed: 07/29/98

Aim: *To study GM-CFU survival against acute irradiation.*

Results: *GM-CFU survival is inversely proportional with acute irradiation (100-250 R)*

Brief Procedure:

- 1) Irradiate the mice acutely with desired doses.
- 2) Sacrifice each mouse by cervical dislocation and sterilize using 70% EtOH and move it into laminar flow hood.
- 3) Remove both femurs carefully using sterile instruments and clean the attached flesh thoroughly.
- 4) Flush the bone marrow with 2% Horse Serum in Dulbecco's Modified Eagles Medium (2% HS-DMEM) using 21G needle and syringe.
- 5) Separate the mononuclear cells by density gradient procedures using Histopaque.
- 6) Plate the desired number of cells (~~400~~ cell suspension) in mixture of ~~60~~ 60% HS-DMEM and ~~40~~ 0.6% bacto agar solution in the presence of ~~100~~ U GM-CFS in 6-well plate
- 7) Keep the plates ~~in petri dishes~~ for 20 min. in laminar flow hood and move them into incubator with 5% CO₂ and 95% air, at 37°C.
- 8) Count the granulocyte macrophage colonies on 7th day.

High Dose Rate Irradiation:

Group#	Attenuator Used	Turn Table Position	# of Sources Used	Dose Rate R/min	Irradiation Time (min)	Total Dose (R)	Comments if any
✓ C1	-	-	-	-	-	0	untreated control
C2							
✓ 1	X-10	#3	2	101.4	0.98	100	Marked on head
✓ 2	X-10	#3	2	101.4	2.46	250	Marked on body
3							
4							
5							
6							

Preparing Media and Agar:

38 Culture Medium (Double Strength): 13.37g (1 pack) of D-MEM powder (Gibco, Cat # 12100-046) + 490 ml deionized water + 16 μ l of L-asparagine (Gibco Cat # 12416-012) at a concentration of 5 μ g/ μ l + 150.4 μ l of DEAE dextran (mol. wt. = 2×10^6 , intrinsic viscosity = 0.7) at a concentration of 1 μ g/ μ l (Sigma Cat # D-9885) + 10 ml of penstrep (Gibco Cat # 600-5070, 5,000 units/ml pen, 5,000 μ g/ml streptomycin) + 3.7 g of NaHCO_3 (Gibco Cat # 11810-025).

Culture medium (2x) with 60% Horse Serum: Add 60% Horse serum in 2x DMEM just before use.

Wash Medium: i) Mix equal amounts of culture medium and sterile deionized water.
ii) Add 2% HS just prior to experiment.

Agar: Prepare 0.6% agar by adding 0.6 g Difco Bacto agar (Difco Cat # 0140-15-4) to 100 ml boiling deionized water. Autoclave on liquid cycle for 20 min.

Comments If any: Horse serum should be added on the previous day of the experiment.

Flushing Bone marrow:

- 1) Remove both femurs and tibias from each mouse and place them in a test tube containing wash medium kept in ice if the femur can not be flushed immediately.
- 2) Flush the marrow from each femur by aspirating 3 ml of Wash Medium through the femur 10 times with a 21G needle/3 ml syringe in a 50 ml conical centrifuge tube. Follow with two flushes with 0.5 ml of fresh Wash Medium.
- 3) Spin the cells at 1200 rpm for 5 minutes at 4°C, decant, break up pellet, resuspend the cells in 5 ml of cold Wash Medium, and vortex the cell suspension.

Comments If any:

Counting the Cells:

Add 10 μ l of cell suspension to 20 ml of Isotone II in a coulter cup and count the cells using coulter counter. Calculate total # of cells in each group.

Coulter Counter Parameters:Current(I)=500 μ A

Full Scale = 1

 T_L = 2.7 T_U = 99.9Manometer Select = 500 μ l

Attenuation= 4

Alarm Threshold = off

Preset Gain = 1

Stirrer control = off

Multiplication Factor to get total # of cells in 5 ml = 20,000 x Coulter count

Group #	Coulter Count without ZG	Avg	Total # of cells	Coulter Count with 2 drops ZG	Avg	Total # of cells
C1	Not performed					
C2						
1						
2						
3						
4						
5						
6						

Comments if any:

Separating Mononuclear cells and washing the cells:

- 1) Transfer 3.5 ml of Histopaque (Sigma Cat #H8889) into fresh 15 ml tubes (1 tube per group).
- 2) Layer the cell suspension carefully on top of Histopaque and centrifuge at 1500 rpm, 4°C, for 30 minutes.
- 3) Using a Pasteur pipette transfer the mononuclear cells into fresh tubes.
- 4) Dilute the cell suspension to 15 ml by adding cold Wash Medium into each tube and spin them at 1200 rpm, 4°C, for 5 min.
- 5) Decant the supernatant, break the pellet, and add 15 ml cold Wash Medium, and spin them again at 1200 rpm, 4°C, for 5 min. Repeat this procedure 2 more times.
- 6) After 3rd wash break the pellet and resuspend in 2 ml Culture medium (2x DMEM) with 60% HS and keep the tubes in dry bath at 37°C.
- 7) Add 20 μ l of cell suspension to 20 ml of Isotone II in a coulter cup and determine total # of cells in each group using coulter counter.

Multiplication Factor to get # of cells/ μ l = $2 \times$ Coulter Count (# of cells/ml = $2000 \times$ Coulter Count)

Group #	Coulter Count without ZG	Avg	# cells per μ l	Coulter Count with 2 drops ZG	Avg	# cells per μ l	cells/ sample
C1	Not performed ↓			2264, 2118, 2205	2195	4,391,333	1,463,777
C2							
1				533, 551, 535*	539	1,079,333	359,777
2				3376, 3344, 3256	3325	6,650,666	2,216,888
3							
4							
5							
6							

Average # of cells per μ l =

DILUTIONS

* This low count may be due to problem associated with density gradient procedure

For dilution — see the attached sheets

Culture Medium: Maintain four 13mm tubes each containing 4.5 ml of Culture medium in dry bath at 37°C.

~~Horse Serum: Maintain five 13mm tubes each containing 4.5 ml of Horse Serum in dry bath at 37°C.~~

~~Agar: Maintain five 16mm tubes each containing 6.5 ml of Agar in dry bath at 37°C.~~

- 1) Warm up dilution tubes (^{three}one per group) to 37°C in dry bath.
- 2) Warm up Agar (~~30 ml~~) and 60% HS in 2x DMEM (~~30 ml~~) in separate tubes to 37°C
- 3) Mark the petri dishes (^{well plate}3 petri dishes ^{3 wells}per group for each dilution). ^{containing 20 ul of GM-CSF} 20 ul
- 4) Mix 1.7ml agar + x ml of 2x DMEM with 60% HS + y ml cell suspension + 0.408 ml GMCSF (x + y = 1.7 ml) in a dilution tube.
- 5) Add 1 ml of mixture 4 to each ^{well plate}petri dish, mix properly and let it gel for about 15 minutes.
- 6) Repeat steps 4 and 5 for each ^{dilution}group.
- 7) Repeat steps 1 to 6 for each ^{dilution}group.
- 8) Incubate the cells in an incubator at 37°C and 5% CO₂, 95% air for 7 days.
- 9) On 8th day of incubation count colonies and determine the survival fractions.

Comments If any:

Agar was solidified at 37°C, no stock bottle was used following occasional heating at microwave oven

Counting the Colonies: (Inverted at 40X or dissecting at 35X)

8/5/98

Group #	Dose (rads)	# of cells plated	# CFU-GM counted	Avg	SF
C1	0	1x10 ⁶	187, 192, 179	186	
C2 C	0	3x10 ⁵	82, 74, 91	82.33	-
C3	0	1x10 ⁵	10, 12, 15	12.33	
2	100	3x10 ⁵	23, 21, 27	23.6	0.28
3	100	1x10 ⁵	7, 8, 5	6.66	
4	250	1x10 ⁶	31, 43, 50	41.33	
5	250	3x10 ⁵	10, 8, 7	8.33	0.10
6	250	1x10 ⁵	4, 3, 2	3	

Cell Count

①

①

couller count for
20ml E zapo-Albumin
2264, 2118, 2205

AVG.
Count
2195.6

Cells/ml
4,391,333

Total # of
cells
8782666

cells/kaur
1,463,777

②

533, 551, 535

539.6

1,079,333

2158666

359777

③

3376, 3344, 3256

3325.3

6,650,666

13301332

2,216,888

Dilution

Group 1 :

(control)

Coulter Count for 20 μ l (MS = 500 μ l) 2264, 2118, 2205	Cells/ml 4,391,333
--	-----------------------

1. 9 x 13 mm glass tube (for dilution)
2. 3 x 13 mm \bar{E} 5.5 ml Agar in each
3. 3 x 13 mm \bar{E} 5.5 2X DMEM
4. 5 x 6-well plate \bar{E} GM-CSF (20 μ l)

Dilution A : 3.4 ml of 1,000,000 cells/ml = 3,400,000 cells

$$\text{Vol. required} = \frac{3400000}{4,391,333} = 0.774 \text{ ml}$$

For dilution, 1.7 ml agar + 774 μ l cells + 925 μ l 2X DMEM

Dilution B : 3.4 ml of 300,000 cells/ml = 1,020,000 cells

$$\text{Vol. required} = \frac{1,020,000}{4,391,333} = 0.232 \text{ ml}$$

For dilution, 1.7 ml agar + 232 μ l cells + 1.46 ^{ml} ~~ml~~ 2X DMEM

Dilution C : 3.4 ml of 100,000 cells/ml = 340,000 cells

$$\text{Vol. required} = \frac{340,000}{4,391,333} = 0.077 = 77 \mu\text{l}$$

For dilution, 1.7 ml agar + 77 μ l cells + 1.62 ^{ml} ~~ml~~ 2X DMEM

Gr 2
(100R)

3.4 ml of 1,000,000 cells/ml = 3,400,000 cells

Dil. A : vol required = $\frac{3400000}{1079333} = 3.15 \text{ ml}$

This dilution was not possible because total the volume of final cell suspension was ~ 2 ml.

Dil. B : vol required = $\frac{1020000}{1079333} = 0.945 \text{ ml}$

1.7 ml agar + 0.945 ml cells + 0.755 ml 2X DMEM

Dil. C : vol required = $\frac{340000}{1079333} = 0.315 \text{ ml}$

1.7 ml agar + 0.315 ml cells + 1.38 ml 2X DMEM

Gr 3 (250R)

Dil. A vol = $\frac{3400000}{6650666} = 0.511 \text{ ml}$

1.7 ml agar + 0.511 ml cells + 1.188 ml 2X DMEM

Dil B vol = 0.153 ml $\left(\frac{1020000}{6650666} = 0.153 \right)$

1.7 ml agar + 0.153 ml cells + 1.546 ml 2X DMEM

Dil C vol = 0.05 ml = 51 $\left(\frac{340000}{6650666} = 0.051 \right)$

1.7 ml agar + 51 μ l cells + 1.649 ml 2X DMEM

10% → 20X
90%

1% → 200X
99%

control

100


200

400

1.0

Survival Fraction

0.01



CITADEL® NO. 652 - SEMI-LOG 2 X 10

Dose (Rad)