

WB Cells
Stacey

WB Protocol

Date: 6/1/00
Exp. Name

Ex.#
Name:

1. Set Rocker-Roller in 37°C incubator with 5% CO₂
-turn on coulter counter or have a hemacytometer ready
-turn centrifuge on so it can start to cool
2. Wash 2-150 cm² flasks with 5ml of PBS twice and then trypsinize the cells using 3 ml of trypsin and place flasks in incubator for 3-4 min. until cells lift when flask in "smacked"
3. Neutralize the trypsin by adding 7 ml of D-Medium and resuspend the cells by pooling into a 50 ml tube (total vol.=10 ml=3 ml trypsin+7 ml D-Medium)
4. Centrifuge the tube at 2000 rpm, 4°C, 10 min.

etc. from and move to the other plate
 2-WBs flask were to D-Medium
 1-WB flask washed to 16ml D-Medium

5. Decant supernatant, gently break pellet by running tube along a test tube rack 3 times
-resuspend pellet in 10 ml of D-Medium and vortex
6. gently syringe cells 5 times using a 5-cc syringe with a 21 gauge needle
7. Do cell count using coulter counter

-100ul of cells in 20 ml of Isotone

-do 3 times:

1.	4478	1. 2387
2.	4541	2. 2375
3.	4428	3. 2395

need - 6 tubes 2mls each = 12mls
 2×10^6 cells/ml
 = 24 mill cells
 But need to overstock
 50 more / 5mls
 = 30 mill cells

background 00010
20002500ul

(if setting is at 50 ul on coulter counter multiply ave. by 4000 or if setting is at 500 ul on coulter counter then multiply ave. by 400)

CELL COUNT:

$4462 \times 400 = 1,784,800$

$2385 \times 400 = 954,000$

8. Dilute cells to 2 X 10⁶ cells/ml in D-Medium

2 x 10⁶ cells/ml in memB.

FINAL VOL: 10.0

avg 49mls

ACTUAL CELL COUNT:

$4 \text{ tubes} = 8 \text{mls} \times 2 \times 10^6 = 16,000,000$
 2mls per tube

9. Transfer 2 ml (4 X 10⁶ cells/ml) of cell suspension into six 14 ml Falcon polypropylene tubes and label tube and cap

place 1 tube in centrifuge, clean
remove pellet & resuspend in 2ml memB.

10. Place tubes in Rocker-Roller for 12 hours at 37°C, 5% CO₂

-syringe, then place rest of tubes in rocker roller

DATE/TIME:

6/1/00
 3:40pm

tubes 1, 3, 4 - look good
 no clumping

WBS

Clive Ct 4462/

$$\begin{aligned} \text{Cell Conc} &= 4462 \times 400 \\ &= 1,784,800 \text{ cell/mL} \end{aligned}$$

Need 15 mL of 2,000,000 cells/mL
= 30,000,000 cells

$$\text{Vol. Required} = \frac{30,000,000}{1,784,800} = 16.8 \text{ mL}$$

we need 9 mL of 2,000,000 cells/mL
= 18,000,000 cells

$$\text{Vol. required} = \frac{18,000,000}{1,784,800} = 10.08$$

1 Dm | 1 - mLms

WBC

ave # 2385

$$\text{cell conc} = 2385 \times 400 = 954,000$$

$$\text{Need} = 6 \text{ tubes}, 2 \text{ ml/s each} = 12 \text{ ml/s} \\ 2 \times 10^6 \text{ cell/ml} = 24 \text{ mil cells}$$

$$\text{Need to make more } 15 \text{ ml/s}, 2 \times 10^6 \text{ cell/ml} \\ = 30 \text{ mil cells}$$

Vol. Req.

$$\frac{30,000,000}{954,000} = 31.4 \text{ ml/s}$$

|| D-men

~~9~~ ml/s

9 ml/s of cells
 $2 \times 10^6 \text{ cell/ml}$

Vol
Req.

$$\frac{18,000,000}{954,000} =$$

Tube 1 & 2 6/2/00

Day 2: condition 1 (wash and plate)

- remove pellets from before roller
- 1. Wash cells with 5 ml of D-Medium and centrifuge at 2000 rpm, 10 min. 4°C
- 2. Decant and break pellet as previously described, vortex and resuspend cells in 2 ml of D-Medium using a 2 ml pipet

Resuspended
COO11
50ul

- 3. Syringe cells 5 times using a 3 cc syringe with a 21 gauge needle

3589

-DO CELL COUNT: Tube 1 2719

Tube 2 2627

- 4. Do 4, 10-fold dilutions 2506

2698

- set up 8 dilution tubes: (each condition must be done in duplicate)
- place 4.5 ml of D-Medium into the 2 sets of 4 tubes and label the tubes 10^{-1} , 10^{-2} , 10^{-3} , and 10^{-4} and also label the dilution tube with the corresponding tube number which contains the cells
- 5. Perform the dilutions by placing 0.5ml of cells into the dilution tubes 10^{-1} .

- vortex the tube.
- then take 0.5ml of the cells from that dilution tube and place it into the next tube, 10^{-2}
- vortex the tube
- then take 0.5ml of the cells from that dilution tubes and place into the next tube, 10^{-3}
- vortex the tube
- finally do the same for the last dilution tube

2ml D-Medium in Dishes

- 6. Seed in triplicate, 1ml of the diluted cells into a 60 mm culture dish, label dish
- Use dilutions- 10^{-3} and 10^{-4} *repeat for next 14ml tube
- 7. Incubate the dishes for 7 days in a 37°C, 5% CO₂ incubator

DATE: 6/2/00 Plate 9:00

- 7 days later:
- 8. Remove culture dishes from incubator
- 9. Remove media and label underside of culture dishes selected for count (50-250 colonies)
- 10. Wash 3 times with normal saline and then 2 times with methanol
- 11. Stain colonies with 0.1% crystal violet and count colonies under fluroscent light

Day 2: condition 2 (clusters) Tube 3

use 3

pellet

1. Remove tubes from the rocker-roller
2. Wash with 5 ml of D-Medium and centrifuge at 2000 rpm, 10 min., 4°C
3. Decant supernatant and break pellet, as described before, vortex and place cell suspension (200ul) into a Helina tube using a 200ul pipet tip
4. Wash 14 ml test tube with 200 ul of D-Medium and place the wash into the Helina tube using a 200ul pipet tip to get a final vol. of 400ul in the Helina tube
5. Centrifuge the Helina tubes at 1000 rpm, 4°C, 5 min.
6. Place the tubes in Rainin pipet rack and then place into refrigerator at 10.5°C for 72 hours

DATE: 6/2/00 10:00

Day 5: 72 hours later: Date:

1. After 72 hour incubation, remove the Helina tubes from the refrigerator
 - set up 2, 14ml tubes with 10 ml of D-Medium in them and label 1 and 2
 - remove supernatant from the Helina tube and place into corresponding 14 ml tube
 - then break pellet in Helina tube and place the cells into the 14 ml tube
 - remove some medium from the top of the 14 ml tube and place it into the Helina tube and wash the Helina tube to collect any remaining cells and place back into the 14 ml tube
2. centrifuge the tubes in a precooled centrifuge for 10 min, 2000 rpm, 4°C
3. While waiting for centrifuge, Label 2 sets of 4 glass dilution tubes with the following dilutions -- 10^{-1} , 10^{-2} , 10^{-3} , and 10^{-4} and labe the dilution tube with the corresponding 14 ml tube #
 - place 4.5 ml of D-Medium into each of the dilution tubes
4. When done centrifuging, decant the supernatant and break the pellet as previously described and then vortex.
 - resuspend the cells in 2 ml of D-Medium with a 2 ml pipet
 - then syringe the cells 5 times using a 3cc syringe with a 21 gauge needle

DO CELL COUNT:

5. Do four-10 fold dilutions

Perform the dilutions by placing 0.5ml of cells into the dilution tube 10^{-1} .

-vortex the tube.

-then take 0.5ml of the cells from that dilution tube and place it into the next tube, 10^{-2}

-vortex the tube

-then take 0.5ml of the cells from that dilution tubes and place into the next tube, 10^{-3}

-vortex the tube

-finally do the same for the last dilution tube

6. Seed in triplicate, 1ml of the diluted cells into a 60 mm culture dish, label dish

-Use dilutions- 10^{-3} and 10^{-4}

-Repeat for next 14 ml tube

7. Incubate the dishes for 7 days in a 37°C , 5% CO_2 incubator

8. Following staining and counting procedure as described above.

tube 4

Day 2, condition 3: incubate 72 hours 10.5°C in rocker-roller

1. Remove cells from rocker-roller(37°C)
2. Wash cells with 5 ml of D-Medium and centrifuge tubes at 2000 rpm, 10 min. 4°C
3. Decant the supernatant and break pellet as previously described, vortex collect = 1407
-resuspend the cells in 2 ml of D-Medium and syringe 3 times using a 3cc syringe with a 21/32 3 gauge needle
4. Place tubes in a rocker-roller which is in a 10.5°C refrigerator for 72 hours.

DATE:

10:08

6/2/00

Day 5: 72 hours later:

1. After the 72 hours, remove the tubes from the rocker-roller
 2. add 5 ml of D-Medium and centrifuge tubes at 2000 rpm, 10 min, 4°C
 3. Decant the supernatant, break pellet, vortex, and resuspend the cells in 2 ml of D-Medium
 4. syringe the cells with a 3cc syringe, 21 gauge needle, 5 times
 5. Do four, 10-fold dilutions
- label 2 sets of 4 glass dilution tubes as follows:
10⁻¹, 10⁻², 10⁻³, and 10⁻⁴ and label the dilution tube with the corresponding 14 ml tube #
 - place 4.5 ml of D-Medium into each of the dilution tubes

*Perform the dilutions by placing 0.5ml of cells into the dilution tube 10⁻¹.

-vortex the tube.

-then take 0.5ml of the cells from that dilution tube and place it into the next tube, 10⁻²

-vortex the tube

-then take 0.5ml of the cells from that dilution tubes and place into the next tube, 10⁻³

-vortex the tube

-finally do the same for the last dilution tube

6. Seed in triplicate, 1ml of the diluted cells into a 60 mm culture dish, label dish

-Use dilutions-10⁻³ and 10⁻⁴

-Repeat for next 14 ml tube

7. Incubate the dishes for 7 days in a 37°C, 5% CO₂ incubator

8. Following staining and counting procedure as described above.

Mediums:

1. MEMA: MEM, 10% heat inactivated FCS, L-Glutamine, pcn/strep
2. MEMB: MEM(Ca free), 10% heat inactivated FCS, L-Glutamine, pcn/strep
 - prevents cells from sticking together
 - for suspension
3. Lindane: a chemical that inhibits formation of gap junctions
4. Triated thymidine: Has short range of ^3H β -particles
 - allows only self-irradiation of labeled cells and no cross-irradiation of unlabeled cells
5. DMSO Dimethyl sulfoxide: "eats-up" free radicals, scavenger
6. Wash MEMA: MEM, 10% calf serum, L-Glutamine, pcn/strep
7. D-medium=5%FBS, 50ug/ml gentamicin

New Cell Lines:

Cell	Gap Junction	HGPRT
WB-r	+	-
WB-s	+	+
WB-aB1	-	-

Split cells:

1. check flask for confluency, 70-80%
2. obtain a waste container
3. remove medium from flask and dispense into waste container
4. wash flask with 10 ml of PBS and remove PBS into waste
5. wash flask with 2 ml of trypsin (use 3ml of trypsin for WB cell line) and put into incubator for 3-6 minutes depending on cell line
6. take out of incubator and "smack" flask to loosen cells and then check under scope to make sure majority of cells are off bottom
7. obtain new culture flasks and label with date, cell type, and initials
8. take 25 ml pipet and fill new flasks with 20 ml of new media
9. add 10 ml of media to flask with cell to neutralize the trypsin and note total vol.=
10. put correct amount of cells into new flasks and leave a small amount in the old flask and add 20 ml of media.
11. incubate

Making Media: MEM_

1. obtain a pouch of media which is good to make 1 L
2. fill round bottom flask with about 300ml of sterile H₂O
3. add pkt of MEM powder
4. wash pkt out with some sterile water to remove all powder
5. then shake flask
6. move to hood and get 25ml pipet and add 30 ml of sodium bicarb (30ml/L)
7. then take flask back to sterile water and fill to 1 liter mark
8. move back to hood and do filtration
9. put filter on new bottle, pour MEM through the filter and vacume will pull the media through.
never let filter get dry turn off vacume then switch to next bottle
10. add L-glut. PCN/strep, and 10% FCS day of use
11. Label bottle MEM_, date, and whatever has been added

L-Glut., pcn/strep, FCS10%:

1. use 17X100 mm snap cap tubes
2. put 10 ml of pcn/strep into 1 tube and 10ml of L-Glut. into another tube
1tube of L-Glut. and 1tube of pcn/strep will go into 1 liter of media

FCS10%= 100ml in 1liter

Freezing cells:

1. prepare freezing mixture
2. remove media from flasks to be frozen
3. wash cells with PBS-5ml and remove PBS
4. wash cells with trypsin-2ml and place in incubator for 3-6 min. depending on cell type
5. "smack" side of flask to loosen cells
6. put 5 ml of media into each flask to neutralize trypsin and then reseed cells into a 50 ml tube and centrifuge at 2000 rpm, 7 min
7. decant the supernatant, break-up pellet and resuspend in correct amount of freezing medium, depends on how many tubes freezing
8. place in cryofreezer apparatus.

for WB cells freezing medium: 50ml

-DMSO 10%=5 ml

-FBS 15%= (10% already in d-media) 5ml

-d-medium=40ml

D-MEDIUM

IN 300 ML STERILE, DISTILLED WATER

1. D-powder=8.97 g (GIBCO ^{Formula} Cat # 78-5470 EF)
2. Na pyruvate(piruv)=10 ml of 0.1M
3. NaCl= 835 g
4. Glucose= 1g

MIX, place magnetic bar in beaker and place on stirring machine

-prepare pH meter by placing in buffer solution that has a pH close to what is eventually desired and turn machine from standby to pH.

-then take pH of the buffer

-once reading is made turn back to standby and rinse probe with dH₂O

-test another pH buffer just for accuracy

-then test sample

-leave probe in beaker for next step

5. add NaOH 1M (~ 0.9ml) to get to pH 6.5, while stirring

6. add 1 g of NaHCO₃ (sod. Bicarb) ***add slowly or will get ppt**leave stirring rod in place

-get a final pH of ~7.17, leave probe in beaker

-add water up to 1 L

7. then filtrate solution

8. then add FBS to a concentration of 5% so, 50ml FBS in 1000ml media or 25 ml FBS in 500ml media

9. gentamicin 25ug/ml add: 2.5 ml to 1000ml of media

Freezing mixture for V79: 10%DMSO, 15%FBS

-for 100 ml

1. filter 10 ml of DMSO

-use 10 cc syringe and fill barrel with DMSO

-then remove needle and replace with sterile filter

-push down plunger on syringe and dispense DMSO into a sterile test tube

2. if using MEMA which already has 10% FBS then just add 5 ml to bring the % up to 15

-if using new MEMA use 15 ml of FBS

*IN ADDITION, if using New MEMA also need to add pcn/strep, and l-glut.

HARVARD UNIVERSITY
SCHOOL OF PUBLIC HEALTH

LABORATORY OF RADIOBIOLOGY
(617) 432-1184

665 HUNTINGTON AVENUE
Boston, Massachusetts 02115

FACSIMILE TRANSMISSION

FAX NUMBER: (617) 432-0107

ATTENTION: Edouard Azarov
DEPARTMENT: Radiology - Radiation Research
FAX NO: 973 - 972 6474
Telephone: 973 - 972 5323
Total #of sheets: 01 (including cover sheet)

FROM: Sonia de Toledo
DEPARTMENT: C. C. B.
TELEPHONE: 617 - 432 1182
Faxed by: Sonia

COMMENTS:

Habit - Good Morning, the car is in for repair today. I may call you later about the car.

D-mediums 1L

D-powder 8.97g

Na piruv. 10 mL of 0.1M

NaCl 0.835g

glucose 1g

add NaOH 1M (≈ 0.9 mL) to pH 6.5

1g NaHCO_3

final pH ≈ 7.17

Sonia

D-Mem -

- 2.5 mL glut
- 5% fbs = 50 mL in 1000 mL

(5 ml / 1000

2000:1

10 µg.

10 mg / ml.

.5

5 mg / ml.

.025 µg

50

5mg

5 | 2000

5 ng in 1000 ml

5000 µg / 1000 ml

5 µg/ml

~~for~~

Order ID C- 10070

LIFE TECHNOLOGIES.

Custom Order Quote

End User Name

Phone FAX

Ship To Address

Shipping Temp

Carrier

Freight Charges

Customer Formula #/LTI Catalog Number Modification: Special Packaging of

Product Info

Minimum Essential Medium *Cat #*

with/without	Component Description	Concentration
w	1.5X Amino Acids	
w	1.5X Vitamins	
w	1X L-Glutamine	
w	2X Non-Essential Amino Acids	
w	Earle's Salts	
wout	D-Glucose	
wout	Phenol Red	
wout	Sodium Bicarbonate (NaHCO3)	

SKU	Suffix	Qty	Unit	Fill	Pkg size	UOM	Container	QCCODE	QCLeadT	UnitPrice	Storage Temp	SubGrp Descripti
78-5470	EF	10	10		1	L		15AR	0	\$104.80	2 to 8C	

*100% 0.5ml / 10
gentamicin*

10.4g/10

Bridget

X-46966

4-6 weeks

100% 20 bottle

100% 500 ml

100%

154.50

3 090.00

Monday, May 22, 2000

Page 1 of 1



Order ID C- 10072

Custom Order Quote

End User Name

Phone FAX

Ship To Address

Shipping Temp

Carrier

Freight Charges

Customer Formula #/LTI Catalog Number Modification: Special Packaging of

Product Info

<u>with/without</u>	<u>Component Description</u>	<u>Concentration</u>
w	1.5X Amino Acids	
w	1.5X Vitamins	
w	1X L-Glutamine	
w	2X Non-Essential Amino Acids	
w	Earle's Salts	
w	Sodium Bicarbonate (NaHCO ₃)	
wout	D-Glucose	
wout	Phenol Red	

SKU	Suffix	Qty	Unit Fill	Pkg size	UOM	Container	QCCODE	QCLeadT	UnitPrice	Storage Temp	SubGrp Descripti
LIQUID	DJ	(20)	0.5	1	L	P1173		0	\$39.55	2 to 8C	
LIQUID	DJ	8	0.5	1	L	P1173		0	\$61.75	2 to 8C	

2 weeks

Experiment Name: V79 mixing
Experiment Performed By: A.B.

Exp.#: 5/26/00, 1
Date:

1. Set rocker-roller in 37°C incubator with 5% CO₂
-set coulter counter
-set centrifuge to cool
 2. Wash 2- 150 cm² flasks with 5 ml of PBS and then Trypsinize the cells using 2 ml of trypsin
-put flasks into incubator for 3 minutes
 3. Neutralize the trypsin and resuspend the cells in 8 ml of MEMB
3000 cells/ml of MEMB (resuspended)
 4. Repool the cells in a ^{Falcon} 50 ml tube and then centrifuge at 1800 rpm, 4°C, 10 min. ²⁰⁰⁰
 5. Decant supernatant, gently break pellet, and resuspend in 10-20 ml of MEMB
1000 to 2000 cells/ml of MEMB used
 6. Gently syringe cells 5 times with a ^{3cc} or ^{5cc} syringe, 21 gauge needle
 7. Perform cell count using coulter counter **Cell count=**
100 ul in 20 ml isobone
 8. Dilute cells to ~4 X 10⁶ cells/ml in MEMB FINAL Vol. = ^{6.5 ml}
- do cell count again to double check- *Cells # 45 x 10⁶ cells* actual ct!
1,072,000 cells/ml
 9. Transfer 2 ml of cell suspension into twenty ^{14 ml} tubes (Falcon plastic test tube-polypropylene) Label each 1-10-each tube will have 2,000,000 cells
25 ul, 10 microliters
 10. Roll tubes in rocker-roller for 3-4 hours at 37°C, 5% CO₂ **Date/Time** *5/26/00; 11-40 a.m.*
start time
 11. After incubation irradiate ¹⁰ of the tubes and keep the other ten in the rocker-roller *end time*
do per tubes
- *keep the tubes to be irradiated cool

Tube #	Total Dose (Rads)	Dose Rate (Rad/Min)	Time (min)	Attenuation
1	0	0	0	0
2	0	0	0	0
3	500	739.0	.67	X-0
4	700	739.0	.95	X-0
5	1000	739.0	1.35	X-0
6	3000	2981.0	1.01	X-0
7	5000	2981.0	1.71	X-0
8	7000	2981.0	2.35	X-0
9	10,000	2981.0	3.35	X-0
10	30,000	2981.0	10.06	X-0

While irradiating keep tubes over ice, *medium* (use at a time, \bar{c})
 "Start" cult in rack (chain) *when you hold tube*

12. When done irradiating wash tubes in 8 ml of cold MEMA wash and centrifuge at 2000 rpm, 4°C, 10 min.

13. Then decant supernatant and break up pellet *as described before*

remove a little from each bottle & place in 4 tubes.
 14. Add the unirradiated cells and media to the irradiated cells, vortex and wash with 5 ml of MEMA wash

15. Centrifuge the cells for 10 min at 2000 rpm 4°C

16. Decant supernatant and break-up pellet and vortex and transfer the cells to a *skinny*, Helina tube (950 μ l) *using 200 μ l pipet tip*

17. Wash test tube with 200 μ l of MEMA and transfer to the skinny tube

18. Centrifuge tubes at 1000 rpm, 4°C for 5 min. - *1st tube in test tubes*

19. Place tubes in a rack and put into refrigerator at 10.5°C for 72 hours

Date/Time: *1/26/00 5:45*

20. After 72 hours, remove tubes 1-10, carefully remove the supernatant using a long-necked pasteur pipet and resuspend pellet in 10 ml of MEMA in a 14 ml tube. Wash the skinny tubes with the MEMA to remove all cells and put into 14 ml tube which corresponds to tubes 1-10.

rewind ←

- put in test tube

remove supernatant from skinny tube Data/Time!
 use to corresponding test tube

21. Centrifuge tubes in a precooled centrifuge for 10 min, 2000 rpm, 4°C

22. Label and prepare tubes for dilution and colony dishes

- load 60 mm petri dish with 4 ml MEMA

16 x 100 mm
 plan to use

- load 2 tubes with 4.5 ml MEMA and label them 1-4, 2-4, 3-4, 4-3, 4-3, 5-3, 5-4

6-3, 6-4, 7-3, 7-4, 8-3, 8-4, 9-3, 9-4, 10-3, 10-4

see tubes on last pg
 10⁻¹, 10⁻², 10⁻³, 10⁻⁴ for each plate

23. Decant supernatant, break-up pellet, and resuspend in 10 ml MEMA wash and centrifuge again

- use 2nd pipet

24. Decant supernatant, break-up pellet, vortex, and resuspend in 2 ml MEMA and then syringe the cells 5 times using a 3 cc syringe, 21 gauge needle. Do Cell Count;

25. Perform four 10-fold dilutions (0.5 ml cells into 4.5 ml MEMA) to get 10⁻¹-10⁻⁴

if will introduce when tubes & vortex, do for rest of tubes
 - use 50 mL pipet - methane

26. Seed 1 ml in triplicate into the 60 mm culture dish

see tubes on last page

27. Then incubate for 7 days in the CO₂ incubator at 37°C

28. Remove dishes from incubator. Remove medium

- label bottom of dish - selected for count - 50-250 colonies

29. Wash each dish 3 times with normal saline and 2 times with methanol

30. Stain with 0.1% crystal violet and count colonies

fluorescent labeled

1-4
 2-4

10⁻¹ 10⁻² 10⁻³ 10⁻⁴
 1-3 1-4

- 1.2
- 2.2
- 3.2, 3.3
- 4.2, 4.3, 4.4
- 5.2, 5.3, 5.4
- 6.2, 6.3, 6.4
- 7.2, 7.3, 7.4
- 8.2, 8.3, 8.4
- 9.3, 9.4, 10.5
- 10.3, 10.4, 10.5

4-2 4-3 4-4

Collect 100ml T done.

cell count
dwell x 4000x
cell growth
00005

① 680, 236, 700, 665

② 690, 672, 693

③ ~~668~~, 715, 666, ~~713~~

④ 670, ~~713~~, 701, 679

⑤ 789, 765, 761, 808 →

index
811 ← 847

⑥ flush → 904, 843, 851

⑦ ~~807~~ 671, 663, 719

⑧ ~~804~~ 774, 736, 758

⑨ ~~803~~, ~~807~~, 732, 748, 711

⑩ ~~807~~, ~~801~~, 804, ~~808~~, 759, ~~813~~ 758
804, 759, 758

V79 COLONY FORMING ASSAY

Experiment Name : ¹³⁷Cs toxicity (acute, cluster, suspension);

Exp. # : ;

Experiment performed by:

Date:

Handwritten notes:
 - 1000 rpm
 - 10 min
 - 1000 rpm
 - 10 min
 - 1000 rpm
 - 10 min

1. Set the rocker-roller at 37°C incubator with 5% CO₂, set the Coulter Counter, wash cells (from two 150 cm² flask, subcultured 1:2, 24h before) with 5 mls PBS, 2 mls trypsinize cells, resuspend in 8 ml MEMB for each flask; pool, vortex, pass five times through 3 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml isotone (Coulter balanced electrolyte solution) *Cell Count = 1449 @ 1374 (3) 1431
2. Dilute to ~4,000,000 cells/ml in MEMB (final volume 11 ml) [Actual count cells/ml] Final vol. = 4.5 ml MEMB + 7.9 stock cells.
3. Transfer 0.5 ml of cell suspension into twenty 14 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall - each tube polypropylene = plastic
4. Roll the tubes for 3-4 h at 37°C, 5% CO₂ UM have 2000 cells 12 ml. Date/Time: 11:40 - 3:30
5. After ~3-4 h incubation period, remove tubes add 8 ml wash MEMA, vortex and centrifuge at 2000 rpm at 4°C for 10 min (precooled centrifuge). Date/Time: 11:40 - 3:30
6. Decant supernatant, click tubes, vortex, resuspend in 3 ml wash MEMA cold - keep cool
7. Centrifuge tubes for 10 min at 2000 rpm, 4°C - then centrifuge
8. Decant supernatant, click tubes, resuspend in 200 ul ice-cold MEMA, - Irradiate 10 of the tubes with the proper dosages - vortex remaining ten tubes of "stock cells" to prevent clumping
9. After all of the tubes have been irradiated, wash the tubes twice with 2 ml MEMA wash. **after first wash add the 200 ul of unirradiated cells to each of the ten 17 X 100mm tubes. 1 2 3 4 - 1 2 3 4 -
10. After last decant, vortex and one tube at a time using a 100 ul sterile micropipet tip - transfer the cells into a "skinny" (Helina) tube - put some mema - wash the skinny tube with 100ul of MEMA - via Test tube - as "wash" -
- *Centrifuge the skinny tubes at 1000 rpm, 4°C, 5min. - wash bottom
11. Transfer tubes at 10°C for 72 h. Date/Time: - 1000 rpm
12. After 72 h, for tubes 1-5, carefully remove the supernatant using long-neck Pasteur pipet, resuspend the pellet in 100 ml

Handwritten notes:
 1000 rpm
 10 min

Handwritten notes:
 1000 rpm
 10 min

Handwritten notes:
 1000 rpm



Counter Counter

- Push power buttons on both mechanics at same time. - Hit COME + LT
- ISO tone II - DC bias-ground check
 - maintain osmolarity
- turn to reset.
- 1. 50ml = 5 sec count
- 50ml = 17 sec count
- then turn to count
- after flush system

- ① open fill knobs
- ② open reset
- ③ close rest
- ④ close fill

clock wise

17. Cells

100 uL of cells in 20mls + subculture

- leave on same setting for flush & count
- mix cells

$$5.672 \times 10^6 \text{ cells/mL} = \frac{5.672 \times 10^6 \text{ cells/mL}}{1 \times 10^6 \text{ con}}$$

$2 \times 10^6 \text{ cells} / 2 \text{ mL}$ need } need 45mls of 1×10^6 per ml
 45 mLs - need } 45×10^6 tot. cells needed.

13. The tubes were irradiated using Mark I irradiator (^{137}Cs gamma-ray), two tube (one tube for pellet and one for the suspension) at a time for a single dose-point, while placing onto a Rainin pipet tip box containing ice as per the Table below

Tube #	Total Dose (Rad)	Dose rate (Rad/min)	Time (min)	Attenuat.
1	0	0	0	0
2	0	0	0	0
3	300	97.3	3.08	X-10
4	600	739.8	0.81	X-0
5	1200	739.8	1.62	X-0
6	0	0	0	0
7	0	0	0	0
8	300	97.3	3.08	X-10
9	600	739.8	0.81	X-0
10	1200	739.8	1.62	X-0

14. After irradiation, carefully remove the supernatant from the top for tubes 6-10, resuspend pellet in 200 ul wash MEMA and transfer the content from all tubes to ten 14 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using pasteur pipet
15. Again add 200 ul wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions in 14 ml tubes
16. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)
17. Labeling and preparation of dilution tubes and colony dishes
- load 60 mm petri dishes with 4 ml MEMA
 - load T-tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.
18. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
19. Centrifuge tubes for 10 min at 2000 rpm, 4°C
20. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times

5/26

Initial cell count \rightarrow 1449, 1374, 1431

Avg cell count = 1418

Cell conc. = 1418×4000

= 5,672,000 cells/ml

We need 45 ml of 1000,000 cells/ml

= 45,000,000 cells

Vol of required = $\frac{45,000,000}{5,672,000}$

= 7.9 ml

did 45 ml
61c need
20 tubes of
2 ml = 40 ml
of media.

Take 8 ml of stock

Add 37 ml of MEM

Avg = 268

275

241

288

Cell conc =

1,072,000 cells/ml

freezing tubes = 1.5 ml's cell & media
per tube

~~add~~

Freezing mixture . prepare 50 ml

add = DMSO - 10% - 5 ml
FBS - 15% - ~~10.5 ml~~ + 5.0 ml
d-medium - 40 ml

media 50 ml

10 ml's FBS 10%
45

5 ml L-glut

5 ml strep/pen