

V79 COLONY FORMING ASSAY

Experiment Name : $^{125}\text{IUdR}$ + MEA; **Exp. # :** 1; **Investigator:** A. Bishayee

Date: 04/09/98

1. Set the rocker-roller at 37°C incubator, set the Coulter Counter, wash cells (from 75 cm² flask, subcultured 1:2, 24h before) with PBS, trypsinize cells, resuspend in 7 ml MEMB, pass five times through 3 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml isotone (Coulter balanced electrolyte solution)
2. Dilute to ~4,00,000 cells/ml in MEMB (final volume 11 ml) [Actual count : 4,54,200 cells/ml]
3. Transfer 1 ml of cell suspension into ten 12 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall
4. Roll the tubes for 3-4 h at 37°C, 5% CO₂ **Date/Time:** 04/09/98 ; 3-15 p.m.
5. Prepare MEMB containing radioactivity in hood
2.6 µl $^{125}\text{IUdR}$ (prepared on 04/09/98) + 5 ml MEMB
6. After 3-4 h, remove test tubes from roller and add MEMB with or without radioactivity according to Table below. **Date/Time:** 4/9/98 ; 7-30 p.m.

Tube #	$^{125}\text{IUdR}$ uCi/ml	Cells in MEMB (ml)	MEMB (ml)	MEMB+ $^{125}\text{IUdR}$ (ml) [1.0 uCi/ml]	MEMA+ MEA 100 ug/ml (1.3 mM)	MEMA	
1	0	1.0	1.0	0	2.0	0	
2	0	1.0	1.0	0	2.0	0	
3	0.01	1.0	0.98	0.02	2.0	0	
4	0.02	1.0	0.96	0.04	2.0	0	
5	0.03	1.0	0.94	0.06	2.0	0	
6	0	1.0	1.0	0	0	2.0	
7	0	1.0	1.0	0	0	2.0	
8	0.01	1.0	0.98	0.02	0	2.0	
9	0.02	1.0	0.96	0.04	0	2.0	
10	0.03	1.0	0.94	0.06	0	2.0	

7. Return test tubes to roller for 12 h, increase the elevation angle of the roller.

Date/Time: 04/09/98 ; 7-45 p.m.

8. While test tubes are rolling label 40 (4x10) gamma-tubes (13 X 100 mm VWR glass test tube)
9. After ~12 h incubation period, remove tubes and centrifuge at 2000 rpm at 4°C for 10 min
(precooled centrifuge) - Date/Time: 04/10/98; 9-00 a.m.
10. During centrifugation, move roller to 10°C and obtain ice
11. Prepare 11 ml of sterile MEA (100 ug/ml) in MEMA, put on ice
12. Remove buckets from centrifuge and carefully remove 100 µl of supernatant and place in pre-labeled gamma-tube.
13. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
14. Centrifuge tubes for 10 min at 2000 rpm, 4°C
15. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
16. Centrifuge tubes for 10 min at 2000 rpm, 4°C
17. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
18. Centrifuge tubes for 10 min at 2000 rpm, 4°C
19. Decant supernatant, click tubes, vortex, resuspend in 2 ml ice cold MEMA containing 0 or 100 ug/ml of MEA as per Table. Keep on ice!
20. Transfer tubes to roller at 10°C for 72 h. Date/Time: 04/10/98; 11-00 a.m.
21. Transfer 10 ul supernatant in three sets of tubes containing small pieces of tissue paper from 100 ul supernatant removed earlier and count them for radioactivity
Date/Time: 04/14/98; 2-00 p.m.
21. After 72 h, remove tubes and place on ice, add 8 ml ice cold wash MEMA.
Date/Time: 04/13/98; 10-00 a.m.
22. Centrifuge tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)
23. Labeling and preparation of dilution tubes and colony dishes
 - load 57 60 mm petri dishes with 4 ml MEMA
 - load 30 T-tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 2.2, 2.3, 2.4, X.2, X.3, X.4, etc.
24. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
25. Centrifuge tubes for 10 min at 2000 rpm, 4°C
26. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
27. Centrifuge tubes for 10 min at 2000 rpm, 4°C
28. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 3 cc syringe with 21 gauge needle
29. Determine cell concentration by transferring 100 µl to Coulter cup
30. Vortex tube, transfer 0.5 ml into dilution tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3, vortex tube X.3 and transfer 0.5 ml to tube X.2 and vortex. Keep tubes on ice.
31. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.

32. Transfer 300 μ l of cell suspension (in triplicate) to gamma tubes for each tube

33. Incubate petridishes for 1 week

34. Count gamma tubes for radioactivity

Date/Time : 04/14/98 ; 3-00 p.m.

35. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol.

Stain colonies with crystal violet

36. Count colonies (50 or more cells). There must be between 25 and 250 colonies for the dish to be a valid data point.

Expt # 1

04/09/98

$$\begin{aligned}\text{Initial Cell Count} &= 10,803; 10,908; 11,072 \\ \text{Avg. Cell Count} &= 10,927 \\ \text{Cell Conc.} &= 10,927 \cdot 6 \times 400 = 4,371,066\end{aligned}$$

For dilution

$$\text{Vol. of Cell suspension required} = \frac{44,000,000}{4,371,066} = 1.00$$

Take 1 ml of cells + 10 ml of MEMB = 11 ml

After dilution,

$$\begin{aligned}\text{Final Count} &= 1200, 1154, 1057, 1131 \\ \text{Avg Cell Count} &= 1135.5 \\ \text{Cell Conc.} &= 1135.5 \times 400 \\ &= 4,54,200 \text{ Cells/ml}\end{aligned}$$

Expt #1

04/08/98

Prepare 5 ml of 1 $\mu\text{Ci}/\text{ml}$ = 5 μCi required.

Stock 1.93 $\mu\text{Ci}/\text{ml}$.

$$\text{Stock required} = \frac{5}{1.93} = 2.59 = 2.6 \mu\text{l}$$

① Take 2.6 μl ^{125}I IUPA

② Add 5 ml of MEMB

Expt # 1

		I-125 window	Open window		
1	1.00	24	244	} 1M	
2	1.00	22	248		
3	1.00	23	289		
4	1.00	23	253	} 2M	
5	1.00	19	240		
6	1.00	27	255		
7	1.00	157	395	} 3M	
8	1.00	168	414		
9	1.00	176	432		
10	1.00	351	586	} 4M	
11	1.00	299	509		
12	1.00	334	586	} 5M	
13	1.00	426	653		
14	1.00	468	705		
15	1.00	468	709	} 6M	
16	1.00	19	230		
17	1.00	21	261		
18	1.00	24	240	} 7M	
19	1.00	20	266		
20	1.00	33	291		
21	1.00	19	224	} 8M	
22	1.00	189	425		
23	1.00	150	389		
24	1.00	178	447	} 9M	
25	1.00	316	544		
26	1.00	335	583		
27	1.00	321	561	} 10M	
28	1.00	448	716		
29	1.00	470	724		
30	1.00	455	482		
<hr/>					
31	1.00	1825	2478	} 9M	
52	1.00	1770	2069		
53	1.00	1642	1940	} 10M	
54	1.00	2095	2397		
55	1.00	2224	2478		
56	1.00	2048	2345		
57	1.00	18	262	} 1C	
58	1.00	24	252		
59	1.00	28	262	} 2C	
60	1.00	33	243		
61	1.00	25	271		
62	1.00	28	272	} 3C	
63	1.00	954	1216		
64	1.00	974	1225		
65	1.00	973	1240	} 4C	
66	1.00	1445	1710		
67	1.00	1582	1842		
68	1.00	1498	1751	} 5C	
69	1.00	2005	2260		
70	1.00	2024	2311		
71	1.00	2152	2426		
72	1.00	27	287	} 6C	
73	1.00	17	247		
74	1.00	21	269	} 7C	
75	1.00	18	248		
76	1.00	28	287		
77	1.00	22	246	} 8C	
78	1.00	1423	1707		
79	1.00	1388	1623		
80	1.00	1392	1668	} 9C	
81	1.00	2255	2539		
82	1.00	2374	2665		
83	1.00	2377	2659	} 10	
84	1.00	2357	2672		
85	1.00	2421	2728		
86	1.00	2360	2557		

Cysteamine
Self-Dose
Experiment
10ML
onto tissues

Cysteamine
Self Dose
Experiment
300ML
washed cells

TABLE-1

Expt. # : A

Date/Time : 04/14/98 ; 2-00 p.m.

Tube #	Medium count for 10 ul (cpm)*	Avg. cpm	dpm $\frac{0.48 \times 1.47}{[cpm/0.142]}$ $cpm/0.7056$	$\mu Ci/ml (A_0)$ on counting [dpm/22200]	$\mu Ci/ml (A_0)$ on addition [A ₀ e ^{-λt}]
1	1, -1, 0	0	0	0	0
2	0, -4, 4	0	0	0	0
3	134, 145, 153	144	204.08	0.00919	0.0097
4	328, 276, 311	305	432.25	0.0194	0.0205
5	403, 445, 445	431	610.82	0.0275	0.0290
6	-4, -2, 1	0	0	0	0
7	-3, 10, -4	0	0	0	0
8	166, 127, 155	149.3	211.59	0.0095	0.0100
9	293, 312, 298	301	426.58	0.0192	0.0203
10	425, 447, 432	434.6	615.9	0.0277	0.0293

04/09/98 ; 7-45

* I-125 window - Empty tube

$$4 \times 24 + 12 + 6.25$$

* Efficiency = 0.48 } I-125
Yield = 1.47 }

= 114.24 h

$$e^{-\lambda t}$$

$$\bar{C} = \frac{0.693 \times 114.24}{1440}$$

$$= 0.9465$$

TABLE-2

Expt. # : 1

Date/Time : 04/14/98 ; 3-00 p.m

Tube #	Radioactivity for 300 ul cell suspension (cpm)	Avg. cpm	dpm $\frac{0.7056}{[cpm/0.142]}$	μ Ci/ml (A) on counting [dpm/666000]	μ Ci/ml (A ₀) after 12 h incubation [A _t e ^{-λt}]
1	-5, 1, 5	0	0	0	0
2	10, 2, 5	0	0	0	0
3	931, 951, 950	944	1337.8	0.0020	0.0021
4	1422, 1559, 1475	1485.3	2105.0	0.0031	0.0033
5	1982, 2001, 2129	2037.3	2887.3	0.0043	0.0045
6	4, -6, -2	0	0	0	0
7	-5, 5, -1	0	0	0	0
8	1400, 1365, 1369	1378	1952.9	0.0029	0.0030
9	2232, 2381, 2354	2312.3	3277.1	0.0049	0.0051
10	2334, 2398, 2287	2323	3292.23	0.0049	0.0051

04/10/98 ; 9-00 a.m.

$$e^{-\lambda t}$$

$$= e^{-\frac{0.693 \times 101}{1440}}$$

$$= e^{-0.0486}$$

$$= 0.9525$$

$$4 \times 24 + 52$$

$$= 101$$

TABLE-3

Expt. # : 1

Date/Time : 04/13/98 ; 11-30 a.m.

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 400]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]
1	951, 893, 903, 863	886.3	354533.3	0
2	416, 402, 478	432	172800	0
3	767, 771, 802	780	312000	0.00673
4	845, 802, 793	813.3	325333.3	0.01383
5	668, 575, 578	607	242800	0.01853
6	874, 876, 858	869.3	347733.3	0
7	874 , —	—	—	—
8	874, 853, 807	844.6	337866.6	0.00887
9	825, 767, 777	789.6	315866.6	0.01614
10	420, 425, 413	421.3	168533.3	0.03026

TABLE-4

Expt # : 1

Date : 04/20/98

Tube.dilution	Colony 1	Colony 2	Colony 3	Avg Colony for X.2	SF
1.2	80	68	60	69.33	
2.2	11	14	14	13*	
3.2	47	49	35	43.66	0.6298
4.3	82	90	78	83.33	0.1201
5.3	61	40	52	51	0.0735
6.2	72	82	95	83	
7.2**	—	—	—	—	—
8.2	20	28	32	26.66	0.3212
9.3	59	53	54	55.3	0.0666
10.4	20	27	29	0.25	0.0030


* not considered

** error during reading
not included

Exp # 1

①—① With MEA
②—② without MEA

DMF = 1.45



1.0

0.1

46 6010

SF

0.00

SEMILOGARITHMIC 4 CYCLES X 70 DIVISIONS
KEUFFEL & ESSER CO. MADE IN U.S.A.

10
9
8
7
6
5
4
3
2
1
9
8
7
6
5
4
3
2
1
9
8
7
6
5
4
3
2
1

0 0.01 0.02 0.03
pci/cell