

14

### COLONY FORMING & HPRT ASSAY FOR V79 CELLS

Experiment Name : 100 % cluster, 3HTdR);  
Investigator: M. Lenarczyk  
21, 2001

Exp.  
Date: May

1. Set the rocker-roller at 37°C incubator with 5% CO<sub>2</sub>, set the Coulter Counter, wash cells (from two 150 cm<sup>2</sup> flasks, sub-cultured 20x10(6) cells/flask 24h before) with PBS, trypsinize cells, each re-suspend in 9 ml MEMB, pass 5x through 10 cc syringe with 21 gauge needle, count the cells using Couter counter (100 ul cell suspension / 20 ml Isotone
  2. Dilute to ~4,000,000 cells/ml in MEMB [Actual count : **3 995 333** cells/ml]
  3. Transfer 1 ml of cell suspension into ten 12 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall
  4. Keep the tubes in the roller for 3-4 h at 37°C, 5% CO<sub>2</sub>
- Date/Time: May 21, 2001/ 18:00

\* - Cell count - 
$$\left[ \frac{(1996 + 2007 + 1990 + 1997)}{4} \right] \times \text{Mode} - 100 \mu\text{l} = 3\ 995\ 333$$

Final count  
= 1996, 2007, 1990  
1997

3,995,333 cells/ml

MS = 100 μl

2982, 2938, 2992  
2970

Cell conc = 5,940,000  
cells/ml

12 ml x 4,000,000  
cells/ml

= 48,000,000

Vol. =  $\frac{48,000,000}{5,940,000}$

= 8 ml + 4 μl

5. Prepare MEMB containing radioactivity in hood

60  $\mu$ l  $^3$ HTdR (Stock : 1  $\mu$ Ci/ $\mu$ l on 4/19/01) + 5 ml MEMB

~~Lot #~~ - NET - 0272 Lot # - 3106-421-

6. After 3-4 h, remove test tubes from roller and add MEMB with or without radioactivity according to Table below.

Date/Time: May 21, 2001 / 20:00

Tube #	$^3$ HTdR uCi/ml	Cells in MEMB (ml)	MEMB (ml)	MEMB+ $^3$ HTdR (ml) [12 uCi/ml]
1	0	1.0	1.0	0
2	0	1.0	1.0	0
3	1	1.0	0.830	0.170
4	1.5	1.0	0.750	0.250
5	2.0	1.0	0.670	0.330
6	2.5	1.0	0.580	0.420
7	3.0	1.0	0.500	0.500
8	4.0	1.0	0.330	0.670
9	5.0	1.0	0.170	0.830
10	6.0	1.0	0	1

7. Return test tubes to roller for 12 h  
20:15

Date/Time: May 21, 2001 /

8. Next day, while test tubes are in roller label 10 gamma-tubes (13 X 100 mm VWR glass test tube)
9. After ~12 h incubation period, remove tubes and centrifuge at 2000 rpm at 4°C for 10 min (precooled centrifuge).

**Date/Time: May 22, 2001 / 11:00**

10. Remove buckets from centrifuge and carefully remove 150 µl of supernatant and place in pre-labeled gamma-tube.
11. Decant supernatant, click tubes, vortex, re-suspend in 10 ml **wash** MEMA
12. Centrifuge tubes for 10 min at 2000 rpm, 4°C
13. Decant supernatant, click tubes, vortex, re-suspend in 10 ml **wash** MEMA
14. Centrifuge tubes for 10 min at 2000 rpm, 4°C regular
15. Decant supernatant, click tubes, vortex, re-suspend in 10 ml **wash** MEMA
16. Centrifuge tubes for 10 min at 2000 rpm, 4°C
17. ~~Decant supernatant, click tubes, vortex, re-suspend in 2 ml of standard culture MEMA~~
19. Decant supernatant, click tubes, vortex, and using 200 µl tips transfer the cells in polypropylene micro-centrifuge tubes (Helena tube, 400 µl)
14. Again add 200 µl ice cold standard culture MEMA, pipet 1-2 times @ transfer the remaining cell in the same polypropylene micro-centrifuge tubes (Total volume ~400 ul)
15. Centrifuge tubes for 5 min at 1000 rpm, 4°C
16. Transfer Hellena tubes at 10.5°C for 72 h

Decant supernatant, click tubes, vortex, transfer the cell suspension in polypropylene micro-centrifuge tubes with attached caps (Helena Plastics, 400 ul) using 200 ul pipet tips

17. Transfer tubes at 10.5°C for 72 h. 05/22/01 **Date/Time: 13:00**
23. Transfer 30 ul supernatant in three sets of 20 ml scintillation vials containing 6 ml liquid scintillation cocktail (EcoLume) from 150 ul supernatant removed earlier (Step 10) and count them for radioactivity **Date/Time:**

24. After 72 h, carefully remove the supernatant from the top, vortex pellet and transfer the cells suspension to ten 12 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using Pasteur pipet

**Date/Time: 05/25/01 13:30**

25. Again add 200 ul wash MEMA in micro-centrifuge tubes, re-suspend and transfer the cell suspensions in 12 ml tubes
26. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)
27. Labeling and preparation of dilution tubes and colony dishes
- load 66, 60 mm petri dishes with 4 ml MEMA
  - load 40 sterile tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3,

HTdR, 100% cluster, V79

05/22/01 / 13:00 72h → 05/25

2.4, 2.5; X.2, X.3, X.4, X.5 etc.

28. Decant supernatant, click tubes, vortex, re-suspend in 10 ml **wash MEMA**
29. Centrifuge tubes for 10 min at 2000 rpm, 4°C
30. Decant supernatant, click tubes, vortex, re-suspend in 10 ml **wash MEMA**
31. Centrifuge tubes for 10 min at 2000 rpm, 4°C
32. Decant supernatant, click tubes, vortex, re-suspend in 10 ml **wash MEMA**
31. Centrifuge tubes for 10 min at 2000 rpm, 4°C
32. Decant supernatant, click tubes, vortex, re-suspend in 2 ml **standard culture MEMA**, pass 5x through 3-5 cc syringe with 21 gauge needle
33. Determine cell concentration by Coulter counting (100 µl cell suspension + 20 ml Isotone II)
34. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5, transfer 0.5 ml into dilution tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3, vortex tube X.3 and transfer 0.5 ml to tube X.2 and vortex. Keep tubes on ice.
35. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.
36. Transfer 100 µl of cell suspension (in triplicate) to 20 ml scintillation vial containing 6 ml cocktail (EcoLume)
37. Incubate P60's for 1 week
38. Count vials for radioactivity **Date/Time :**
39. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol. Stain colonies with 0.05% crystal violet
40. Count colonies. There must be between 25 and 250 colonies for the dish to be a valid data point.

Exp 05/21/01  
 Y79, 100% cluster, HDR

Medicum count (1st)

PAGE: 1

USER: 6 ID:H3 HOWELL PRESET TIME: 1.00 THU 24 MAY 2001 14:45  
 SAMPLE REPEAT: 1 CYCLE REPEAT: 1 SCR:N RS232:N  
 H#: 1 ABC:N QCF:N RCM:N  
 CHANNEL 1-LL: 0 UL: 400 ZSIGMA: 2.00 BKG SUB: 0.00 BKG ZSIG: 0.00 LSR: 0  
 DATA CALC: CPM, UNKNOWN REPLICATES: 1 NORM FACTOR: 0 1.00000  
 HALF LIFE(DAYS):N

SAM	POS	CH	CPM	ZSIG%	TIME	EL TIME	AVG H#	ERR
1	29- 1	1	14.00	53.45	1.00	1.42	76.0	
2	29- 2	1	5.00	89.44	1.00	3.00	77.0	
3	29- 3	1	11.00	60.30	1.00	4.63	77.0	
4	29- 4	1	10.00	63.25	1.00	6.25	77.0	
5	29- 5	1	10.00	63.25	1.00	7.82	76.0	
6	29- 6	1	12.00	57.74	1.00	9.38	77.0	
7	29- 7	1	22800.00	1.97	0.45	10.40	76.0	
8	29- 8	1	22504.44	1.99	0.45	11.41	76.0	
9	29- 9					11.88		101
10	29-10	1	32180.00	1.88	0.35	12.77	77.0	
11	29-11	1	33826.66	1.99	0.30	13.63	77.0	
12	29-12	1	32825.71	1.87	0.35	14.53	75.0	
13	29-13	1	43968.00	1.91	0.25	15.33	77.0	
14	29-14	1	39446.66	1.84	0.30	16.24	78.0	
15	29-15	1	44848.00	1.89	0.25	17.04	79.0	
16	29-16	1	55820.00	1.89	0.20	17.79	77.0	
17	29-17	1	55850.00	1.89	0.20	18.54	78.0	
18	29-18	1	55470.00	1.90	0.20	19.35	76.0	
19	**- 1	1	70800.00	1.94	0.15	20.10	78.0	
20	**- 2	1	72926.66	1.91	0.15	20.80	77.0	
21	**- 3	1	72206.66	1.92	0.15	21.50	78.0	
22	**- 4	1	96320.00	1.66	0.15	22.26	77.0	
23	**- 5	1	95260.00	1.67	0.15	22.96	78.0	
24	**- 6	1	97546.66	1.65	0.15	23.66	79.0	
25	**- 7	1	124119.99	1.47	0.15	24.36	80.0	
26	**- 8	1	118100.00	1.84	0.10	25.07	79.0	
27	**- 9	1	114119.99	1.53	0.15	25.77	77.0	
28	**-10	1	146066.66	1.35	0.15	26.48	79.0	
29	**-11	1	135526.66	1.40	0.15	27.19	79.0	
30	**-12	1	136560.00	1.71	0.10	27.90	79.0	
31	**-13	1	29508.57	1.97	0.35	28.85	0.0	

Standard  
 DPM-9820

Exp:  
05/21/01

V79, 100% cluster, HDR

Medium count (2nd)

PAGE: 1

USER: 6 ID:H3 HOWELL PRESET TIME: 1.00 THU 24 MAY 2001 16:51  
SAMPLE REPEAT: 1 CYCLE REPEAT: 1 SCR:N RS232:N  
H#: 1 ABC:N QCF:N RCM:N  
CHANNEL 1-LL: 0 UL: 400 2SIGMA: 2.00 BKG SUB: 0.00 BKG 2SIG: 0.00 LSR: 0  
DATA CALC: CPM, UNKNOWN REPLICATES: 1 NORM FACTOR: 1.00000  
HALF LIFE(DAYS):N

SAM	POS	CH	CPM	2SIG%	TIME	EL TIME	AVG H#	ERR
1	29- 1	1	13.00	55.47	1.00	1.42	76.0	
2	29- 2	1	8.00	70.71	1.00	3.00	75.0	
3	29- 3	1	14.00	53.45	1.00	4.58	77.0	
4	29- 4	1	10.00	63.25	1.00	6.14	77.0	
5	29- 5	1	12.00	57.74	1.00	7.71	77.0	
6	29- 6	1	10.00	63.25	1.00	9.33	77.0	
7	29- 7	1	22695.55	1.98	0.45	10.33	76.0	
8	29- 8	1	22082.00	1.90	0.50	11.38	76.0	
9	29- 9					11.85		101
10	29-10	1	32151.43	1.89	0.35	12.80	78.0	
11	29-11	1	33220.00	1.85	0.35	13.71	77.0	
12	29-12	1	33496.66	2.00	0.30	14.56	77.0	
13	29-13	1	43672.00	1.91	0.25	15.36	77.0	
14	29-14	1	40036.00	2.00	0.25	16.21	78.0	
15	29-15	1	44144.00	1.90	0.25	17.08	80.0	
16	29-16	1	55540.00	1.90	0.20	17.83	77.0	
17	29-17	1	56410.00	1.88	0.20	18.58	77.0	
18	29-18	1	55240.00	1.90	0.20	19.33	76.0	
19	**- 1	1	71400.00	1.93	0.15	20.08	77.0	
20	**- 2	1	73706.66	1.90	0.15	20.78	78.0	
21	**- 3	1	72853.33	1.91	0.15	21.48	78.0	
22	**- 4	1	95473.33	1.67	0.15	22.18	77.0	
23	**- 5	1	96426.66	1.66	0.15	22.88	78.0	
24	**- 6	1	97420.00	1.65	0.15	23.58	80.0	
25	**- 7	1	123160.00	1.80	0.10	24.28	79.0	
26	**- 8	1	117779.99	1.50	0.15	24.98	78.0	
27	**- 9	1	115286.66	1.52	0.15	25.69	78.0	
28	**-10	1	146950.00	1.65	0.10	26.40	79.0	
29	**-11	1	138186.66	1.39	0.15	27.10	78.0	
30	**-12	1	134060.00	1.41	0.15	27.81	79.0	
31	**-13	1	29797.14	1.96	0.35	28.78	0.0	

Stavola  
DRA 98200

Plate 200,000 for MN  
1,000,000 for 6 HTRT

MS = 500  $\mu$ l; Background = 14

↓  
DAY 0 of mutant expression

TABLE-3

Expt. #: V79, 100%, HTRT

Date/Time: 05/25/04; 3-00 pm

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 <sup>6</sup> Cells/ml]
1	6134, 6219, 6161	6171	2,468,533	
2	<del>6789, 6628, 6961</del> 7301, 7004, 7688	6792	2,717,066	
3	6934, 6950, 6751	6878	2,751,333	
4	7188, 7166, 7101	7151	2,860,666	
5	7464, 7568, 7366	7466	2,986,400	
6	<del>7036, 6720, 6864, 6800</del>	6794	2,717,866	
7	6437, 6516, 6442	6465	2,586,000	
8	6548, 6789, 6661	6666	2,666,400	
9	6044, 6415, 6176	6211	2,484,666	
10	6515, 6607, 6469	6530	2,612,133	

4/6 →  
4/6 →

Mode - 500  $\mu$ l  
Bcign - see above

For survival:

- 1.2
- 2.2
- 3.2
- 4.2
- 5.2 5.3
- 6.2 6.3
- 7.2 7.3
- 8.2 8.2
- 9.2 9.3
- 10.2 10.3

Cell activity (1st)  
 V79, 100%, HicR  
 May 21, 2001

PAGE: 1

USER: 6 ID:H3 HOWELL PRESET TIME: 1.00 FRI 25 MAY 2001 17:01  
 SAMPLE REPEAT: 1 CYCLE REPEAT: 1 SCR:N RS232:N  
 H#: 1 ABC:N GCF:N RCM:N  
 CHANNEL 1-LL: 0 UL: 400 2SIGMA: 2.00 BKG SUB: 0.00 BKG 2SIG: 0.00 LSR: 0  
 DATA CALC: CPM, UNKNOWN REPLICATES: 1 NORM FACTOR: 1.00000  
 HALF LIFE(DAYS):N

SAM	POS	CH	CPM	2SIG%	TIME	EL TIME	AVG H#	ERR
1	**	1	11.00	60.30	1.00	1.42	88.0	
2	**	2	9.00	66.67	1.00	3.00	86.0	
3	**	3	13.00	55.47	1.00	4.58	87.0	
4	**	4	11.00	60.30	1.00	6.14	89.0	
5	**	5	10.00	63.25	1.00	7.71	89.0	
6	**	6	12.00	57.74	1.00	9.28	90.0	
7	**	7	2450.00	4.04	1.00	10.89	88.0	
8	**	8	2448.00	4.04	1.00	12.46	87.0	
9	**	9	2965.00	3.67	1.00	14.03	89.0	
10	**	10	7151.00	2.37	1.00	15.61	88.0	
11	**	11	7461.00	2.32	1.00	17.18	90.0	
12	**	12	7668.00	2.28	1.00	18.75	90.0	
13	**	13	8304.00	2.52	1.00	20.33	91.0	
14	**	14	4966.00	2.84	1.00	21.90	88.0	
15	**	15	4902.00	2.86	1.00	23.53	89.0	
16	**	16	4678.00	2.92	1.00	25.10	90.0	
17	**	17	3811.00	3.24	1.00	26.67	89.0	
18	**	18	4373.00	3.02	1.00	28.23	89.0	
19	29	1	9558.00	2.05	1.00	29.87	89.0	
20	29	2	8447.00	2.18	1.00	31.43	88.0	
21	29	3	9503.00	2.05	1.00	33.01	90.0	
22	29	4	8970.00	2.11	1.00	34.57	87.0	
23	29	5	9266.00	2.08	1.00	36.14	87.0	
24	29	6	10171.00	1.98	1.00	37.71	89.0	
25	29	7	13630.67	1.98	0.75	39.07	93.0	
26	29	8	12164.71	1.97	0.85	40.49	88.0	
27	29	9	12307.06	1.96	0.85	41.91	88.0	
28	29	10	12753.75	1.98	0.80	43.27	87.0	
29	29	11	15354.29	1.93	0.70	44.54	88.0	
30	29	12	15958.46	1.96	0.65	45.75	89.0	
31	29	13	29491.43	1.97	0.35	46.72	0.0	



Cell activity (2uol

U79, HialR, 100%

May 21, 2001.

PAGE: 1

USER: 6 ID: H3 HOWELL      PRESET TIME: 1.00      TUE 29 MAY 2001 10:11  
SAMPLE REPEAT: 1 CYCLE REPEAT: 1 SCR: N      R5232: N  
H#: 1 ADC: N QCF: N RCM: N  
CHANNEL 1-LL: 0 UL: 400 2SIGMA: 2.00 BKG SUB: 0.00 BKG 2SIG: 0.00 LSR: 0  
DATA CALC: CPM, UNKNOWN REPLICATES: 1      NORM FACTOR: 1.00000  
HALF LIFE(DAYS): N

SAM	POS	CH	CPM	2SIG%	TIME	EL TIME	AVG H#	ERR
1	**	1	11.00	60.30	1.00	1.42	88.0	
2	**	2	11.00	60.30	1.00	3.00	87.0	
3	**	3	11.00	60.30	1.00	4.58	88.0	
4	**	4	11.00	60.30	1.00	6.14	89.0	
5	**	5	10.00	63.25	1.00	7.71	89.0	
6	**	6	8.00	70.71	1.00	9.28	90.0	
7	**	7	1909.00	4.58	1.00	10.84	89.0	
8	**	8	1942.00	4.54	1.00	12.40	87.0	
9	**	9	2410.00	4.07	1.00	13.97	89.0	
10	**	10	5644.00	2.66	1.00	15.53	87.0	
11	**	11	5737.00	2.64	1.00	17.10	89.0	
12	**	12	5998.00	2.58	1.00	18.67	90.0	
13	**	13	5173.00	2.78	1.00	20.24	90.0	
14	**	14	3979.00	3.17	1.00	21.82	89.0	
15	**	15	3831.00	3.23	1.00	23.40	90.0	
16	**	16	3678.00	3.30	1.00	24.97	90.0	
17	**	17	2968.00	3.67	1.00	26.60	90.0	
18	**	18	3406.00	3.43	1.00	28.17	90.0	
19	29-	1	7419.00	2.32	1.00	29.79	89.0	
20	29-	2	6589.00	2.46	1.00	31.36	88.0	
21	29-	3	7556.00	2.30	1.00	32.92	90.0	
22	29-	4	6643.00	2.45	1.00	34.50	86.0	
23	29-	5	7143.00	2.37	1.00	36.07	88.0	
24	29-	6	7841.00	2.26	1.00	37.64	88.0	
25	29-	7	10642.11	1.99	0.95	39.16	93.0	
26	29-	8	9280.00	2.08	1.00	40.73	88.0	
27	29-	9	9368.00	2.07	1.00	42.32	89.0	
28	29-	10	10225.00	1.98	1.00	43.89	86.0	
29	29-	11	11546.67	1.96	0.90	45.42	89.0	
30	29-	12	12587.50	1.99	0.80	46.78	88.0	
31	29-	13	29540.00	1.97	0.35	47.74	0.0	

Standard  
DPM - 98200

TABLE-4

Expt # V79, 100% HCR

Date :

38. Count vials for radioactivity

Date/Time :

3  
}

Tube.dilution	Colony 1	Colony 2	Colony 3	Avg Colony	SF
1.2	86	78	72		
2.2	50	76	63		
3.2	71	81	68		
4.2	68	64	66		
5.2	62	55	55		
6.2	81	66	67		
7.2	29(?)	38	46		
8.2	59	53	51		
9.2	38	46	38		
10.2	40	36	36		

2  $\mu\text{g}/\mu\text{l}$ .

3  $\mu\text{g}/\text{ml}$ .

CYTOCHALASIN-B

~~10 x 3 ml~~

10 sample x ~~3 ml~~ + 3 ml = 30 ml MEMB

~~10 x~~ 10 sample  $\rightarrow$  30 ml

35 ml.

3  $\mu\text{g}/\text{ml}$

x - ~~35~~ ml  
35

105  $\mu\text{g}$   
||

105  $\mu\text{g}$  - y  $\rightarrow$  52.5  $\mu\text{l}$  of stock CytB  
2  $\mu\text{g}$  -  $\mu\text{l}$  (2  $\mu\text{g}/\mu\text{l}$ )

Take 52.5  $\mu\text{l}$  of 2  $\mu\text{g}/\mu\text{l}$  CytB

Add ~~it~~ into the tube with 35 ml MEMB

Mix  $\rightarrow$  (final conc. of CytB ~~52.5  $\mu\text{g}/35\text{ml}$~~ )

$$= 105 \mu\text{g} / 35 \text{ml} = \underline{\underline{3 \mu\text{g}/\text{ml}}}$$

↑  
Micrococci

DAY 0 = May 25, 2001.

HPRT mutant expression - DAY 3

May 28, 2001

Exp: V79, 100%, HTdK

Sample:	Coulter Count			# of cells per ml. $\times 10^6$	Total # of cells $\times 10^6$	Volume for $10^6$ cells (ml)
1	8643	8671	8699	3.45	20.7	0.29
2	9465	9199	9250	3.7	22.2	0.27
3	7846	7945	7764	3.1	21.9	0.32
4	6119	6208	5965	2.4	14.6	0.41
5	Contamination			←—————		
6	9556	9329	9476	3.8	22.6	0.27
7	9301	9238	9357	3.8	22.2	0.27
8	8744	8582	8382	3.4	20.5	0.29
9	8121	8066	8227	3.2	19.4	0.31
10	7639	7846	7535	3.1	18.3	0.33

Background - 36

Mode - 500  $\mu$ l.

Procedure:

1. Wash 1x with PBS-A
2. Trypsinize (2ml trypsin/T75, 2-3min RT) & harvest the cells (+8ml MEM)
3. Spin the cells down (5 min, 2K rpm, 4°C)
4. Syringe the cells 3x
5. Count the cells (100  $\mu$ l cell susp. + 20 ml Isotone)
6. Plate  $10^6$  cells / T75 flask + 15ml MEM

KIPET mutant expression DAY 6

May 31, 2009

Exp: V79, 100%, HTdR,

Sample	Coulter count		# of cells per ml $\times 10^6$	Total # of cells per in 7ml.	Volume for 10 <sup>6</sup> cells.	
1	6426	<del>5162</del> 6762	6720	2.7	18.5	0.37
2	7890	7715	7858	2.2	15.2	0.46
3	9414	9836	9784	3.9	27.1	0.26
4	6843	6773	6778	2.7	18.9	0.37
5	<hr/>					
6	<del>8571</del>	7052	6875	2.8	19.5	0.36
7	9266	9352	9420	3.7	26.1	0.27
8	<del>8028</del> <del>6825</del>	7731 <del>7652</del>	<u>7731</u>	3.1	21.9	0.32
9	10083	10308	10302	4.1	28.6	0.24
100	9767	9590	9710	3.9	27.1	0.26

Beckpr - 17

Mode - 500  $\mu$ l

HPRT mutant expression, DAY 9

, June 3, 2001

Exp: V79, 100%, HDR,

Sample	Coulter count			# of cells per ml.	Total # of cells in 6 ml.	10 <sup>6</sup> cells (Vol.)
1	9864	9779	10066	4.0	23.7	0.25
2	11462	11509	11161	4.5	27.3	0.22
3	10164	9996	10152	4.0	24.2	0.25
4	9967	10586	10756	4.2	25.0	0.24
5	<hr/>				<hr/>	<hr/>
6	9956	9819	9979	4.0	23.7	0.25
7	9471	9533	9522	3.8	22.7	0.33
8	10066	10401	10343	4.1	24.6	0.24
9	9493	9510	9830	3.8	23.0	0.26
10	11261	11208	11261	4.5	26.9	0.22

Belign - 17

Mode - 500 µl

HPET<sup>-</sup> expression DAY 11  
 Exp: V79, 100%, H1dR,

June 5, 2001

Sample	Couter count			# of cells per ml $\times 10^6$	Total # of cells (in 6 ml)	$10^6$ cells (ml)
• 1	1668	1672	1642	0.66	3.9	1.52
• 2	1076	1131	1098	0.43	2.6	2.32
3	4557	4716	4518	1.83	<del>11.00</del>	0.55
4	3971	4078	3996	1.60	9.6	0.63
5	<hr/>					
6	3919	3886	3878	1.55	9.3	0.65
• 7	1985	1949	1962	0.80	4.7	1.25
8	3083	3019	3061	1.22	7.3	0.82
9	3789	3614	3772	1.48	8.9	0.68
10	4934	4986	4924	1.97	11.8	0.51

BeCupr-16

11.0e - 500µl

NOTE : T75 #1, 2 & 7 was capped too tight!

HPRi - Challenge

June, 2001

Exp: V79, 100%, UTR,

Sample	Colter Count	# cells/ml	Total # (# cells/ml)	$0.2 \times 10^6$
1	5379, 5556 5437	2178533	13.02	92
2	4488, 4359 4494	1774400	10.6	111
3	4429 4526 4466	1785067	10.7	111
4	4349 4142 4271	1697200	10.2	118
? 6	5306 5198 5084	2074000	12.9	95
? 7	5098 5141 5032	2025067	12.2	100
? 8	5509 5395 5390	2168133	13.0	91
9	5239 4778 5110	2012533	12.1	100
10	5110 5136 5012	2030000	12.2	100

Belogr- 11

Mode - 500

STOCK  $200000/5ml = 40000/ml$

$200000/5ml = 40000/ml$

↓  
10x  $0.5ml + 4.5ml \text{ MEMA} = 4000/ml$

$20000/5ml = 4000$

↓  
100x  $0.5ml + 4.5ml \text{ MEMA} = 400/ml$

$20000/5ml = 4000/ml$

$4000/5ml =$

↓  
for PE  $0.5ml / P60' \times 3$



Plating efficiency for HPBI Challenge

13 July, 2001.

Exp. 100%, HADR, V79

1	100	138	140	13	15	17
2	135	128	129			
3	157	180	160	17	21	
4	173	155	180	18	11	17
5	<hr/>					
6	132	150	124	20	22	
7	129	128	119	17	21	24
8	130	137	135	25	17	
9	119	100	110			20
10	115	132	127	20	17	19

HPET - 100% HTDR, 179

Plated for 61, 62 counts on June 7

Fixed & stained on June 18

1 0, 0, 0, 0, 0, 0, 1, C  
2 0, 6, 1, 2, 1, 1, 1, C, C,  
3 2, 1, 4, 4, 4, 1, 1, 3, 4,  
4 0, 0, 1, 4, 2, 3, 3, 1, 3, C,  
5 -----  
6 0, 0, 0, 0, 2, 2, 1, 2, 1, 1-0  
7 0, 7, 7, 4, 7, 10, 5, 8, C, C  
8 0, 0, 0, 2, 0, 0, 1, ①, 1, C  
9 0, 0, 2, 1, 2-4, 2, 1, 1, 6, 3,  
10 0, 2, 3, 1, 3, 3, 3, 3, C, C

**COLONY FORMING ASSAY – V79 cells / 50 % cluster/ 3HTdR**

**Experiment Name :** <sup>3</sup>HTdR toxicity (cluster, 50% labeling);

**Exp. June 07, 2001**

**Experiment performed by:** M. Lenarczyk

**Date:** 06/07/01

1. Set the rocker-roller at 37°C incubator with 5% CO<sub>2</sub>, set the Coulter Counter, wash cells (from 4x 175 cm<sup>2</sup> flask, sub-cultured 1:2, 24h before) with DPBS, trypsinize cells (2 ml/flask), re-suspend cells in 15 ml MEMB/flask, pool 2 flask in 50mm conical tube, spin the cells down (4 o C, 2000 rpm, 10"), aspirate supernatant, and wash the cells in 15 ml DPBS, aspirate supernatant, pass five times through 10 cc syringe with 21 gauge needle, perform cell count by transfer 100 µl in Coulter cup containing 20 ml IsoLume
2. Dilute to ~2,000,000 cells/ml in MEMB <sup>4528, 4603, 4687</sup> [actual count : 1.840 000 cells/ml]
3. Transfer 1 ml of cell suspension into 20 12 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall
4. Keep the tubes in the roller for 3-4 h at 37°C, 5% CO<sub>2</sub> **Date/Time: June., 07, 2001 / 16:00**
5. Prepare MEMB containing radioactivity in hood  
800 µl <sup>3</sup>HTdR (Stock : 1 µCi/µl on ., ., 2001 ) + 4200 ml MEMB
6. After 3-4 h, remove first set of 10 test tubes from roller and add MEMB with or without radioactivity according to Table below. **Date/Time: June., 07, 2001 / 19:00**

Tube #	<sup>3</sup> HTdR µCi/ml	Cells in MEMB (ml)	MEMB (ml)	MEMB+ <sup>3</sup> HTdR Working sol. Activity [160µCi/ml] (ml)
1	0	1.0	1.0	0
2	0	1.0	1.0	0
3	10	1.0	875	0.125
4	20	1.0	0.750	0.250
5	30	1.0	0.625	0.375
6	40	1.0	0.500	0.500
7	50	1.0	0.375	0.625
8	60	1.0	0.250	0.750
9	70	1.0	0.125	0.875
10	80	1.0	0	1

7. Return test tubes to roller for 12 h .

**Date/Time: June, 07, 2001 / 19:30**

1.  $1.8 \times 10^6 / \text{ml} \rightarrow$
2. Roll start  $\rightarrow 16:00 \rightarrow$   $\rightarrow 19:00$  (3 hrs).
3. ~~19:30~~  
Add H<sub>2</sub>O<sub>2</sub>

	↓ 19:30 (July 7, 2001)	→ 10:45 (July 8, 2001)	(15 hrs)	H-3 + MEMB (160 μCi/ml)
		MEMB		
1	0	100	1	0
2	0	100	1	0
3	10		0.875	0.125
4	20		0.750	0.250
5	30		0.650	0.375
6	40		0.500	0.500
7	50		0.375	0.625
8	60		0.250	0.750
9	70		0.125	0.875
10	80		0	1.000
				<u>4.500</u>

3 ml MEMB + 800 μCi

4.2 ml + 800 μCi (160 μCi/ml)

4. Take 30 μl for medium activity counting

5. Wash 3x in Wash MEMA  $\rightarrow$

6. Emerge

6-TO

1.  $2.5 \text{ mg} + 10 \text{ ml NEHA}$   
 $2.5 \text{ mg} + 10 \text{ ml NEHA}$  } dissolve  $\rightarrow 5 \text{ mg} / 20 \text{ ml}$  (STOCK SOL).

2. ~~Add~~ Add STOCK SOL. (sterilize ~~using~~) into  $255 \text{ ml NEHA}$   $\xrightarrow{\text{mix}}$  (WORKING SOL)

3. Take 3 ml of Working sol & add into P100' with 7 ml.  
It gives you final concentration of 6-TO

$$5 \text{ mg} / 275 \text{ ml} = 5000 \mu\text{g} / 275 \text{ ml}$$

$$\times \quad - \quad 3 \text{ ml}$$

$$\downarrow$$
$$54.55 \mu\text{g}$$

$$54.55 \mu\text{g} / 3 \text{ ml} \xrightarrow{\text{add}} + 7 \text{ ml} = 54.55 \mu\text{g} / 10 \text{ ml}$$

$$\downarrow \downarrow$$
$$\boxed{5.45 \mu\text{g} / \text{ml}}$$

July 7, 2001

---

8. Next day, while test tubes are in roller label 10 gamma-tubes (13 X 100 mm VWR glass test tube)

9. After ~12 h incubation period, remove tubes and centrifuge at 2000 rpm at 4°C for 10 min (precooled centrifuge). **Date/Time: June, 07, 2001 / 10:45 ( hrs )**

NOTE: Set of non-labeled cells were not rolled. The roller is broken. !!!!

10. Remove buckets from centrifuge and carefully remove 150 µl of supernatant from labeled tubes and place in pre-labeled gamma-tube.

11. Decant supernatant, click tubes, vortex, re-suspend in 10 ml **wash** MEMA.

Steps 11-16 are common for labeled as well as non-labeled tubes.

12. Centrifuge tubes for 10 min at 2000 rpm, 4°C

13. Decant supernatant, click tubes, vortex, re-suspend in 10 ml **wash** MEMA

14. Centrifuge tubes for 10 min at 2000 rpm, 4°C

15. Decant supernatant, click tubes, vortex, re-suspend in 5 ml of **regular** MEMA

16. Centrifuge tubes for 10 min at 2000 rpm, 4°C