




## V79 HPRT Mutants/Cell Hypoxic versus Aerobic Clusters



## V79 HPRT Mutants/Cell Hypoxic versus Aerobic Clusters



colonies	av col	std col	mutants	av muts	std mut	mut/cell	std mut/cell	Survival	PE	S/S0	StdS/S0									
1	150	166	149	155	9.53939	1	2	1	1	1	0.8	0.632456	5.16129E-06	4.08036E-06	140	166	152	0.7	1	0.0752
2	137	129	152	139.33	11.6762	0	1	0	1	0					135	149	157			
3	121	137	145	134.33	12.2202	12	17	19	15	17	16	2.645751	0.000119107	1.96954E-05	52	63	73	0.3	0.4182	0.0701
4	152	130	119	133.67	16.8028	38	45	29	50	39	40.2	7.918333	0.000300748	5.92394E-05	36	41	31	0.2	0.2403	0.03337
5	117	125	139	127	11.1355	60	62	59	41	58	56	6.514693	0.000440945	6.70448E-05	5.3	6	6.8	0	0.0403	0.00501
6	165	147	155	155.67	9.0185	2	0	0	1	0	0.7	0.823273	4.49679E-06	5.28869E-06	121	115	139	0.7	1	0.09369
7	141	159	139	148.33	11.0151	2	1	1	0	0					137	152	144			
8	129	125	138	130.67	6.65833	11	6	5	7	14	8.6	3.781534	6.58163E-05	2.89403E-05	69	75	82	0.4	0.5594	0.04342
9	147	152	118	139	18.3578	29	31	19	37	20	27.2	7.828892	0.000195683	5.48841E-05	47	56	68	0.3	0.4183	0.06343
10	167	145	149	153.67	11.7189	36	39	47	40	48	42	5.244044	0.000273319	3.41261E-05	10	11	12	0.1	0.0817	0.00867

Date	Cell Doublings	Counts	Vol(ml)			
9/24/99	0	589	598	571	2	4688000
10	0	611	627	631		4984000
3	0	541	559	561		4429333
6	0	629	642	629		5066667
12	0	667	656	672		5320000
0	0	542	581	559		4432000
0	0	620	635	642		5058667
3	0	529	549	557		4360000
6	0	607	598	622		4872000
12	0	511	509	507		4072000

483.5 1.0  
0.91  
1.04  
1.100  
474.5 0.92  
1.02  
0.86

Date	Doublings	Sample	Total Doublings							
9/29/99	411	431	435	6	1E+07	3.35276	1.05	1.00	1	11.12182
7	471	481	459		1.1E+07	3.47612			2	10.91181
day 3	389	362	372		8984000	3.16736	0.90	0.86	3	10.93538
	332	321	341		7952000	2.99132	0.80	0.76	4	10.62649
	441	456	465		1.1E+07	3.44573	1.10	1.05	5	11.37441
	432	444	456		1.1E+07	3.41359			6	11.02331
	409	422	436		1E+07	3.34142	1.05	1.00	7	10.95364
	381	392	401		9392000	3.23143	0.95	0.89	8	10.74782
	356	365	369		8720000	3.12433	0.87	0.83	9	10.93045
	403	372	385		9280000	3.21412	0.93	0.89	10	11.03198

10/1/99	772	761	756	6	1.8E+07	4.19472
	686	655	677		1.6E+07	3.99856
	701	711	722		1.7E+07	4.09356
	656	631	634		1.5E+07	3.84186
	732	745	739		1.8E+07	4.14796
	741	756	762		1.8E+07	4.17568
	635	659	662		1.6E+07	3.96791
	672	657	659		1.6E+07	3.99132
	713	732	742		1.7E+07	4.12695
	699	710	729		1.7E+07	4.09626

10/4/99	499	488	502	6	1.2E+07	3.57434
	436	456	462		1.1E+07	3.43723
	522	532	542		1.3E+07	3.67448
	536	542	539		1.3E+07	3.69332
	561	572	585		1.4E+07	3.78073
	437	452	462		1.1E+07	3.43403
	501	533	529		1.3E+07	3.64432
	490	471	478		1.2E+07	3.52507
	522	535	542		1.3E+07	3.87717
	531	555	563		1.3E+07	3.72159

## V79 COLONY FORMING ASSAY

Experiment Name :  $^{137}\text{Cs}$  toxicity (acute, cluster, suspension);  
Experiment performed by: A. Bishayee

Exp. #: 2;  
Date: 09/20/99

1. Set the rocker-roller at 37°C incubator with 5% CO<sub>2</sub>, set the Coulter Counter, wash cells (from two 150 cm<sup>2</sup> flusk, subcultured 1:2, 24h before) with PBS, trypsinize cells, resuspend in 7 ml MEMB for each flusk, pool, vortex, pass five times through 3 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml isotone (Coulter balanced electrolyte solution)
2. Dilute to ~4,000,000 cells/ml in MEMB (final volume 11 ml) [Actual count : cells/ml]
3. Transfer 1 ml of cell suspension into ten 14 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall
4. Roll the tubes for 16 h at 37°C, 5% CO<sub>2</sub>      Date/Time: 09/20/99 ; 6-00 PM
5. After ~16 h incubation period, remove tubes, add 8 ml wash MEMA, vortex and centrifuge at 2000 rpm at 4°C for 10 min (precooled centrifuge). Date/Time: 09/21/99 ; 9-30 a.m.
6. Decant supernatant, click tubes, vortex, resuspend in 3 ml wash MEMA
7. Centrifuge tubes for 10 min at 2000 rpm, 4°C
8. Decant supernatant, click tubes, resuspend in 200 ul ice cold MEMA, transfer the cell suspension in polypropylene microcentrifuge tubes with attached caps (Helena Plastics, 400 ul) using pipet tips
9. Again add 200 ul ice cold MEMA, resuspend and transfer the cell suspensions in the same polypropylene microcentrifuge tubes (Total volume ~400 ul)
10. Centrifuge tubes for 5 min at 1000 rpm, 4°C
11. Transfer tubes at 10°C for 72 h.      Date/Time: 09/21/99 ; 11-00 a.m.
12. After 72 h, for tubes 1-5, carefully remove the supernatant, resuspend the pellet in 400 ml MEMA and place all tubes on the perforated plate of Rainin pipet tip box containing ice (to maintain ~10.5°C)

~~Done~~

13. The tubes were irradiated using Mark I irradiator (<sup>137</sup>Cs gamma-ray), two tube (one tube for pellet and one for the suspension) at a time for a single dose-point, while placing onto a Rainin pipet tip box containing ice as per the Table below

Date/Time: 09/24/99;

0020156

Tube #	Total Dose (R)	Dose rate (Rad/min)	Time (min)	Attenuat.
1	0	0	0	0
2	0	0	0	0
3	300	97.3	3.08	X-10
4	600	<del>387.59.8</del>	0.81155	X-0
5	1200	<del>387.59.8</del>	1.62309	X-0
6	0	0	0	0
7	0	0	0	0
8	300	97.3	3.08	X-10
9	600	<del>387.59.8</del>	0.81155	X-0
10	1200	<del>387.59.8</del>	1.62309	X-0

14. After irradiation, carefully remove the supernatant from the top for tubes 6-10, resuspend pellet in 200 ul wash MEMA and transfer the content from all tubes to ten 14 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using pasteur pipet

15. Again add 200 ul wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions in 14 ml tubes

16. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)

17. Labeling and preparation of dilution tubes and colony dishes

- load 60 mm petri dishes with 4 ml MEMA

- load T-tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.

18. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

19. Centrifuge tubes for 10 min at 2000 rpm, 4°C

20. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 3 cc syringe with 21 gauge needle

1. Determine cell concentration by transferring 100  $\mu$ l to Coulter cup
2. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5 and transfer 0.5 ml to tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3 and vortex tube X.3 and transfer 0.5 ml to tube X.2. Keep tubes on ice.
3. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.
4. Incubate petridishes for 1 week
5. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol. Stain colonies with 0.05% crystal violet
6. Count colonies. There must be between 25 and 250 colonies for the flask to be a valid data point.

mutation assay

Date	Day	Procedure
24/99	0	Plate $10^6$ cells from each tube in P100 E 10ml MEM10
1/99	3	Plate $10^6$ cells from each plate in P100 E 10ml MEM10
11/99	7	Plate $10^6$ cells from each plate in P100 E 10 MEM10
4/99	10	i) count and plate 200,000 cells from each P100 to another P100 X5 E 10ml MEM5 E 10 $\mu$ M 5Gua
		ii) Plate 200 cells in <del>100</del> P60 X3 in MEM5
1/99	17	Colony count

$$600 / .05 = 12,000 / \text{ml}$$

$$\frac{200}{2400000}$$

$$4.8 \times 10^6$$

09/24/99

Cells suspended in 2ml MEMA,  
 For counting, take 100ul in 20 ml

- 589, 598, 571 in 50 ul
- 611, 627, 631
- 541, 559, 561
- 629, 642, 629
- 667, 656, 672
- 542, 561, 559
- 620, 635, 642
- 529, 549, 557
- 607, 598, 622
- 511, 509, 507

Cells suspended in 6 ml MEMA

09/24/99

- 0 411, 431, 435 426 ± 12.9 } 444.5 = 1.00
- 0 471, 461, 459 463 ± 6.4 }
- 3 389, 362, 372 374 ± 13.7 0.84
- 6 332, 321, 341 331 ± 10.0 0.74
- 12 441, 456, 465 454 ± 12.1 1.02
- 0 432, 444, 456 444 ± 12.0 } 433 = 1.00
- 0 409, 422, 436 422 ± 13.5 }
- 3 381, 392, 401 391 ± 10.0 0.90
- 6 356, 365, 369 363 ± 6.66 0.84
- 12 403, 372, 385 387 ± 15.0 0.89

10/01/99

Cells suspended in Gme MEXA

- 772, 761, 756
- 666, 655, 677
- 701, 711, 722
- 656, 631, 634
- 732, 745, 739
- 741, 756, 762 ✓
- 635, 659, 662 ✓
- 672, 657, 659
- 713, 732, 742 ✓
- 699, 710, 729

10/04/99

Cells suspended in Gme of MEXA

- 499, 488, 502
- 436, 456, 462 ✓
- 522, 532, 542 ✓
- 536, 542, 539
- 561, 572, 585
- 437, 452, 462 ✓
- 501, 533, 529
- 490, 471, 478
- 522, 535, 542 ✓
- 531, 555, 563



Plating Efficiency  
200 cells were plated for each tube

10/11/99

Tube #	# of colonies	Avg # of colonies
1	150, 166, 149	147.16
2	137, 129, 152	
3	121, 137, 145	134.33
4	152, 130, 119	133.66
5	117, 125, 129	127
6	165, 147, 155	151
7	141, 159, 139	
8	129, 125, 138	130.6
9	147, 152, 118	139.0
10	167, 145, 149	153.6

Mutant colonies

200,000 cells were plated for each tube

Tube #	# of colonies	Avg. # of colonies	Mutants/Cell
1	1, 2, 1, 1, 1	0.8	0.00000544
2	0, 1, 0, 1, 0		
3	12, 17, 19, 15, 17	16	0.000119
4	38, 45, 29, 50, 39	40.2	0.000300
5	60, 62, 59, 41, 58	54	0.000425
6	2, 0, 0, 1, 0	0.7	0.0000046
7	2, 1, 1, 0, 0		
8	11, 6, 5, 7, 14	8.6	0.000065
9	29, 31, 19, 37, 20	27.2	0.000195
10	36, 39, 47, 40, 48	42	0.000273