

TABLE-3

Expt. #: 3

Date/Time: 06/19/00; 3:00 p.m.

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]	kBq/Cluster [pCi/cell x 148]
1	592, 569, 575	578			
2	601, 623, 617	613			
3	509, 523, 535	522	2089333	0.00118	0.175
4	542, 567, 572	560	2241333	0.003	0.449
5	511, 532, 542	528	2113333	0.009	1.33
6	635, 652, 666	651	2604000	0.015	2.22
7	562, 589, 555	568	2274666	0.0501	7.42
8	611, 589, 599	599	2398666	0.0945	13.99
9	581, 592, 571	581	2325333	0.1309	19.37
10	610, 639, 621	623	2493333	0.1571	23.24