

TABLE-3

Expt. # : 1

Date/Time : 07/20/99 ; 1-30 p.m.

Tube #	Coulter count for 100 ul cell suspension <i>mono. select > 50 µl</i>	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]
1	475, 469, 443	462.3	1849333	0
2	432, 426, 456	438	1752000	0
3	565, 578, 557	566.6	2266666	0.0015
4	479, 469, 441	463	1852000	0.0034
5	469, 482, 486	479	1916000	0.0052
6	130, 118, 118*	122	488000	0.0239
7	462, 457, 463	460.6	1842666	0.0204
8	475, 476, 476	475.6	1902666	0.0371
9	481, 471, 491	481	1924000	0.0632
10	476, 483, 495	484.6	1938666	0.0804

* Cells must have been lost during washing