

9/27 Exp wib5001

9/27 6:45pm Retention pt 2, 75 cm<sup>2</sup> flask from 10:20am 9/25  
 Remove medium, wash 2X 15ml PBS PS  
 Add 1ml trypsin, 37°C, 4min  
 Resusp. in 1.5 ml MEMB \*

Syringe 5X 3cc ZIG  
 Coulter count 100µl, 20µl Isotone, 500µl marker  
 1000 978 968 Bkg 18

Transfer 3x 600µl to 7ml vials & add 5ml Ecoline, vortex

9/28 4:40pm Retention pt 3, 75 cm<sup>2</sup> flask from 10:20am 9/25  
 Remove medium, wash 1X w/ 15ml PBS PS  
 Add 1ml trypsin, 37°C, 4min  
 Resusp. in 1.5 ml MEMB

Syringe 5X 3cc ZIG  
 Coulter count 100µl, 20µl Isotone, 500µl marker  
 3100 2959 2903

Transfer 3x 600µl to 7ml vials & add 5ml Ecoline  
 vortex