

TABLE-3

Expt. # 2 ✓

Date/Time : 02/15/99

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]
1	711, 745, 735			
2	852, 837, 860			
3	811, 799, 789			
4	735, 745, 760			
5	809, 811, 835			
6	799, 805, 811	805	3220000	0.1604
7	760, 741, 752	751	3004000	0.1581
8	715, 721, 726	721	2882666	0.1488
9	760, 775, 745	760	3040000	0.1625
10	801, 812, 823	812	3248000	0.1657

0.1591