

TABLE-3

Expt. # : 3

Date/Time : 05/22/98; 2-15 p.m

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 <sup>6</sup> Cells/ml]
1	321, 317, 316	318	1272000	0
2	354, 320, 355	343	1372000	0
3	356, 352, 346	351	1405333	0.000078
4	370, 350, 355	358	1433333	0.00021
5	324, 322, 316	320	1282666	0.00049
6	326, 330, 318	324	1298666	0.00083
7	328, 318, 312	319	1277333	0.00205
8	309, 327, 323	319	1278666	0.00409
9	327, 320, 319	322	1288000	0.00823
10	314, 311, 330	318	1273333	0.01218