

TABLE-3

Expt. # : 4

Date/Time : 02/20/98 ; 11-10 a.m.

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell ^{ml} [uCi/cell x 10 ⁶ Cells/ml]
1	647, 597, 594, 633, 621	633.66	253466.6	0
2	693, 621, 619, 651	630.33	252133.3	0
3	727, 700, 684, 645	676.33	270533.3	0.0321
4	686, 610, 605, 611	608.66	243466.6	0.0627
5	602, 617, 633	617.33	246933.3	0.1082
6	703, 733, 739	725	290000	0
7	640, 628, 664	644	257600	0.0383 ⁰
8	729, 753, 735, 721	736.33	294533.3	0.0699 ^{0.0335}
9	749, 738, 698	728.33	291333.3	0.0706
10	748, 675, 669, 645	663	265200	0.1342